

LARYVOX® PRO HME

FILTER CASSETTE



LARYVOX® PRO
HME



LARYVOX® PRO
HME HIGHFLOW

FAHL

LARYVOX® PRO HME

FILTER CASSETTE



DISCRETE DESIGN &

EASY SPEAKING OPTION

FAHL

Andreas Fahl
Medizintechnik-Vertrieb GmbH
August-Horch-Str. 4a
51149 Köln - Germany
Phone +49 (0) 22 03 / 29 80-0
Fax +49 (0) 22 03 / 29 80-100
email info@fahl.de
www.fahl.de

LARYVOX® PRO HME - filter cassette and speaking option

LARYVOX® PRO HME (Heat-Moisture Exchanger) provides all functions of a conventional filter cassette - it filters, warms and humidifies the respiratory air. In addition, LARYVOX® PRO HME is equipped with an easy-to-operate speaking option for patients with a voice prosthesis and has an inconspicuous design.

FILTERING - WARMING - HUMIDIFYING

The filter medium integrated in the LARYVOX® PRO HME housing serves for exchanging heat and moisture and filters the respiratory air. Due to the loss of nasal breathing, filtration, warming and humidification of the respiratory air must be ensured in order to prevent irritation of the airways due to dry air, dust or harmful substances contained in the respiratory air.

After surgery, tracheotomised / laryngectomised patients often suffer from increased production of mucous and irritation of the throat associated with an urge to cough. These symptoms can impair breathing and voice rehabilitation. Regular use of HME filter cassettes alleviates the symptoms and enhances the quality of life.

LARYVOX® PRO HME is available in two versions: the HighFlow version with low breathing resistance, which is thus particularly suitable in the case of significant physical strain and for new HME users, and furthermore LARYVOX® PRO HME with normal breathing resistance.

The filter cassette is intended for single use only and must be replaced at the latest after 24 hours. LARYVOX® PRO HME can be worn both during the day and at night.

DISCRETE DESIGN AND HYGIENIC APPLICATION

LARYVOX® PRO HME consists of a transparent plastic housing and a beige closing lid. The housing is equipped with a protective device below the filter insert. This prevents on the one hand that secretions come into contact with the HME filter. On the other hand, the protective device enhances hygienic application by preventing secretions from getting onto the patient's fingers or clothing.

Thanks to the flat profile of LARYVOX® PRO HME and the beige closing lid, the HME filter cassette is very inconspicuous. The simple handling of the filter cassettes facilitates its affixing exactly over the stoma, expectoration, and care of the stoma. Lateral breathing apertures enable unhindered air flow even if the cassette is covered with clothing/tracheostoma protection material.



LARYVOX® PRO HME




LARYVOX® PRO HME HighFlow

The HME filter cassette can be conveniently fixed in place with the aid of any base plate (e.g. LARYVOX® Tape) or tracheostomy tube (e.g. LARYNGOTEC®) with a 22 mm combi-adapter. The LARYVOX® PRO HME filter cassette can easily be removed from the 22 mm holding ring at any time.


EASY SPEAKING OPTION

LARYVOX® PRO HME serves laryngectomised patients with voice prosthesis to achieve a good voice result.

By exerting slight finger pressure on the closing lid, the patient can accomplish a very precise closure of the tracheostoma. Even if the housing lid is not operated exactly on-centre, complete closure can still be easily achieved, thereby avoiding unwanted loss of air which laterally bypasses the HME and reducing unpleasant side noises. The device can be operated even through the patient's clothing.



Application of the speaking option (closing)



Application of the speaking option (opening)



The front ① and side ② openings for the airflow facilitate unhindered and easy breathing!

ORDER INFORMATION:

▶ LARYVOX® PRO HME
pack of 30

REF 49841

▶ LARYVOX® PRO HME HighFlow
pack of 30

REF 49840