LARYVOX® HME FILTERKASSETTE

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PICTROGRAM LEGEND

EN

 REF	Catalogue number		Use by
 LOT	Batch code		Do not reuse
	Consult instructions for use		Content (in pieces)
	Keep away from sunlight		Manufacturer
	Keep dry		

LARYVOX® HME FILTER CASSETTE

I. FOREWORD

These instructions for use are valid for LARYVOX® HME filter cassettes. The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling. **Please read the instructions for use carefully before using the product for the first time!**

II. INTENDED USE

LARYVOX® HME filter cassettes are particle filters consisting of a foam filter and a plastic housing. The LARYVOX® HME filter cassettes are affixed over the tracheostoma either on a special base plate with holding ring or on a tracheostomy tube with combi-adaptor or 22 mm connector. LARYVOX® HME filter cassettes filter the breathing air, thereby reducing intrusion of particles into the patient's airways. The filter cassettes reduce the formation of viscous secretions in the lungs. They store the moisture and warmth of the exhaled air in the filter medium and emit these to the inhaled air again during inspiration.

When used for the first time the products must only be selected, used and inserted by clinicians or medical professionals trained in their use.

III. PRODUCT DESCRIPTION

The LARYVOX® HME filter cassette is available in three different versions:

1. LARYVOX® HME filter cassette (REF 49800)

- Normal respiratory resistance
- For patients directly after laryngectomy
- For situations with normal physical exertion
- Plastic housing with outer diameter of 22 mm
- Compatible with all commonly used 22 mm receptacle systems (such as DURA-VENT® KOMBI tracheostomy tube or LARYVOX® Tape base plate)

2. LARYVOX® HME HighFlow filter cassette (REF 49810)

- Low respiratory resistance
- For new users of HMEs
- For situations with strong physical exertion
- Plastic housing with outer diameter of 22 mm
- Compatible with all commonly used 22 mm receptacle systems (such as DURA-VENT® KOMBI tracheostomy tube or LARYVOX® Tape base plate)

3. LARYVOX® O₂ HME filter cassettes (REF 49802)

- Normal respiratory resistance
- With O₂ connection (the oxygen connecting piece fitted at the side of the housing (5mm ISO plug-in connection) enables a connecting tube to be connected up, thereby allowing the filter cassette to be combined with an oxygen source).
- Ideal for postoperative application
- Plastic housing with an outer diameter of 22 mm and an inner diameter of 15 mm
- Compatible with all commonly used 22 mm receptacle systems (such as DURA-VENT® KOMBI tracheostomy tube or LARYVOX® Tape base plate) as well as with a 15 mm standard connector

IV. WARNING

LARYVOX® HME filter cassettes are single-patient products and are intended for single use only. They must not be cleaned or disinfected, as use of disinfectant can render the filter unusable. Sterilisation and reuse can impair the function and is therefore not permissible either. The LARYVOX® HME filter cassette must under no circumstances be inserted directly into the tracheostoma!

V. COMPLICATIONS

If the particle filter is clogged or if signs of clogging appear, the LARYVOX® HME filter cassette must immediately be removed and disposed of.

VI. PRECAUTIONS

Care must be taken to ensure that the valve opening is not blocked by clothing or the like in order to enable unhindered breathing (risk of respiratory distress).

VII. CONTRAINDICATIONS

- The LARYVOX® HME should be removed if respiratory resistance is too high.
- Use in patients with strong exsiccosis (dehydration), or excessive secretion in the lungs and airways, or strong air leakage flow (exhaled air not flowing through the filter cassette) is not permissible.
- Not for use in patients with impaired consciousness and/or in non-independent persons requiring the help of others, who are not able to remove the LARYVOX® HME filter cassette by themselves in emergency situations.

VIII. INSTRUCTIONS

- Check whether the packaging and the product are in a faultless condition.
- Check the use-by or expiry date. Do not use after this date.
- Fit ❶ the LARYVOX® HME filter cassettes vertically onto the 22- or 15-mm connection of the tracheostomy tube or the base plate, applying gentle pressure.
- Check the respiratory resistance.
- Remove the LARYVOX® HME filter cassettes after a ❷ wearing period of 24 hours at most, or if the filter cassette is contaminated with secretions for instance. To do this, hold the base plate or the tracheostomy tube in position with two fingers ❸ and remove the filter cassette by tilting it downward slightly.

IX. SERVICE LIFE

- The maximum period of use / wearing period is 24 hours.
- The filter must be changed even several times a day if required, for instance if the respiratory resistance is in danger of increasing due to accumulation of secretions.

X. GENERAL TERMS AND CONDITIONS

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the filter cassette or by repairs if these modifications or repairs were not carried out by the manufacturer himself. This applies both to damages to the filter cassette itself caused thereby and to any consequential damages caused thereby.

If the filter cassette is used beyond the period of use specified under Section IX, and/or if the filter cassette is used, maintained (cleaned, disinfected) or stored in non-compliance with the instructions and specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

Product specifications described herein are subject to change without notice.

LARYVOX® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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