
FAHL

## LARYNGOTEC® PRO TRACHEALKANÜLE

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---

No inner cannula

Outer cannula including 22 mm combi-adapter (KOMBI)

Several small holes in a specific arrangement in the center of the tube, is referred to as „sieve“ (LINGO)

Tube holder included

Length short

MRT suitable

Angle

Manufacturer

Use by

Consult instructions for use

Batch code

Catalogue number

Content (in pieces)

Single patient use

Keep dry

Keep away from sunlight

# LARYNGOTEC® PRO TRACHEOSTOMY TUBES

EN

## I. FOREWORD

These instructions for use are valid for Fah!® tracheostomy tubes. The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling of Fah!® tracheostomy tubes.

**Please read the instructions for use carefully before using the product for the first time!**

Keep them in an easily accessible place for future reference.

## II. INTENDED USE

Fah!® tracheostomy tubes are intended to stabilise the tracheostoma following laryngectomy or tracheostomy.

The tracheostomy tube is designed to keep the tracheostoma open.

When used for the first time the products must only be selected, used and inserted by clinicians or medical professionals trained in their use.

**Users must be trained by medical professionals on how to use the Fah!® tracheostomy tubes safely and appropriately.**

Fah!® LINGO tube variants are intended for tracheotomised patients who retain their larynx or for laryngectomees who use a shunt valve (voice prosthesis users).

## III. WARNINGS

**In case of an unstable tracheostoma, the airways must always be secured before removing the tracheostomy tube and a replacement tube must be kept ready for insertion. The replacement tube must be inserted immediately, even before cleaning and disinfecting the changed tube.**

**Patients must be trained by medical professionals on how to use the Fah!® tracheostomy tubes safely.**

**Fah!® tracheostomy tubes must never be blocked, e.g. by mucus or encrustation. Blockage can cause asphyxiation!**

**Mucus in the trachea can be aspirated (removed by suction) through the tracheostomy tube with a tracheal suction catheter.**

**Damaged tracheostomy tubes must not be used and must be discarded immediately. Use of a damaged tube may result in airway compromise.**

**Tracheostomy tubes are single-patient products and intended for use by one patient only.**

**Do not re-use with another patient. Do not reprocess for re-use by another patient.**

**Do not use tracheostomy tubes during laser treatment (laser therapy) or with electrosurgical devices! Contact with a laser beam may damage the tube.**

### **CAUTION!**

**Tracheostomy tubes with speech function are only recommended for tracheotomised patients with normal secretion and normal mucous tissue.**

**If there is strong secretion, a tendency to form granulation tissue, if the patient is undergoing radiation therapy or if scabs have formed, a perforated cannula version is only recommended under regular surveillance by a doctor and provided that shorter replacement intervals are observed (as a rule once a week), because the perforations in the outer tube can increase the formation of granulation tissue.**

## IV. COMPLICATIONS

The following complications could emerge when using this device:

Contamination of the stoma may make it necessary to remove the tube. Contaminations can furthermore cause infections that necessitate the use of antibiotics.

Unintentional inhalation of a tube/button that was not sized correctly may need to be removed by a doctor/physician. If mucus blocks the tube/button it should be removed and cleaned.

## V. CONTRAINDICATIONS

**Do not use if the patient is allergic to the material.**

**Do not use with ventilated patients.**

## VI. PRECAUTIONS

The correct tracheostomy tube size should be selected by the treating physician or by trained medical professionals.

When using Fah!® tracheostomy tubes, the products may become infected for instance with yeast (Candida), bacteria or the like, which may lead to the material becoming cloudy and its service life being shortened. Should this occur, the tube must be replaced immediately.

Use only auxiliary equipment with a 22 mm connection with Fah!® KOMBI adapters to prevent accidental loss of accessories or damage to the tube.

We highly recommend always keeping at least two replacement tubes on hand to ensure continuous supply.

Do not perform any type of adjustment, repair or alteration on a tracheostomy tube or compatible accessories. If a product is damaged, it must be discarded immediately.

## VII. PRODUCT DESCRIPTION

The LARYNGOTEC® PRO tracheostomy tubes are products made of medically pure silicone.

We supply the Fah!® tracheostomy tubes in various sizes and lengths.

Fah!® tracheostomy tubes are reusable medical products intended for single-patient use.

Fah!® tracheostomy tubes are only to be used by the same patient and not by another patient.

The tube tip is rounded to prevent irritation of the mucus membranes in the trachea.

The standard version of the silicone tube is recommended if a normal tracheostomy tube is no longer required, but a laryngectomy tube is required to stabilize the tracheostoma.

The KOMBI versions of the Fah!® tracheostomy tubes have a combi-adapter on the neck flange. It has an opening with a diameter of 22 mm and is thereby compatible with standard filter and valve systems with a 22 mm connector. This allows the attachment of special filters for filtering, humidifying and heating the air.

It may be advisable to use tubes of different lengths alternately so the tip of the tube does not always touch the same spot in the trachea and cause irritation. This prevents pressure points or the formation of granulation tissue. Always discuss the exact procedure with your doctor/physician.

## MRT (MAGNETIC RESONANCE TOMOGRAPHY) INFORMATION

**The Fah!® silicone tube can be worn by the patient even while undergoing magnetic resonance tomography (MRT).**

### Neck flange

A particular feature of the Fah!® tracheostomy tube is the specially shaped neck flange, which is designed to fit the anatomy of the neck.

The neck flange of the tracheostomy tubes has two side eyelets for inserting a tube holder. (see picture 1) Please follow the instructions provided with the tube holder carefully when attaching and removing it from your tracheostomy tube.

Ensure that the Fah!® tracheostomy tube is placed in the tracheostoma without tension and the position is not changed when fastening the tube holder.

## VIII. DIRECTIONS FOR TUBE INSERTION AND REMOVAL

### For the Doctor/Physician

The appropriate tracheostomy tube must be selected by a doctor/physician or trained medical professionals.

Select a tube that fits the patient's anatomy to optimise comfort and ventilation (breathing in and out).

### For the Patient

#### CAUTION!

**Check the use-by or expiry date. Do not use after this date.**

It is advisable to use sterile disposable gloves.

Carefully examine the tube before first use to make sure that it is not damaged and that there are no loose parts.

Should you notice any anomaly or anything unusual, DO NOT use the tube. Return the tube to the manufacturer for inspection.

Tracheostomy tubes must be cleaned thoroughly every time before being inserted. Cleaning is also recommended prior to first use if the tracheostomy tube is not supplied as a sterile product!

The tube must always be cleaned and, if necessary, disinfected as follows before re-inserting according to the instructions provided below.

The neck flange respectively the funnel-shaped housing (retainer ring) must not be inserted into the tracheostoma. Make sure that it is always outside the tracheostoma (see picture 2).

If secretions collect in the lumen of the Fah!® tracheostomy tube and these cannot be removed by coughing or suctioning, the tube or stoma button should be removed and cleaned.

After cleaning and/or disinfection, carefully examine the Fah!® tracheostomy tube for sharp edges, cracks, or other signs of damage, since these may impair function and/or injure the mucus membranes in the airways.

Discard all damaged tubes.

## 1. Insertion of the tube

### Step-by-step instructions to insert Fah!® tracheostomy tubes.

Before application, users should clean their hands (see picture 3).

Remove tube from the package (see picture 4).

If you are inserting the tube yourself, use a mirror to make insertion of the Fah!® tracheostomy tube easier.

When inserting the Fah!® tracheostomy tube, hold it by the neck flange with one hand (see picture 5).

Pull the tracheostoma slightly apart with your free hand to allow the tip of the tube to fit into the tracheostoma more easily.

Now carefully insert the tube into the tracheostoma during the inspiration phase (while breathing in) while tilting your head slightly back (see picture 5).

Advance the tube into the trachea.

Straighten your head once the tube has been inserted further into the trachea.

To facilitate insertion of the tracheostomy tube, it is advisable to lubricate the outer tube by wiping with FAHL® OPTIFLUID® Lubricant Gel (REF 31600, 36105).

Please observe the instructions for use for the respective product.

### CAUTION!

**Only use water-based lubricants!**

**Oil-based lubricants can damage or destroy the product.**

## 2. Removing the tube

### CAUTION!

**Accessories such as a tracheostoma valve or HME (Heat Moisture Exchanger) must be removed first before proceeding to remove the Fah!® tracheostomy tube.**

Proceed very carefully to avoid injury to the mucus membranes.

### Step-by-step instructions to remove the Fah!® tracheostomy tubes:

The tracheostomy tubes should be removed with the head slightly tilted back. Grip the tube or button at the side by the neck flange or the housing (see picture 5).

Carefully remove the tracheostomy tubes.

## IX. CLEANING AND DISINFECTION

### CAUTION!

**For reasons of hygiene and to avoid the risk of infection, you should thoroughly clean the Fah!® tracheostomy tube at least twice a day, more often during heavier secretion production.**

**In case of an unstable tracheostoma, the airways must always be secured before removing the tracheostomy tube and a replacement tube must be kept ready for insertion. The replacement tube must be inserted immediately, even before cleaning and disinfecting the changed tube.**

Keep this in mind, your personal cleaning schedule, including additional disinfection if required, must always be determined in consultation with your doctor/physician and according to your individual needs.

**General cleaning may be insufficient to meet the special hygiene requirements to prevent infections for patients with specific problems (e.g. MRSA, ORSA etc.), who are subject to an increased danger of re-infection. We recommend chemical disinfection of the tubes according to the instructions provided below. Please consult your doctor/physician.**

### CAUTION!

**the tracheostomy tubes must not be cleaned using a dishwasher, a steam cooker, a microwave oven, a washing machine or similar appliances!**

### CAUTION!

**Cleaning agent and disinfectant residues on the tracheostomy tube can lead to irritations of the mucus membranes or other health impairments.**

Tracheostomy tubes are considered as belonging to the category of instruments with hollow cavities. Special care must therefore be taken when performing cleaning or disinfection that the tube is fully wetted with the solution used and is unobstructed (see picture 6).

**The solutions used for cleaning and disinfection must be freshly prepared prior to each use.**

### 1. Cleaning

Fah!® tracheostomy tubes must be cleaned / replaced regularly in accordance with the individual patient's needs.

Only use cleaning agents if the tube is outside the tracheostoma.

A mild, pH-neutral washing lotion can be used for cleaning the tube. We recommend using the special cannula-cleaning powder (REF 31110) as directed by the manufacturer.

Never clean Fah!® tracheostomy tubes with cleaning agents that are not approved by the manufacturer. Never use aggressive household cleaners, detergents, high-concentration alcohol or denture cleaners.

This may be very hazardous to health. The tube may also be destroyed or damaged.

Depending on the individual needs of the patient, which must be determined by the treating physician, it may be necessary to immerse the tracheostomy tube in boiling water (212 degrees Fahrenheit/100 degrees Celsius) for at least 5 minutes but not more than 10 minutes.

#### **WARNING**

**Let the tube/button cool before reinsertion to avoid scalds and burns.**

#### **Cleaning Steps**

Remove any additional items before cleaning.

First, rinse the tube thoroughly under running water (see picture 7).

Use only lukewarm water to prepare the cleaning solution and observe the instructions for use of the cleaning agent.

We recommend the use of the cannula-cleaning tub with sieve insert (REF 31200) to facilitate the cleaning procedure.

To do this, hold the sieve insert by the upper edge to prevent contact with and contamination of the cleaning solution (see picture 8).

Place only one tracheostomy tube at a time in the sieve insert of the cannula-cleaning tub. If more than one tube is cleaned at a time, they may be compressed too strongly and damaged by excessive pressure.

After the immersion time has expired (see instructions for use of the cannula cleaning powder), wash the tracheostomy tube thoroughly several times with lukewarm, clear water (see picture 7). There must be no residues of the cleaning agent on the tube when it is inserted into the tracheostoma.

If necessary, for example, if persistent secretion residues cannot be removed by the cleaning bath, additional cleaning with a special cannula-cleaning brush (OPTIBRUSH®, REF 31850 or OPTIBRUSH® Plus with fibre top, REF 31855) may help. Only use the cleaning brush, if the tube is removed and already outside the tracheostoma.

Always insert the cannula cleaning brush into the cannula from the cannula tip (see picture 9).

Use the brush as directed and proceed with great care to avoid damage to the soft tube material.

Thoroughly rinse the tube or button under lukewarm running water or sterile saline solution (0.9% NaCl solution).

Dry the tube thoroughly with a clean and lint-free cloth after wet cleaning.

Never use tracheostomy tubes with impaired functionality or with damages such as sharp edges or cracks, because this may lead to injuries of the mucus membranes in the trachea. If damages are found, discard the tube.

#### **2. Chemical Disinfection Instructions**

The Fahl® tracheostomy tube can be disinfected by cold disinfection with special chemical disinfectants.

Disinfection should always be done if ordered by a doctor/physician due to specific health concerns caused by disease, infection, or your specific situation.

Disinfecting is generally indicated to prevent cross-infection and in inpatient situations (e.g. hospitals, nursing homes, and/or other health care facilities) to limit infections.

#### **CAUTION!**

**Always clean according to the cleaning procedure described above before proceeding to disinfection (if applicable).**

**Disinfectants that release chlorine or that contain strong alkalis or phenol derivatives must under no circumstances be used. This may seriously damage or even destroy the tube or button.**

For this purpose, OPTICIT® tube disinfectant (REF 31180) should be used in accordance with the manufacturer's instructions.

As alternative, we recommend a disinfectant based on glutaric aldehyde as active ingredient. Always observe the area of application and spectrum of activity specified for the disinfectant by the manufacturer. Follow the instructions for the disinfectant.

Dry the tube thoroughly with a clean and lint-free cloth after wet cleaning.

#### **X. STORAGE**

Currently not used and cleaned tubes should be stored in a clean plastic container, dry environment away from dust, sunlight and heat.

We highly recommend keeping at least two replacement tubes on hand to ensure continuous supply.

#### **XI. SERVICE LIFE**

These tracheostomy tubes are products intended for single-patient use.

When used and looked after correctly, in particular observing the instructions for use, cleaning, disinfection and maintenance provided in the instructions for use, these tracheostomy tubes have an expected service life of three to six months.

The maximum period of use should not exceed 6 months.

The service life of a tube is influenced by many factors. The composition of the secretions, the quality of cleaning and other aspects, for example, are very important.

Damaged tubes must be replaced immediately.

EN

## **XII. LEGAL NOTICES**

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the tracheostomy tube, especially due to shortening of the same or perforations, if these modifications or repairs were not carried out by the manufacturer himself. This applies both to damages to the tubes themselves caused thereby and to any consequential damages caused thereby.

If the tracheostomy tube is used beyond the period of use specified under Section XI, and/or if the tracheostomy tube is used, maintained (cleaned, disinfected) or stored in non-compliance with the instructions and specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious adverse event occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, this is to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

Product specifications described herein are subject to change without notice.

LARYNGOTEC® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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