
HUMIDOPHONE®

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PICTOGRAM LEGEND

en

Pictogrammes listed below you'll find on the product packaging, if applicable.

Catalogue number

Batch code

Use by

Date of manufacture

Content (in pieces)

Consult instructions for use

Keep away from sunlight

Keep dry

Do not reuse

Non-sterile

Medical devices

CE marking

CE marking identification number of the notified body

Manufacturer

HUMIDOPHONE® SPEAKING VALVE

en

1. FOREWORD

These instructions for use are valid for **HUMIDOPHONE®** speaking valves with filter function. When used for the first time, the products must be selected, used and inserted by appropriately trained physicians or medical professionals/medical device consultants.

These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

Please read the instructions for use carefully before using the product for the first time!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the product.

It contains important information on the product!

2. SAFETY INFORMATION

⚠ WARNING!

The HUMIDOPHONE® must only be used by patients with partially or fully preserved larynx!

⚠ WARNING!

The HUMIDOPHONE® is a single-patient product and is intended for single use only.

It must not be cleaned or disinfected, as use of cleaning agent can render the filter unusable. Sterilisation and reuse impair the function and therefore are also not permissible.

⚠ WARNING!

The HUMIDOPHONE® must under no circumstances be inserted directly into the tracheostoma (risk of suffocation!).

⚠ WARNING!

Care must be taken to ensure that the valve opening is not blocked by clothing or the like in order to enable unhindered breathing (risk of respiratory distress).

⚠ WARNING!

The HUMIDOPHONE® must only be used by patients with a clear state of consciousness and spontaneous breathing!

The latter does not apply if the patient is in an institution under nursing supervision for 24 hours or in home care.

⚠ WARNING!

The HUMIDOPHONE® may not be used in conjunction with a nebuliser or atomiser.

The HUMIDOPHONE® may not be used in conjunction with a nebuliser or atomiser.

This would invariably entail the risk of increased respiratory resistance and thus obstruction of breathing.

⚠ WARNING!

If the filter is clogged or if signs of clogging appear, the HUMIDOPHONE® speaking valves must immediately be removed and disposed of.

3. INTENDED PURPOSE

The speaking valves enable spontaneously breathing tracheotomised patients with a partially or fully preserved larynx to speak without having to use their fingers in the presence of a tracheostomy tube with voice option or a tracheostomy patch.

It filters the breathing air, thereby reducing intrusion of particles into the patient's airways.

Speaking valves with an O₂ connection also allow the supply of additional oxygen via the integrated O₂ connection.

4. INDICATION

Speaking valves are indicated for spontaneously breathing tracheotomised patients with a partially or fully preserved larynx.

5. CONTRAINDICATION

⚠ WARNING!

Tracheostomy tubes with a HUMIDOPHONE® must not be used on any account by laryngectomised patients, as this could lead to serious complications, even suffocation!

Use in patients with strong exsiccosis (dehydration), or excessive secretion and water accumulation in the lungs and airways (left-ventricular heart failure, pulmonary oedema), or strong air leakage flow (exhaled air not flowing through the filter cassette) is not permissible.

6. COMPLICATIONS AND ADVERSE REACTIONS

⚠ WARNING!

The HUMIDOPHONE® must be removed before going to sleep to prevent potential obstruction of the airways!

⚠ WARNING!

Prior to using a HUMIDOPHONE®, the cuff of the tracheostomy tube must be completely free of air.

7. PRODUCT DESCRIPTION

The HUMIDOPHONE® is a speaking valve as well as a particle filter, consisting of a foam filter and a plastic housing as well as of a lid with valve function.

The plastic housing is equipped on the patient side with a standardized opening, thereby allowing for connection to a fenestrated/sieved tracheostomy tube.

During inhalation, the valve diaphragm of the HUMIDOPHONE® opens and the inhaled air flows through the foam filter element.

After inhalation, the valve closes, thereby enabling the user to speak.

The pressure rise during expiration supports the valve closure process.

HUMIDOPHONE® - On the patient side, the plastic housing is equipped with a central opening with 22 mm external diameter, thereby allowing for connection to a 22 mm receptacle system.

HUMIDOPHONE® TYP 15 - On the patient side, the plastic housing is equipped with a central opening with 15 mm internal diameter, thereby allowing for connection to a 15 mm standard connector.

HUMIDOPHONE® PLUS - On the patient side, the plastic housing is equipped with a central opening with 15 mm internal diameter, thereby allowing for connection to a 15 mm standard connector.

The external diameter of the opening is 22 mm, and is thus also compatible with all commonly used 22 mm receptacle systems.

HUMIDOPHONE® PLUS O₂- On the patient side, the plastic housing is equipped with a central opening with 15 mm internal diameter, thereby allowing for connection to a 15 mm standard connector.

The external diameter of the opening is 22 mm, and is thus also compatible with all commonly used 22 mm receptacle systems.

An oxygen connecting piece fitted to the side of the housing (5 mm plug-in connection) enables a connecting tube to be connected up, thereby allowing combination with an oxygen source.

8. APPLICATION

Check whether the packaging and the product are in a faultless condition.

Check the use-by or expiry date.

Do not use the product after this date.

Use with a 15 mm connector: mount the **HUMIDOPHONE®** with the central 15 mm opening by means of light pressure on the 15 mm standard connector of the tracheostomy tube ❶.

Use with a 22 mm connector: mount the **HUMIDOPHONE®** by means of light pressure on the housing ring of the tracheostomy tube ❷, until the housing snaps in audibly through the outer fastening rim.

Remove the **HUMIDOPHONE®** by holding the tracheostomy tube firmly with two fingers in order to extract the filter cassette by tugging lightly ❸/❹.

9. HYGIENE INSTRUCTIONS

The **HUMIDOPHONE®** is a single-patient product and is intended for single use only.

They must not be cleaned or disinfected, as use of disinfectant can render the filter unusable.

In addition, bacterial or fungicidal colonisation of the filter material can lead to an increased risk of infection.

Cleaning, disinfection or (re-)sterilisation as well as use longer than 24 hours can impair the safety and function of the product and are therefore not permissible!

10. STORAGE

This product should be stored in a dry environment away from sunlight and/or heat.

11. SERVICE LIFE

The maximum period of use / wearing period is 24 hours ❺.

If necessary, the filter must be changed several times a day, for instance if the respiratory resistance is in danger of increasing due to accumulation of secretions.

12. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

13. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages (in particular for functional deficiencies, injuries, infections, and/or other complications or adverse events) caused by unauthorised product alterations, by repairs not performed by the manufacturer or by improper handling, care (cleaning/disinfection) and/or storage of the products in violation of the provisions of these instructions for use.

This applies – to the extent permitted by law – both to damages to the products themselves caused thereby as well as to any consequential damages caused thereby.

When using the product beyond the period of use specified under Section 11, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.

HUMIDOPHONE® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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