











FAHL

**DURAVENT®**

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
# PICTOGRAM LEGEND

EN

Pictogrammes listed below you'll find on the product packaging, if applicable.

|   |   |   |   |
|---|---|---|---|
|    | No inner cannula                                  |  | Date of manufacture                                   |
|    | With one inner cannula                            |  | Manufacturer  |
|    | With two inner cannulas                           |  | Use by  |
|    | With three inner cannulas                         |  | Consult instructions for use                          |
|    | Inner cannula with 22 mm-connector (KOMBI)        |  | CE marking identification number of the notified body |
|    | Inner cannula with 15 mm-connector (UNI)          |  | Batch code  |
|    | Inner cannula with 15 mm-swivel connector (VARIO) |  | Catalogue number                                      |
|    | Inner cannula with low profile                    |  | Content (in pieces)                                   |
|    | With oxygen support                               |  | Single patient use                                    |
|    | Decannulation plug                                |  | Do not use if package is damaged                      |
|    | Sieving   |  | Keep dry  |
|    | Tube holder included                              |  | Keep away from sunlight                               |
|    | Specially for children                            |  | Medical devices                                       |
|    | Inner cannula fenestrated                         |   |   |
|   | Length short                                      |   |   |
|  | Length medium                                     |   |   |
|  | Length long                                       |   |   |
|  | Length XL   |   |   |
|  | MRT suitable                                      |   |   |
|  | Angle   |   |   |

# DURAVENT<sup>®</sup> TRACHEOSTOMY TUBES

## I. FOREWORD

These instructions for use are valid for Fahl<sup>®</sup> tracheostomy tubes. The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling of Fahl<sup>®</sup> tracheostomy tubes.

**Please read the instructions for use carefully before using the product for the first time!**

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the tracheostomy tube. It contains important information on the product!

## II. INTENDED USE

Fahl<sup>®</sup> tracheostomy tubes are intended to stabilise the tracheostoma following laryngectomy or tracheostomy.

The tracheostomy tube is designed to keep the tracheostoma open.

When used for the first time the products must only be selected, used and inserted by clinicians or medical professionals trained in their use.

Overlong tubes can be indicated particularly in cases of deep-seated tracheal stenoses.

Fahl<sup>®</sup> LINGO tube variants are intended for tracheotomised patients who retain their larynx or for laryngectomees who use a shunt valve (voice prosthesis users).

## III. WARNINGS

**Patients must be trained by medical professionals on how to use the Fahl<sup>®</sup> tracheostomy tubes safely.**

**Fahl<sup>®</sup> tracheostomy tubes must never be blocked, e.g. by mucus or encrustation. Blockage can cause asphyxiation!**

**Mucus in the trachea can be aspirated (removed by suction) through the tracheostomy tube with a tracheal suction catheter.**

**Damaged tracheostomy tubes must not be used and must be discarded immediately. Use of a damaged tube may result in airway compromise.**

**When the tracheostomy tubes are inserted or removed, this can lead to irritations, coughing or slight bleeding. If bleeding persists, consult your doctor immediately!**

**Do not use tracheostomy tubes during laser treatment (laser therapy) or with electrosurgical devices! Contact with a laser beam may damage the tube.**

### CAUTION!

Tracheostomy tubes with speech function are only recommended for tracheotomised patients with normal secretion and normal mucous tissue.

Tracheostomy tubes that contain metal parts must under no circumstances be used during a radiation therapy (radiotherapy), as this can lead to serious skin damage! If you need to wear a tracheostomy tube during a radiation therapy, then only use plastic tracheostomy tubes without any metal parts. In the case of plastic tracheostomy tubes with silver speaking valve, the valve can for instance be removed completely from the tube, including the safety chain, by removing the inner cannula from the outer cannula before starting a radiation therapy.

### CAUTION!

If there is strong secretion, a tendency to form granulation tissue, if the patient is undergoing radiation therapy or if scabs have formed, a perforated cannula version is only recommended under regular surveillance by a doctor and provided that shorter replacement intervals are observed (as a rule once a week), because the perforations in the outer tube can increase the formation of granulation tissue.

## IV. COMPLICATIONS

The following complications could emerge when using this device:

Contamination of the stoma may make it necessary to remove the tube. Contaminations can furthermore cause infections that necessitate the use of antibiotics.

Unintentional inhalation of a tube that was not sized correctly may need to be removed by a doctor/physician. If mucus blocks the tube it should be removed and cleaned.

## V. CONTRAINDICATIONS

**Do not use if the patient is allergic to the material.**

### CAUTION!

**Never under any circumstances use cannula variants without cuff for patients on mechanical ventilation!**

**CAUTION!**

Only use perforated/fenestrated cannula variants during artificial respiration after consultation with the treating doctor/physician.

**CAUTION!**

Tracheostomy tubes with speaking valves must not be used on any account by laryngectomised patients, as this could lead to serious complications, even suffocation!

**VI. PRECAUTIONS**

The correct tracheostomy tube size should be selected by the treating physician or by trained medical professionals.

Use only auxiliary equipment with a 22 mm connection with Fah!® KOMBI adapters to prevent accidental loss of accessories or damage to the tube.

We highly recommend always keeping at least two replacement tubes on hand to ensure continuous supply.

Do not perform any type of adjustment, repair or alteration on a tracheostomy tube or compatible accessories. If a product is damaged, it must be discarded immediately.

**VII. PRODUCT DESCRIPTION**

The DURAVENT® tracheostomy tubes are products that are made of medical-grade polyvinyl chloride.

The tracheostomy tubes are made of thermosensitive medical synthetics that develop their optimal product properties at body temperature.

We supply the Fah!® tracheostomy tubes in various sizes and lengths.

Please see the size table in the appendix.

Fah!® tracheostomy tubes are reusable medical products intended for single-patient use.

Fah!® tracheostomy tubes are only to be used by the same patient and not by another patient.

It may be advisable to use tubes of different lengths alternately so the tip of the tube does not always touch the same spot in the trachea and cause irritation. This prevents pressure points or the formation of granulation tissue. Always discuss the exact procedure with your doctor/physician.

**MRT (MAGNETIC RESONANCE TOMOGRAPHY) INFORMATION**

These tracheostomy tubes can even be worn by the patient during magnetic resonance tomography (MRT).

This only applies to tracheostomy tubes without metal parts though.

In the case of tracheostomy tubes with silver speaking valve, the valve, including the safety chain, can for instance be removed completely from the tube by removing the inner cannula with speaking valve from the outer cannula before being exposed to radiation.

**1. Neck flange**

A particular feature of the Fah!® tracheostomy tube is the specially shaped neck flange, which is designed to fit the anatomy of the neck.

The size specifications are printed on the neck flange.

The neck flange of the tracheostomy tubes has two side eyelets for inserting a tube holder.

All Fah!® tracheostomy tubes with fastening eyelets include a disposable tube holder. The tube holder secures the tracheostomy tube to the neck.

Please follow the instructions provided with the tube holder carefully when attaching and removing it from your tracheostomy tube.

Ensure that the Fah!® tracheostomy tube is placed in the tracheostoma without tension and the position is not changed when fastening the tube holder.

**2. Connectors and adapters**

Connectors and adapters are used to connect compatible tube accessories.

The options available to a patient depend on the clinical condition, such as status post laryngectomy or tracheostomy.

Connectors/adapters are normally permanently attached to the inner cannula. The connector in this case is the universal (15 mm) connector with the aid of which so-called artificial noses or heat and moisture exchangers (HMEs, filters for exchanging heat and moisture) can be attached.

This connector is also available in a special version as 15 mm swivel connector. The swivel version of the 15 mm connector is suitable, for instance when using a respirator tube system, for absorbing the torsional forces that arise in such situations and for relieving strain on the tracheostomy tube and stabilising it in position in such a way that irritations of the mucous membrane in the trachea can thereby be avoided.

The 22 mm combi adapter enables attaching compatible filters and valve systems with a 22 mm connector.

### 3. Cannula tube

The cannula tube borders directly on the neck flange and directs the airflow into the trachea.

The Easy Lock hexagonal lock serves to securely attach inner cannulas and compatible accessories.

The X-ray contrast strip in the tube that runs along the side allows the tube to be depicted radiographically and the position of the tube to be checked.

### 4. Inner cannula

The convenient Easy Lock hexagonal lock serves to securely attach the inner cannula and compatible accessories.

The inner cannulas are, depending on specification, either equipped with a standard attachment or permanently connected to particular adapters/connectors, or they can be connected to detachable accessories such as for instance speaking valves.

The inner cannulas can easily be removed from the outer cannula, thus allowing the air supply to be increased quickly if necessary (for instance in case of dyspnoea).

The connection between inner and outer cannula must first be released by slight lifting (unbuttoning with thumb and index finger). This causes holding prongs/engaging elements of the inner cannula to be pulled off the Easy Lock hexagonal connection of the outer cannula.

Inner cannulas must never be used without outer cannula but must always be affixed to the outer cannula.

#### 4.1 Speaking valves

Tracheostomy tubes with speaking valve (LINGO/PHON) are used after tracheotomy with a complete or partially retained larynx and enable the user to speak.

In the case of tracheostomy tubes with silicone speaking valve, the valve is directly attached to the inner cannula and is included in delivery of these tube versions.

**In tracheostomy tubes with silver speaking valve, the speaking valve can be detached from the inner cannula by pushing it out.**

**In tracheostomy tubes with silicone speaking valve, the speaking valve can be detached from the inner cannula by pulling it off.**

### 5. Decannulation plug

The decannulation plug is included in delivery of the tracheostomy tubes with speaking valve and may only be used by tracheotomised patients who have retained their larynx. It may only be inserted under supervision by a doctor/physician. It allows the tracheostomy tube to be sealed and the air supply to be interrupted for a short time and helps the patient train to regain control of breathing via mouth/nose.

#### CAUTION!

**Contraindicated in laryngectomised patients and patients with chronic obstructive pulmonary disease (COPD)! In such cases, the decannulation plug must not under any circumstances be inserted!**

**The decannulation plug must never under any circumstances be used when the tracheostomy tube is in sealed condition! The decannulation plug may only be used with perforated outer cannula without inner cannula.**

#### CAUTION!

**A decannulation plug is enclosed with the tracheostomy tubes with speaking valve for the purpose of preparing for possible decannulation following temporary tracheotomy. The plug can be used to block the air supply via the tracheostomy tube for a short time so that the patient can get used to breathing via mouth/nose again. Decannulation may only be carried out under medical supervision. The plug may only be inserted upon instruction by the doctor/physician. There is a risk of suffocation! Please also make sure you observe the description of indications provided with the respective product variants/specifications!**

#### CAUTION!

**Make sure to read all accompanying product information, instructions for use, indications and contraindications. Discuss use of the product with your doctor/physician before first use.**

## VIII. DIRECTIONS FOR TUBE INSERTION AND REMOVAL

### For the Doctor/Physician

The appropriate tracheostomy tube must be selected by a doctor/physician or trained medical professionals.

Select a tube that fits the patient's anatomy to optimise comfort and ventilation (breathing in and out).

The inner cannula can be removed at any time to increase air supply or for cleaning. This can for instance be necessary if the cannula is clogged up with secretion residues which cannot be removed by coughing or because no equipment for suctioning off the secretions is available.

## For the Patient

### CAUTION!

**Check the use-by or expiry date. Do not use after this date.**

Carefully examine the tube before first use to make sure that it is not damaged and that there are no loose parts.

Should you notice any anomaly or anything unusual, DO NOT use the tube. Return the tube to the manufacturer for inspection.

Tracheostomy tubes must be cleaned thoroughly every time before being inserted. Cleaning is also recommended prior to first use if the tracheostomy tube is not supplied as a sterile product!

The tube must always be cleaned and, if necessary, disinfected as follows before re-inserting according to the instructions provided below.

If secretion collects in the lumen of the Fah!® tracheostomy tube or stoma button and cannot be removed by coughing or aspiration, the tube should be removed and cleaned.

After cleaning and/or disinfection, carefully examine the Fah!® tracheostomy tube for sharp edges, cracks, or other signs of damage, since these may impair function and/or injure the mucus membranes in the airways.

Never under any circumstances continue using damaged tracheostomy tubes.

### 1. Insertion of the tube

#### Step-by-step instructions to insert Fah!® tracheostomy tubes.

Before application, users should clean their hands (see picture 3).

Remove tube from the package (see picture 4).

Then a tracheal compress is pushed onto the tracheostomy tube.

To increase the lubricity of the tracheostomy tube and thus facilitate insertion into the trachea, it is recommended to wipe the outer tube with an OPTIFLUID® stoma oil wipe (REF 31550), which allows even distribution of the stoma oil on the tube (see Figs. 4a and 4b) or FAHL® OPTIFLUID® Lubricant Gel 20g tube (REF 36100) or FAHL® OPTIFLUID® Lubricant Gel 3g sachet (Ref 36105).

If you are inserting the tube yourself, use a mirror to make insertion of the Fah!® tracheostomy tube easier.

When inserting the Fah!® tracheostomy tube, hold it by the neck flange with one hand (see picture 5).

Pull the tracheostoma slightly apart with your free hand to allow the tip of the tube to fit into the tracheostoma more easily.

Special aids for dilating the tracheostoma are also available. These allow the tracheostoma to be dilated gently and evenly, for instance also in emergency situations with collapsing tracheostoma (see picture 6).

Make sure that the tube is not damaged by friction when using an instrument for assistance.

Now carefully insert the tube into the tracheostoma during the inspiration phase (while breathing in) while tilting your head slightly back (see picture 7).

Advance the tube into the trachea.

Straighten your head once the tube has been inserted further into the trachea.

The tracheostomy tubes should always be attached with a special tube holder. This stabilises the tube and thus ensures that the tracheostomy tube is securely seated in the tracheostoma (see picture 1).

### 2. Removing the tube

#### CAUTION!

**Accessories such as a tracheostoma valve or HME (Heat Moisture Exchanger) must be removed first before proceeding to remove the Fah!® tracheostomy tube.**

#### CAUTION!

If the tracheostoma is unstable, or in emergency situations (puncture/dilation tracheostomy), the tracheostoma can collapse after withdrawal of the tracheostomy tube, thereby impairing air supply. A fresh tracheostomy tube must be kept ready for use in such cases and must be quickly inserted if necessary. A tracheal dilator (REF 35500) can be used for temporarily securing the air supply.

Proceed very carefully to avoid injury to the mucus membranes.

#### Step-by-step instructions to remove the Fah!® tracheostomy tubes:

The tracheostomy tubes should be removed with the head slightly tilted back. Grip the tube or button at the side by the neck flange or the housing (see picture 7).

Carefully remove the tracheostomy tubes.

In the case of the DURAVENT® tracheostomy tubes with silver speaking valve, the speaking valve attachment has a pre-formed rail inside, into which a flat valve plate made of sterling silver with a round valve flap is inserted. The inner cannula can be removed from the outer cannula together with the adapted silver speaking valve. The valve plate does not have to be removed beforehand for this. A chain that is attached to the speaking valve receptacle safeguards the valve against falling down.

For this purpose, the connection between inner and outer cannula must be released by slight lifting (unbuttoning with thumb and index finger). This causes holding prongs/engaging elements of the inner cannula to be pulled off the Easy Lock hexagonal connection of the outer cannula.

For reinsertion, proceed in reverse sequence as described above to reinsert the inner cannula and attach it to the Easy Lock hexagonal connection of the outer cannula.

## IX. CLEANING AND DISINFECTION

### CAUTION!

**For reasons of hygiene and to avoid the risk of infection, you should thoroughly clean the Fahl® tracheostomy tube at least twice a day, more often during heavier secretion production.**

**In case of an unstable tracheostoma, the airways must always be secured before removing the tracheostomy tube and a replacement tube must be kept ready for insertion. The replacement tube must be inserted immediately, even before cleaning and disinfecting the changed tube.**

### CAUTION!

**For reasons of hygiene and to avoid the risk of infection, you should thoroughly clean the Fahl® tracheostomy tube at least twice a day, more often during heavier secretion production.**

Keep this in mind, your personal cleaning schedule, including additional disinfection if required, must always be determined in consultation with your doctor/physician and according to your individual needs.

Regular disinfection is only required if this is indicated from the medical point of view as ordered by the doctor/physician. The reason for this is that the upper airways are not free of microbes in healthy patients either.

**General cleaning may be insufficient to meet the special hygiene requirements to prevent infections for patients with specific problems (e.g. MRSA, ORSA etc.), who are subject to an increased danger of re-infection. We recommend chemical disinfection of the tubes according to the instructions provided below. Please consult your doctor/physician.**

### CAUTION!

**Cleaning agent and disinfectant residues on the tracheostomy tube can lead to irritations of the mucus membranes or other health impairments.**

Tracheostomy tubes are considered as belonging to the category of instruments with hollow cavities. Special care must therefore be taken when performing cleaning or disinfection that the tube is fully wetted with the solution used and is unobstructed (see picture 8).

The solutions used for cleaning and disinfection must be freshly prepared every day.

### 1. Cleaning

Fahl® tracheostomy tubes must be cleaned / replaced regularly in accordance with the individual patient's needs.

Only use cleaning agents if the tube is outside the tracheostoma.

A mild, pH-neutral washing lotion can be used for cleaning the tube. We recommend using the special cannula-cleaning powder (REF 31110) as directed by the manufacturer.

Never clean Fahl® tracheostomy tubes with cleaning agents that are not approved by the manufacturer. Never use aggressive household cleaners, detergents, high-concentration alcohol or denture cleaners.

This may be very hazardous to health. The tube may also be destroyed or damaged.

Tubes can also be cleaned by thermal disinfection at not more than 65 °C. Use clean water with a temperature of not more than 65 °C for this purpose. Make sure that the temperature is kept constant (use a thermometer to monitor the temperature) and never under any circumstances expose the tubes to boiling water for disinfection. This could seriously damage the tracheostomy tube.

### Cleaning Steps

Remove any additional items before cleaning.

The inner cannula must also be removed from the outer cannula.

The inner cannula is to be cleaned in the same manner as the outer cannula.

First, rinse the inner and outer cannula thoroughly under running water (see picture 9).

Use only lukewarm water to prepare the cleaning solution and observe the instructions for use of the cleaning agent.

We recommend the use of the cannula-cleaning tub with sieve insert (REF 31200) to facilitate the cleaning procedure.

To do this, hold the sieve insert by the upper edge to prevent contact with and contamination of the cleaning solution (see picture 10).

Place only one tracheostomy tube at a time in the sieve insert of the cannula-cleaning tub. If more than one tube is cleaned at a time, they may be compressed too strongly and damaged by excessive pressure. You can place the inner and outer cannula side by side for this procedure.

Immerse the sieve insert with the tracheostomy tube components in the prepared cleaning solution.

After the immersion time has expired (see instructions for use of the cannula cleaning powder), wash the tracheostomy tube thoroughly several times with lukewarm, clear water (see picture 9). There must be no residues of the cleaning agent on the tube when it is inserted into the tracheostoma.

If necessary, for example, if persistent secretion residues cannot be removed by the cleaning bath, additional cleaning with a special cannula-cleaning brush (OPTIBRUSH®, REF 31850 or OPTIBRUSH® Plus with fibre top, REF 31855) may help. Only use the cleaning brush, if the tube is removed and already outside the tracheostoma.

Always insert the cannula cleaning brush into the cannula from the cannula tip (see picture 11).

Use the brush as directed and proceed with great care to avoid damage to the soft tube material.

In tracheostomy tubes with speaking valve, the speaking valve must first be detached from the inner cannula. The valve itself must not be cleaned with the brush since it could otherwise be damaged or break off.

Thoroughly rinse the tube or button under lukewarm running water or sterile saline solution (0.9% NaCl solution).

Dry the tube thoroughly with a clean and lint-free cloth after wet cleaning.

Never use tracheostomy tubes with impaired functionality or with damages such as sharp edges or cracks, because this may lead to injuries of the mucus membranes in the trachea. If damages are found, the tube must not be used under any circumstances.

## 2. Chemical Disinfection Instructions

The Fah!® tracheostomy tube can be disinfected by cold disinfection with special chemical disinfectants.

Disinfection should always be done if ordered by a doctor/physician due to specific health concerns caused by disease, infection, or your specific situation.

Disinfection is generally indicated to prevent cross-infection and in inpatient situations (e.g. hospitals, nursing homes, and/or other health care facilities) to limit infections.

### CAUTION!

**Always clean according to the cleaning procedure described above before proceeding to disinfection (if applicable).**

**Disinfectants that release chlorine or that contain strong alkalis or phenol derivatives must under no circumstances be used. This could badly damage or even destroy the tracheostomy tube.**

#### Disinfection steps

For this purpose, OPTICIT® tube disinfectant (REF 31180) should be used in accordance with the manufacturer's instructions.

As alternative, we recommend a disinfectant based on glutaric aldehyde as active ingredient. Always observe the area of application and spectrum of activity specified for the disinfectant by the manufacturer.

After disinfection the tube must be rinsed inside and out very thoroughly with sterile saline solution (NaCl 0.9%) and then dried.

Dry the tube thoroughly with a clean and lint-free cloth after wet cleaning.

### CAUTION!

**Heating to over 65°C, disinfection with boiling water or steam sterilisation are not allowed and lead to damage of the tracheostomy tube.**

## X. STORAGE/CARE

Currently not used and cleaned tubes should be stored in a clean plastic container, dry environment away from dust, sunlight and heat.

**After the inner cannula has been cleaned, disinfected if necessary and then dried, the outer surface of the inner cannula should be lubricated by applying stoma oil (OPTIFLUID® Stoma Oil, 25 ml bottle REF 31525/ Stoma oil wipe REF 31550) or lubricant gel (FAHL® OPTIFLUID® Lubricant Gel 20 g tube REF 36100/ 3 g sachet REF 36105).**

We highly recommend keeping at least two replacement tubes on hand to ensure continuous supply.

## XI. SERVICE LIFE

These tracheostomy tubes are products intended for single-patient use.

When used and looked after correctly, in particular observing the instructions for use, cleaning, disinfection and maintenance provided in the instructions for use, these tracheostomy tubes have an expected service life of three to six months.

The service life of a tube is influenced by many factors. The composition of the secretions, the quality of cleaning and other aspects, for example, are very important.

The maximum service life is not extended by using the tube intermittently (e.g. alternately with other tubes at normal cleaning intervals).

The tracheostomy tubes must be replaced at the latest after 6 months of use or 6 months after the packaging has been opened - if damaged or the like, they must be replaced immediately though - the tracheostomy tubes must not be used any longer and must be disposed of correctly.

#### **CAUTION!**

**Modifications to the tracheostomy tube of any kind, in particular perforations as well as repairs to the tube, may only be carried out by the manufacturer himself or by companies expressly authorised for this in writing by the manufacturer! Inexpertly performed work on tracheostomy tubes can lead to serious injuries.**

## **XII. LEGAL NOTICES**

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the tracheostomy tube, especially due to shortening of the same or perforations, if these modifications or repairs were not carried out by the manufacturer himself. This applies both to damages to the tubes themselves caused thereby and to any consequential damages caused thereby.

If the tracheostomy tube is used beyond the period of use specified under Section XI, and/or if the tracheostomy tube is used, maintained (cleaned, disinfected) or stored in non-compliance with the instructions and specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious adverse event occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, this is to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.


Product specifications described herein are subject to change without notice.

DURAVENT® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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