
## PHON VALVE

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---

en

Pictogrammes listed below you'll find on the product packaging, if applicable.

Speaking valve (PHON)

MRT suitable

Catalogue number

Batch code

Use by

Date of manufacture

Content (in pieces)

Consult instructions for use

Keep away from sunlight

Keep dry

Single patient use

Non-sterile

Medical devices

CE marking

Manufacturer

# PHON VALVE SPEAKING VALVE

en

## 1. FOREWORD

These instructions for use apply to **PHON VALVE** speaking valves.

These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

### **⚠ WARNING!**

**When used for the first time, the products must be selected, used and inserted by appropriately trained physicians or medical professionals/medical device consultants.**

Please read the instructions for use carefully before using the product!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the product.

It contains important information on the product!

## 2. SAFETY INFORMATION

### **⚠ WARNING!**

**Speaking valves are single-patient products and are thus intended exclusively for use by one patient only. They can be used several times by the same patient.**

Re-use of any kind in other patients and thus also reconditioning for re-use in other patients is not allowed.

### **⚠ WARNING!**

**Care must be taken to ensure that the valve opening is not blocked by clothing or the like in order to enable unhindered breathing (risk of respiratory distress).**

### **⚠ WARNING!**

**Under no circumstances may speaking valves be inserted directly into the tracheostoma.**

## 3. INTENDED PURPOSE

The speaking valves enable spontaneously breathing tracheotomised patients with a partially or fully preserved larynx to speak without having to use their fingers in the presence of a phonation cannula in combination with a fenestrated inner cannula or a tracheostomy patch.

## 4. INDICATION

Speaking valves are indicated for spontaneously breathing tracheotomised patients with a partially or fully preserved larynx.

## 5. CONTRAINDICATION

The **PHON VALVE** speaking valve must only be used by patients with partially or fully preserved larynx!

### **⚠ WARNING!**

**Tracheostomy tubes with a PHON VALVE speaking valve must not be used on any account by laryngectomised patients, as this could lead to serious complications, even suffocation!**

The **PHON VALVE** speaking valve must only be used by patients with a clear state of consciousness and spontaneous breathing!

The speaking valve should be removed if respiratory resistance is too high.

The **PHON VALVE** speaking valve must not be used in case of increased secretion formation (hypersecretion).

The **PHON VALVE** speaking valve must not be used in case of pronounced vocal cord paralysis (recurrent paresis).

## 6. COMPLICATIONS AND ADVERSE REACTIONS

### **⚠ WARNING!**

**The PHON VALVE speaking valve must be removed before going to sleep to prevent potential obstruction of the airways!**

### **⚠ WARNING!**

**Prior to using a PHON VALVE speaking valve, the cuff of the tracheostomy tube must be completely free of air.**

Voice formation is only possible when using a phonation cannula in combination with a fenestrated inner cannula.

## 7. PRODUCT DESCRIPTION

When breathing in, the valve membrane of the **PHON VALVE** speaking valve opens.

After inhalation, the valve closes, thereby enabling the user to speak.

The pressure rise during expiration supports the valve closure process.

The **PHON VALVE** speaking valve is available in a number of different versions:

### **PHON VALVE speaking valve (REF 27140)**

The **PHON VALVE** speaking valve with valve flap for tracheostomy tubes consists of a plastic housing.

The silicone valve flap is attached at the centre of the lid and closes the speaking valve in resting condition.

The plastic housing is equipped on the patient side with a central opening with 15 mm internal diameter, thereby enabling connection to a 15 mm standard connector.

### **TRACHLINE® PHON VALVE speaking valve (REF 27141)**

The **TRACHLINE® PHON VALVE** speaking valve consists of a plastic housing as well as a lid with valve function.

The valve flap made of silicone is attached at the centre of the lid and closes the speaking valve in resting condition.

It is equipped on the patient side with a bayonet closure and thereby ensures connection to a **TRACHLINE®** tracheostomy tube.

The speaking valve is intended for use only and exclusively in combination with a phonation cannula.

## 8. APPLICATION

Check the use-by or expiry date.

Do not use the product after this date.

Take the valve carefully out of the packaging, clean and dry it in accordance with the instructions (9. Hygiene Instructions) prior to first use.

The user must check the **PHON VALVE** speaking valve for completeness, absence of damage and proper functioning each time before inserting the speaking valve again.

If the **PHON VALVE** speaking valve is incomplete, damaged or no longer functioning properly, it must not be used anymore and must be disposed of immediately.

**⚠ WARNING!**

**Before using the speaking valve, attach the enclosed warning label to the blocking tube of the tracheostomy tube!**

**CAUTION!**

To avoid injury, it is advisable to first connect the speaking valve to the tracheostomy tube (❶, ❷) and only then to insert the tracheostomy tube with attached speaking valve into the tracheostoma.

Otherwise, excessive pressure could be exerted on the mucous membranes via the cannula tube when the speaking valve is put into place.

**PLEASE NOTE!**

Do not use excessive force when connecting the speaking valve to the tracheostomy tube to avoid damage to the housing.

**TRACHLINE® PHON VALVE (REF 27141)**

Place the speaking valve with the bayonet lock on the 15 mm standard connector of the tracheostomy tube ❷.

Fixate the speaking valve in the tracheostomy tube with a slight anticlockwise rotation (on the patient's side when inserted), thus closing the bayonet lock ❸.

**9. HYGIENE INSTRUCTIONS**

The **PHON VALVE** speaking valve must be cleaned thoroughly at least twice a day, and in situations with strong secretion production correspondingly more often, to prevent functional impairment of the speaking valve by adhering residues of secretions and to minimise the risk of infection.

**Cleaning Steps**

Carefully remove the speaking valve by twisting and at the same time pulling it out of the adapter of the tracheostomy tube.

**REF 27141 TRACHLINE® PHON VALVE**

To remove the speaking valve, loosen the bayonet lock with a slight clockwise rotation (on the patient's side when inserted)❹.

Place the speaking valve to soak in clear, lukewarm water for approx. 10 minutes for cleaning, or use a mild, neutral-pH washing lotion.

Then rinse under flowing water and allow to dry in air.

**PLEASE NOTE!**

Do not use any sharp-edged cleaning utensils!

**PLEASE NOTE!**

Only clean the **PHON VALVE** speaking valve with cleaning agents approved by the manufacturer.

**PLEASE NOTE!**

Disinfection and sterilisation of the **PHON VALVE** speaking valve is not permitted.

**PLEASE NOTE!**

Heating, disinfection with boiling water or steam sterilisation are not permitted and lead to damage of the speaking valve.

## 10. STORAGE

Speaking valves should be stored in a dry environment away from sunlight and/or heat.

## 11. SERVICE LIFE

The maximum period of use is 29 days.

The service life of the **PHON VALVE** is influenced by many factors.

Thus, the composition of secretions, conscientious cleaning and storage as well as other aspects can be of decisive significance.

A speaking valve showing damage of any kind, even if only minor, must be replaced immediately.

## 12. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

## 13. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies both to damage caused to the products themselves as a result thereof and to all consequential damage caused thereby.

When using the product beyond the period of use specified under Section 11, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.

**TRACHLINE®** is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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