LARYVOX® PLUG

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PIKTOGRAMM-LEGENDE

EN

	Catalogue number
	Batch code
	Use by
	Single patient use
	Consult instructions for use
	Content (in pieces)
	Keep dry
	Keep away from sunlight
	Manufacturer

I. FOREWORD

These instructions for use are valid for the LARYVOX® PLUG sealing plug. These instructions for use serve to inform the physician, medical professional staff and user in order to ensure technically correct handling.

Please read the instructions for use carefully before using the product for the first time!

Keep them in an easily accessible place for future reference.

II. INTENDED USE

The LARYVOX® PLUG sealing plug serves for the transient use with a defective shunt valve to prevent transprosthetic fluid leakage (centrally through the shunt valve).

The LARYVOX® PLUG (①a) is supplied non-sterile with an insertion rod (①b).

When used for the first time the products must only be selected, used and inserted by clinicians or medical professionals trained in their use.

III. WARNINGS

Patients must be trained by medical professionals on how to use the LARYVOX® PLUG safely.

The LARYVOX® PLUG is not intended for permanent use and is only used for the transient sealing of a defective shunt valve. It is not a permanent replacement for a functioning shunt valve. Please consult a specialist as soon as possible to have the shunt valve professionally examined and replaced if necessary. LARYVOX® PLUG is a single-patient product and thus intended exclusively for use in a single patient.

Re-use of any kind in other patients and thus also reconditioning for re-use in other patients is not allowed.

Only clean the LARYVOX® PLUG with cleaning agents approved by the manufacturer.

The LARYVOX® PLUG must be cleaned after each use according to the cleaning instructions (see Section IX. CLEANING).

Sterilisation of the LARYVOX® PLUG is not possible.

The safety strap on the LARYVOX® PLUG must not be changed or removed under any circumstances.

Never force the LARYVOX® PLUG or insert it too deeply into the shunt valve. There is a risk that the shunt valve may be pushed/displaced into the oesophagus. In this case, a physician must be consulted immediately.

Please also refer to the instructions for use of the shunt valve you are currently using.

IV. COMPLICATIONS

If there are signs of shunt valve leakage or LARYVOX® PLUG malfunction, a physician should be consulted as soon as possible to check the sealing plug and the shunt valve.

V. CONTRAINDICATIONS

Do not use if the patient is allergic to the material.

The LARYVOX® PLUG must not be used with a 4mm size shunt valve.

CAUTION!

No tracheostoma accessories (e.g. HME filter cassettes, base plates, tracheostomy tubes) may be used while using the LARYVOX® PLUG. Use could result in an accidental change of position of the shunt valve/LARYVOX® PLUG and cause aspiration (inhalation) into the lungs.

VI. PRECAUTIONS

Do not perform any type of adjustment or repair on the LARYVOX® PLUG. In case of damage, the product may not be used and must be disposed of immediately.

VII. PRODUCT DESCRIPTION

The LARYVOX® PLUG sealing plug is a product made of medically pure silicone for use with a shunt valve. The LARYVOX® PLUG is used to transient seal a shunt valve to prevent fluid leakage through the shunt valve until the shunt valve has been properly inspected and replaced if necessary.

The LARYVOX® PLUG consists of a safety strap ②a and a sealing plug ②b.

CAUTION!

Using the LARYVOX® PLUG in the shunt valve makes speech temporarily impossible.

VIII. INSTRUCTIONS

CAUTION!

Never use the LARYVOX® PLUG if it displays visible deficiencies or faults. Contact the manufacturer in such cases.

The LARYVOX® PLUG must be thoroughly cleaned before each use, also before first use.

CAUTION

Check the use-by or expiry date. Do not use the product after this date.

Always clean your hands with soap and water before using the LARYVOX® PLUG.

Insertion of the LARYVOX® PLUG

Start use of the LARYVOX® PLUG in front of a mirror under optimal lighting conditions so that the shunt valve is clearly visible. Then insert the LARYVOX® PLUG as follows:

Guide the tip of the insertion rod into the rear opening of the LARYVOX® PLUG ③.

Fasten the safety strap to the insertion rod (④ / ⑤ a and b).

Hold the insertion rod with the safety strap facing downwards ⑥ and carefully insert the LARYVOX® PLUG into the shunt valve until it is fixated correctly ⑦.

CAUTION!

Never force the LARYVOX® PLUG or insert it too deeply into the shunt valve. There is a risk that the shunt valve may be pushed into the oesophagus.

Check the secure placement of the LARYVOX® PLUG by slightly pulling the safety strap ⑧.

Now loosen the safety strap from the insertion rod ⑨.

Slowly and carefully remove the insertion rod from the opening of the LARYVOX® PLUG ⑩ while rotating it slightly.

CAUTION!

The insertion rod must be removed with special care to avoid dislocation (displacement) / aspiration (inhalation) of the shunt valve into the lungs.

Medical assistance must be sought immediately in the event of aspiration.

If necessary, secure the safety strap with medical adhesive tape ⑪.

Removing the LARYVOX® PLUG

Always clean your hands with soap and water before removing the LARYVOX® PLUG.

Remove the LARYVOX® PLUG in front of a mirror under optimal lighting conditions so that the shunt valve and the sealing plug are clearly visible. Then remove the LARYVOX® PLUG as follows:

Carefully pull the safety strap to remove the LARYVOX® PLUG from the shunt valve ⑫.

IX. CLEANING AND DISINFECTION

Clean the LARYVOX® PLUG after each use.

Keep this in mind, your personal cleaning schedule, including additional disinfection if required, must always be determined in consultation with your doctor/physician and according to your individual needs.

General cleaning may be insufficient to meet the special hygiene requirements to prevent infections for patients with specific problems (e.g. MRSA, ORSA etc.), who are subject to an increased danger of re-infection. We recommend chemical disinfection of the tubes according to the instructions provided below. Please consult your doctor/physician.

CAUTION!

Cleaning agent and disinfectant residues on the product can lead to irritations of the mucous membranes or other health impairments.

1. Cleaning

Under no circumstances clean the LARYVOX® PLUG with cleaning agents not approved by the manufacturer. Never use aggressive household cleaning agents, highly concentrated alcohol or denture cleaning agents.

This may be very hazardous to health! The LARYVOX® PLUG may also be destroyed or damaged.

Cleaning Steps

Place LARYVOX® PLUG to soak in clear, lukewarm water for approx. 2-3 minutes for cleaning, or use a mild, neutral-pH washing lotion. Then rinse with plenty of clear water and allow to dry in air.

Do not use any sharp-edged cleaning utensils!

Depending on the individual needs of the patient, which must be determined by the treating physician, it may be necessary to immerse the LARYVOX® PLUG in boiling water (100 degrees Celsius resp. 212 degrees Fahrenheit) for at least 5 minutes but not more than 10 minutes.

Dry the product thoroughly with a clean and lint-free cloth after wet cleaning.

2. Chemical Disinfection Instructions

Cold-disinfection of the LARYVOX® PLUG with special chemical disinfectants is possible.

Disinfection should always be done if ordered by a doctor/physician due to specific health concerns caused by disease, infection, or your specific situation.

Disinfecting is generally indicated to prevent cross-infection and in inpatient situations (e.g. hospitals, nursing homes, and/or other health care facilities) to limit infections.

CAUTION

Always clean according to the cleaning procedure described above before proceeding to disinfection (if applicable).

Disinfectants that release chlorine or that contain strong alkalis or phenol derivatives must under no circumstances be used. This could badly damage or even destroy the LARYVOX® PLUG.

Disinfection steps

For this purpose, OPTICIT® tube disinfectant (REF 31180) should be used in accordance with the manufacturer's instructions. As alternative, we recommend a disinfectant based on glutaric aldehyde as active ingredient. Always observe the area of application and spectrum of activity specified for the disinfectant by the manufacturer.

After disinfection, the LARYVOX® PLUG must be rinsed very thoroughly with sterile saline solution (NaCl 0.9%) and then dried.

Follow the instructions for the disinfectant.

Dry the product thoroughly with a clean and lint-free cloth after wet cleaning.

X. STORAGE

Cleaned products that are not currently in use should be stored in a dry environment in a clean plastic bag and protected from dust, sunlight and/or heat.

XI. SERVICE LIFE

The LARYVOX® PLUG is a single-patient product.

When used and looked after correctly, in particular observing the instructions for use, cleaning, disinfection and maintenance provided in the instructions for use, the LARYVOX® PLUG has an expected service life of up to 12 months.

The maximum period of use should not exceed 12 months.

XII. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or other adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the LARYVOX® PLUG or by repairs, if these modifications or repairs were not carried out by the manufacturer himself. This applies both to damages to the LARYVOX® PLUG caused thereby and to any consequential damages caused thereby.

When used beyond the period of use specified under Section V, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the LARYVOX® PLUG in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice. LARYVOX® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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