
COMBIPHON®

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PICTROGRAM LEGEND

en

Pictogrammes listed below you'll find on the product packaging, if applicable.

	With oxygen support
	Speaking valve (PHON)
	MRT suitable
	Catalogue number
	Batch code
	Use by
	Date of manufacture
	Content (in pieces)
	Consult instructions for use
	Keep away from sunlight
	Keep dry
	Single patient use
	Non-sterile
	Medical devices
	Do not use if package is damaged
	CE marking
	CE marking identification number of the notified body
	Manufacturer

COMBIPHON® SPEAKING VALVE

1. FOREWORD

These instructions for use are valid for **COMBIPHON®** speaking valves.

These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

⚠ WARNING!

When used for the first time, the products must be selected, used and inserted by appropriately trained physicians or medical professionals/medical device consultants.

Please read the instructions for use carefully before using the product!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the product.

It contains important information on the product!

2. SAFETY INFORMATION

⚠ WARNING!

Speaking valves are single-patient products and are thus intended for multiple use by one patient only.

Re-use of any kind in other patients and thus also reconditioning for re-use in other patients is not allowed.

⚠ WARNING!

Care must be taken to ensure that the airflow is not blocked by clothing or the like (risk of respiratory distress) to enable unhindered breathing.

⚠ WARNING!

There is a fire hazard in connection with oxygen devices and oxygen therapies.

Do not use in the vicinity of sparks or open flames.

⚠ WARNING!

Do not grease connections, tubes or other accessories in order to avoid fire hazard.

⚠ WARNING!

Smoking during oxygen therapy is dangerous and can lead to severe fire-induced injury.

⚠ WARNING!

Only use lotions or ointments labelled as oxygen-compatible to avoid the risk of fire or burns.

⚠ WARNING!

The product must never under any circumstances be inserted directly into the tracheostoma!

3. INTENDED PURPOSE

The speaking valves enable spontaneously breathing tracheotomised patients with a partially or fully preserved larynx to speak without having to use their fingers if a phonation cannula in combination with a fenestrated inner cannula or a tracheostomy patch is present.

Speaking valves with an O₂ connection also allow for the supply of additional oxygen via the integrated O₂ connection.

4. INDICATION

Speaking valves are indicated for spontaneously breathing tracheotomised patients with a partially or fully preserved larynx.

5. CONTRAINDICATION

The COMBIPHON® speaking valve must only be used by patients with partially or fully preserved larynx!

⚠ WARNING!

Tracheostomy tubes with COMBIPHON® speaking valve must never under any circumstances be used by laryngectomised patients (patients without larynx), as this could lead to serious complications, and even to suffocation!

The COMBIPHON® speaking valve must only be used by awake and alert patients with a spontaneous breathing!

The speaking valve should be removed if respiratory resistance is too high.

The **COMBIPHON® SLIM** speaking valve must not be used in case of increased secretion formation (hypersecretion).

The **COMBIPHON® SLIM** speaking valve must not be used in case of pronounced vocal cord paralysis (recurrent laryngeal nerve paralysis).

6. COMPLICATIONS AND ADVERSE REACTIONS

⚠ WARNING!

The COMBIPHON® speaking valve must be removed before going to sleep to prevent potential obstruction of the airways!

⚠ WARNING!

Before using a COMBIPHON® speaking valve, the cuff of the tracheostomy tube must be completely deflated and void of air.

Voice formation is only possible when using a phonation cannula in combination with a fenestrated inner cannula.

7. PRODUCT DESCRIPTION

When breathing in, the valve membrane of the **COMBIPHON®** speaking valve opens.

After breathing in, the valve closes, thereby enabling the user to speak.

The pressure rise during expiration supports the valve closure process.

The **COMBIPHON®** speaking valve is available in a number of different versions:

COMBIPHON® SLIM speaking valve (REF 27108)

The **COMBIPHON® SLIM** speaking valve with valve flap for tracheostomy tubes consists of a small, transparent plastic housing.

The valve flap is attached at the side of the lid and closes the speaking valve in resting condition.

The plastic housing is equipped on the patient side with a central opening with 15 mm internal diameter, thereby allowing for connection to a 15 mm standard connector.

COMBIPHON® SLIM O₂ speaking valve (REF 27109)

The **COMBIPHON® SLIM O₂** speaking valve is additionally equipped with an oxygen connecting piece at the side of the housing (5 mm plug-in connection) to which a connecting tube can be connected up, thereby allowing for combination with an oxygen source.

COMBIPHON® speaking valve with multi-adapter (REF 27131)

The **COMBIPHON®** speaking valve with multi-adapter is a speaking valve with silicone membrane.

It consists of a plastic housing and a lid with valve function.

The valve flap is attached at the centre of the lid and closes the speaking valve in resting condition.

The plastic housing is equipped on the patient side with a central opening with 15 mm internal diameter, thereby allowing for connection to a 15 mm standard connector.

The external diameter of the end piece/connector is 22 mm, and is thus also compatible with all commonly used 22 mm receptacle systems.

COMBIPHON® O₂ speaking valve with multi-adapter (REF 27132)

The **COMBIPHON® O₂** speaking valve with multi-adapter is additionally equipped with an oxygen connecting piece at the side of the housing (5 mm ISO plug-in connection) to which a connecting tube can be connected up, thereby allowing for combination with an oxygen source.

The speaking valve is intended for use only and exclusively in combination with a phonation cannula.

8. APPLICATION

Check the use-by/expiry date.

Do not use the product after this date.

Carefully take the valve out of the packaging; clean and dry it prior to first use in accordance with the instructions (9. Hygiene Instructions).

The user must check the **COMBIPHON**[®] speaking valve for completeness, absence of damage and proper functioning each time before inserting the speaking valve again.

If the **COMBIPHON**[®] speaking valve is incomplete, damaged or no longer functioning properly, it must not be used anymore and must be disposed of immediately.

⚠ WARNING!

Before using the speaking valve, attach the enclosed warning label to the blocking tube of the tracheostomy tube!

CAUTION!

To avoid injury, it is advisable to first connect the speaking valve to the tracheostomy tube (❶, ❷) and only then to insert the tracheostomy tube with attached speaking valve into the tracheostoma.

Otherwise, excessive pressure could be exerted on the mucous membranes via the cannula tube when the speaking valve is snapped into place.

PLEASE NOTE!

Do not use excessive force when connecting the speaking valve to the tracheostomy tube to avoid damage to the housing.

For additional administration of oxygen, a connecting tube can be connected to the O₂ connection.

To this end, push the connecting tube onto the O₂ connection.

9. HYGIENE INSTRUCTIONS

The **COMBIPHON**[®] speaking valve must be cleaned thoroughly at least twice a day, and in situations with strong secretion production correspondingly more often, to prevent functional impairment of the speaking valve by adhering residues of secretions and to minimise the risk of infection.

Cleaning Steps

Carefully remove the speaking valve by twisting and at the same time pulling it out of the adapter of the tracheostomy tube.

Place the speaking valve to soak in clear, lukewarm water for approx. 10 minutes for cleaning, or use a mild, neutral-pH cleaning lotion.

Then rinse thoroughly under flowing water while performing kneading movements and allow to dry in air.

PLEASE NOTE!

Do not use any sharp-edged cleaning utensils!

PLEASE NOTE!

Only clean the **COMBIPHON**[®] speaking valve with cleaning agents approved by the manufacturer.

PLEASE NOTE!

Disinfection and sterilisation of the **COMBIPHON**[®] speaking valve is impermissible.

PLEASE NOTE!

Heating, disinfection with boiling water or steam sterilisation are not permitted and lead to damage of the speaking valve.

10. STORAGE

This product should be stored in a dry environment away from sunlight and/or heat.

11. SERVICE LIFE

The maximum period of use is 29 days.

The service life of the **COMBIPHON**[®] speaking valve is influenced by many factors.

Thus, the composition of secretions, conscientious cleaning and storage as well as other aspects can be of decisive significance.

A speaking valve showing damage of any kind, even if only minor, must be replaced immediately.

12. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

13. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies – to the extent permitted by law – both to damages to the products themselves caused thereby as well as to any consequential damages caused thereby.

When using the product beyond the period of use specified under Section 11, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.


Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.

COMBIPHON[®] is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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