





FAHL® MULTI PLUG

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I. FOREWORD

These instructions for use are valid for FAHL® MULTI PLUG. The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling. **Please read the instructions for use carefully before using the product for the first time!**

II. INTENDED USE

FAHL® MULTI PLUG is used for sealing a tracheostomy tube/base plate and enables safe decannulation.

Decannulation is in many cases the therapy objective in the scope of treatment of patients with tracheostomy tube.

If the patient is weaned gradually from the tracheostomy tube without any complications, decannulation, with the result of final closure of the tracheostoma, can be planned.

Decannulation can be carried out if the primary cause for creating the tracheostoma no longer persists.

III. WARNING

FAHL® MULTI PLUG may only be inserted upon instruction by the doctor/physician. There is a risk of suffocation! FAHL® MULTI PLUG is a single-patient product and is thus intended exclusively for one-time or repeated use in one individual patient.

FAHL® MULTI PLUG must under no circumstances be inserted directly into the tracheostoma!

FAHL® MULTI PLUG cannot be used in conjunction with TRACHEOTEC®, TRACHEOTEC® PRO/SILC and TRACHEOSILC®!

IV. CONTRAINDICATIONS

FAHL® MULTI PLUG must only be used by patients with partially or fully preserved larynx!

Classic contraindications of decannulation include tracheal stenoses and tracheomalacia.

The decannulation plug FAHL® MULTI PLUG must under no circumstances be used for laryngectomised patients and patients with chronic obstructive pulmonary disease (COPD)!

The decannulation plug must never under any circumstances be used when the tracheostomy tube is in blocked condition!

The FAHL® MULTI PLUG decannulation plug may only and exclusively be used with a sieved outer cannula combined with a fenestrated inner cannula.

V. COMPLICATIONS

The following complications can occur prior to, during or after decannulation and can impede decannulation or render decannulation impossible: granulations, aspiration, infections as well as strong hypersalivation and atelectases.

VI. PRODUCT DESCRIPTION

The plastic housing is equipped on the patient side with a central opening with 15 mm internal diameter, thereby allowing for connection to a tracheostomy tube with a 15 mm standard connector ③.

The external diameter of the connector is 22 mm, and is thus also compatible with all commonly used 22 mm receptacle systems ④.

If Laryvox® Seal is used with a base plate, the area around the tracheostoma is sealed air-tight. If Laryvox® Seal is used with a tracheostomy tube, the tracheostomy tube is sealed air-tight.

FAHL® MULTI PLUG has a bar for easier handling during insertion and removal in the tracheostoma.

An opening integrated in the bar allows the attachment of a strap to protect against loss of the FAHL® MULTI PLUG, e.g. in case of a coughing fit.

Compatible with all commonly used 15 mm - and 22 mm receptacle systems (such as DURAVENT® KOMBI tracheostomy tube or LARYVOX® Tape base plate).

VII. INSTRUCTIONS

If the relevant indications are fulfilled, the decannulation process can be started.

CAUTION!

Continuous monitoring of the vital parameters, particularly in the initial phase of the decannulation process, is required.

1. At first, the tracheostomy tube is de-blocked under supervision for short de-blocking phases (15 - 30 minutes). To do this, the air is removed from the cuff of the tracheostomy tube using a syringe ①. A suctioning device should always be available, ready for use. De-blocking in the scope of logopedic therapy would be recommendable.

CAUTION!

Decannulation may only be carried out under medical supervision.

2. Check by manually sealing the tracheostomy tube, whether it is possible to breathe out and breathe in through mouth and nose. An endoscopic check of the act of swallowing is advisable at this point.
3. If, after several de-blocking cycles, the new air flow control is not associated with any complications, i.e. there are no signs of aspiration of saliva, then the de-blocking phases are increased to up to 24 hours and a speaking valve (COMBIPHON® O2) is used ②. If necessary, an inner cannula with an appropriate adapter must be inserted into the tracheostomy tube to allow the speaking valve to be used. After several swallowing attempts, the risk of aspiration must be further assessed.

CAUTION!

If complications occur, a new decannulation attempt can only be started afresh after treatment of these complications (Point 1).

4. If there are no signs of aspiration, the tracheostomy tube with cuff can be replaced by a tracheostomy tube without cuff or a tracheostomy tube with phonation device.
5. If there is no infection, and if the spontaneous swallowing rate is satisfactory, the cough reflex and/or deliberate coughing is preserved and the respiratory situation is stable, plugging/sealing of the tracheostomy tube with FAHL® MULTI PLUG can be started ④. Before the tracheostomy tube is sealed, an endoscopic examination of the trachea is mandatory in order to rule out stenoses or granulations.

CAUTION!

The FAHL® MULTI PLUG decannulation plug may only be inserted upon instruction by the doctor/physician.

There is a risk of suffocation!

If the sealing cap is used, it must urgently be ensured that this is only used if a sieved tracheostomy tube with fenestrated inner cannula has been inserted and the low-pressure cuff of the tube has been de-blocked!

6. If no problems arise when the tracheostomy tube is plugged/sealed, the patient can be decannulated. To do this, the tracheostomy tube is removed from the tracheostoma after ruling out the indication ④ and the Laryvox® Tape Extra Fine base plate is fixed in place on the tracheostoma as described in the instructions and is then sealed with FAHL® MULTI PLUG ⑤. Surveillance of the decannulation process should at first be performed with 24-hour monitoring of the vital parameters.
7. The duration of closure is patient-specific and must be assessed and determined by the treating physician/therapist. If the patient tolerates decannulation over this period without any problems, the base plate with FAHL® MULTI PLUG can be removed and the tracheostoma closed individually (surgically/spontaneously).

VIII. SERVICE LIFE

The service life depends on the individual frequency of use and on individual handling.

The products must be replaced immediately when damaged.

IX. CLEANING

Place FAHL® MULTI PLUG to soak in clear, lukewarm water for approx. 2-3 minutes for cleaning, or use a mild, neutral-pH washing lotion.

Then rinse with plenty of clear water and allow to dry in air.

Do not use any sharp-edged cleaning utensils!

X. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH accepts no liability for damage caused by unauthorised product modifications, repairs not carried out by the manufacturer or by improper handling, care (cleaning/disinfection) and/or storage of the products contrary to the provisions of these instructions for use.

Should a serious adverse event occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, this is to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Product specifications described herein are subject to change without notice.

FAHL® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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