


















## DURACUFF® TRACHEALKANÜLE

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# PICTROGRAM LEGEND

Pictogrammes listed below you'll find on the product packaging, if applicable.

	No inner cannula		Angle
	With one inner cannula		Outer diameter O. D.
	With two inner cannulas		Inner diameter I. D.
	Outer cannula including 15 mm-swivel connector (VARIO)		Date of manufacture
	Inner cannula with 22 mm-connector (KOMBI)		Manufacturer
	Inner cannula with 15 mm-connector (UNI)		Use by
	Inner cannula with 15 mm-swivel connector (VARIO)		Consult instructions for use
	With suction line (SUCTION)		CE marking identification number of the notified body
	With oxygen support		Batch code
	With low-pressure cuff (CUFF)		Catalogue number
	Several small holes in a specific arrangement in the center of the tube, is referred to as sieve (LINGO)		Sterilized using ethylene oxid
	Speaking valve (PHON)		Do not resterilize
	Obturator (insertion aid)		Content (in pieces)
	Tube holder included		Single patient use
	Outer cannula fenestrated		Do not use if package is damaged
	Inner cannula fenestrated		Keep dry
	Length medium		Keep away from sunlight
	Outer cannula fenestrated		Decannulation plug
	For ventilation		Medical devices
	Length long		

# DURACUFF® TRACHEOSTOMY TUBES

EN

## I. FOREWORD

These instructions for use are valid for Fahl® tracheostomy tubes. The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling of Fahl® tracheostomy tubes.

**Please read the instructions for use carefully before using the product for the first time!**

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the tracheostomy tube. It contains important information on the product!

## II. INTENDED USE

Fahl® tracheostomy tubes are intended to stabilise the tracheostoma following laryngectomy or tracheostomy.

The tracheostomy tube is designed to keep the tracheostoma open.

Tracheostomy tubes with low pressure cuff are always indicated after tracheostomy has been performed for any cause if the cannula needs to be sealed against the wall of the trachea.

When used for the first time the products must only be selected, used and inserted by clinicians or medical professionals trained in their use.

**Users must be trained by medical professionals on how to use the Fahl® tracheostomy tubes safely and appropriately.**

Overlong tubes can be indicated particularly in cases of deep-seated tracheal stenoses.

Fahl® LINGO tube variants are intended for tracheotomised patients who retain their larynx or for laryngectomees who use a shunt valve (voice prosthesis users).

## III. WARNINGS

**In case of an unstable tracheostoma, the airways must always be secured before removing the tracheostomy tube and a replacement tube must be kept ready for insertion. The replacement tube must be inserted immediately, even before cleaning and disinfecting the changed tube.**

**Patients must be trained by medical professionals on how to use the Fahl® tracheostomy tubes safely.**

**Fahl® tracheostomy tubes must never be blocked, e.g. by mucus or encrustation. Blockage can cause asphyxiation!**

**Mucus in the trachea can be aspirated (removed by suction) through the tracheostomy tube with a tracheal suction catheter.**

**Damaged tracheostomy tubes must not be used and must be discarded immediately. Use of a damaged tube may result in airway compromise.**

**When the tracheostomy tubes are inserted or removed, this can lead to irritations, coughing or slight bleeding. If bleeding persists, consult your doctor immediately!**

**Tracheostomy tubes are single-patient products and intended for use by one patient only.**

**Do not re-use with another patient. Do not reprocess for re-use by another patient.**

**Do not use tracheostomy tubes during laser treatment (laser therapy) or with electrosurgical devices! Contact with a laser beam may damage the tube.**

### CAUTION!

**Tracheostomy tubes with speech function are only recommended for tracheotomised patients with normal secretion and normal mucous tissue.**

**Tracheostomy tubes that contain metal parts must under no circumstances be used during a radiation therapy (radiotherapy), as this can lead to serious skin damage! If you need to wear a tracheostomy tube during a radiation therapy, then only use plastic tracheostomy tubes without any metal parts. In the case of plastic tracheostomy tubes with silver speaking valve, the valve can for instance be removed completely from the tube, including the safety chain, by removing the inner cannula from the outer cannula before starting a radiation therapy.**

### CAUTION!

**If there is strong secretion, a tendency to form granulation tissue, if the patient is undergoing radiation therapy or if scabs have formed, a perforated cannula version is only recommended under regular surveillance by a doctor and provided that shorter replacement intervals are observed (as a rule once a week), because the perforations in the outer tube can increase the formation of granulation tissue.**

## IV. COMPLICATIONS

The following complications could emerge when using this device:

Contamination of the stoma may make it necessary to remove the tube. Contaminations can furthermore cause infections that necessitate the use of antibiotics.

Unintentional Inhalation of a tube that was not sized correctly may need to be removed by a doctor/physician. If mucus blocks the tube it should be removed and cleaned.

## V. CONTRAINDICATIONS

Do not use if the patient is allergic to the material.

### CAUTION!

Never under any circumstances use cannula variants without low-pressure cuff for patients on mechanical respiration!

### CAUTION!

Only use perforated/fenestrated cannula variants during artificial respiration after consultation with the treating doctor/physician.

### CAUTION!

For patients on ventilation, tube versions with outer and inner cannula must only be used in combination with DURACUFF® CLIP with 15 mm swivel connector for arresting the inner cannula. Otherwise, the patient's life is at risk in the event of disconnection (detachment of the inner cannula from the outer cannula)!

### CAUTION!

Tracheostomy tubes with speaking valves must not be used on any account by laryngectomised patients, as this could lead to serious complications, even suffocation!

## VI. PRECAUTIONS

The correct tracheostomy tube size should be selected by the treating physician or by trained medical professionals.

Use only auxiliary equipment with a 22 mm connection with Fahl® KOMBI adapters to prevent accidental loss of accessories or damage to the tube.

We highly recommend always keeping at least two replacement tubes on hand to ensure continuous supply.

### CAUTION!

When increased tensile forces occur during mechanical ventilation, e.g. due to stiff swivel connectors connected to the tracheostomy tube or due to uncontrolled movements of the patient, the inner cannula can be accidentally twisted out of the outer cannula. The patient must therefore be monitored, or the tracheostomy tube and/or the tube system must be replaced if necessary.

Do not perform any type of adjustment, repair or alteration on a tracheostomy tube or compatible accessories. If a product is damaged, it must be discarded immediately.

## VII. PRODUCT DESCRIPTION

The DURACUFF® tracheostomy tubes are products that are made of medical-grade polyvinyl chloride.

The tracheostomy tubes are made of thermosensitive medical synthetics that develop their optimal product properties at body temperature.

We supply the Fahl® tracheostomy tubes in various sizes and lengths.

The fenestrated inner cannulas have a red ring that indicates that they are fenestrated, thus enabling the user to discern whether the inner cannula is fenestrated even when the inner cannula is inserted.

Fahl® tracheostomy tubes are reusable medical products intended for single-patient use.

Fahl® tracheostomy tubes are only to be used by the same patient and not by another patient.

The package contains 1 tube, supplied sterile by ethylene oxide sterilization (EO).

Please see the size table in the appendix.

A particular feature of the Fahl® tracheostomy tubes is the specially shaped neck flange.

The tube tip is rounded to prevent irritation of the mucus membranes in the trachea.

There are two lateral fastening eyelets on the neck flange for attaching a tube holder.

It may be advisable to use tubes of different lengths alternately so the tip of the tube does not always touch the same spot in the trachea and cause irritation. This prevents pressure points or the formation of granulation tissue. Always discuss the exact procedure with your doctor/physician.

## MRT (MAGNETIC RESONANCE TOMOGRAPHY) INFORMATION

### CAUTION!

As the tracheostomy tubes with cuff have a small metal spring in the non-return valve of the pilot balloon with inflation tube, the tracheostomy tube with cuff should not be used while an MRT (magnetic resonance tomography) examination is being performed.

MRT is a diagnostic imaging technique in which inner organs, tissues and joints are depicted with the help of magnetic fields and radio waves. Metallic objects can be drawn into the magnetic field and can induce changes by their acceleration. Although the metal spring is extremely small and lightweight, interactions that could give rise to health impairments or to malfunction or damage of the technical equipment used or of the tracheostomy tube itself nevertheless cannot be ruled out in this context. If wearing a tracheostomy tube is indicated for keeping the tracheostoma open, we recommend use of a metal-free tracheostomy tube instead of the tracheostomy tube with cuff while an MRT is being performed.

#### **CAUTION!**

**Do not use a tracheostomy tube with low-pressure cuff while undergoing a nuclear magnetic resonance tomography/MRT!**

#### **1. Neck flange**

A particular feature of the Fah!® tracheostomy tube is the specially shaped neck flange, which is designed to fit the anatomy of the neck.

The size specifications are printed on the neck flange.

The neck flange of the tracheostomy tubes has two side eyelets for inserting a tube holder.

All Fah!® tracheostomy tubes with fastening eyelets include a disposable tube holder. The tube holder secures the tracheostomy tube to the neck.

The medium-length cannulas (see size table in the appendix) are colour-coded by means of a yellow ring on the inner cannula and the yellow imprint on the neck flange of the cannula as well as on the pilot balloon.

The insertion aid (obturator) included in delivery facilitates insertion of the tracheostomy tube.

Please follow the instructions provided with the tube holder carefully when attaching and removing it from your tracheostomy tube.

Ensure that the Fah!® tracheostomy tube is placed in the tracheostoma without tension and the position is not changed when fastening the tube holder.

#### **2. Connectors and adapters**

Connectors and adapters are used to connect compatible tube accessories.

The options available to a patient depend on the clinical condition, such as status post laryngectomy or tracheostomy.

Connectors/adapters are normally permanently attached to the inner cannula. The connector in this case is the universal (15 mm) connector with the aid of which so-called artificial noses or heat and moisture exchangers (HMEs, filters for exchanging heat and moisture) can be attached.

This connector is also available in a special version as 15 mm swivel connector. The swivel version of the 15 mm connector is suitable, for instance when using a respirator tube system, for absorbing the torsional forces which arise in such situations and to relieve strain on the tracheostomy tube and stabilise it in position so that irritations of the mucous membrane in the trachea are minimized.

The 22 mm combi-adapter can be used to fasten compatible filter and valve systems with 22 mm attachments – such as, but not limited to the HUMIDOPHONE® speaking valve with filter function (REF 46480), HUMIDOTWIN® Heat and Moisture Exchanger (HME) (REF 46460), COMBIPHON® speaking valve (REF 27131), LARYVOX® HME filter cassettes (REF 49800).

#### **3. Cannula tube**

The cannula tube borders directly on the neck flange and directs the airflow into the trachea.

The Easy Lock hexagonal lock serves to securely attach inner cannulas and compatible accessories.

The X-ray contrast strip in the tube that runs along the side allows the tube to be depicted radiographically and the position of the tube to be checked.

#### **3.1 Low-pressure cuff**

In the product versions with low-pressure cuff, the very thin-walled, high-volume low-pressure cuff adapts well to the trachea and ensures reliable sealing if inflated correctly. The low-pressure cuff can be inflated like a balloon. The small pilot balloon on the inflation tube indicates whether the tracheostomy tube is in sealed (inflated) or non-sealed condition.

The low pressure cuff is inflated via a pilot line with one-way valve and pilot balloon.

##### **3.1.1 Leakproofness test of the cannula and the low-pressure cuff (if present)**

The cannula and the low pressure cuff must be tested for leakages directly before and after every insertion and thereafter at regular intervals. For this purpose, inflate the low-pressure cuff to 15 to 22 mmHg (1 mmHg corresponds to 1.35951 cmH<sub>2</sub>O) and then observe whether a spontaneous drop in pressure occurs (recommended for inflating and testing:

Cuff pressure gauge control inflator, REF 19500). There should be no substantial drop in pressure during the observation period. This leakproofness test must also be performed prior to every renewed insertion (e.g. after cleaning the tracheostomy tube) (see picture 7c).

**The following are possible signs indicating leakages in the cuff (balloon):**

- Externally visible damage to the balloon (holes, fissures etc.)

- Audible and/or tactile perception of air escaping from the balloon
- Water in inlet tubes of the tracheostomy tube (after cleaning!)
- Water in the cuff (after cleaning!)
- Water in the pilot balloon (after cleaning!)
- No cough stimulus when pressure is applied to the pilot balloon

#### **CAUTION!**

When testing the balloon or when inserting, removing or cleaning the tracheostomy tube, never under any circumstances use sharp or pointed objects, such as for instance forceps or clamps, since these can damage or even destroy the balloon. Should one of the above-named signs of leakage be detectable, the tracheostomy tube must under no circumstances be used any more since it is no longer functioning properly!

#### **3.2 Obturator**

**Please check whether the obturator can easily be removed from the cannula before inserting the tracheostomy tube!**

**After checking the obturator for ease of movement, reinsert the obturator into the cannula for inserting the tracheostomy tube.**

#### **3.3 Suction opening (only in Suction tracheostomy tube versions)**

Secretions that have accumulated above the inflated cuff can be removed via the suction opening in the outer tube of the Suction tracheostomy tubes.

The suction tube leading to the outside can be connected to a syringe or a suction device. The user must decide in consultation with the treating physician and after appropriate risk analysis, which of these suction variants to use. The individual disease condition of the patient should be considered here. Suction using a suction device may in any case only be performed if the suction device is equipped with vacuum adjustment. The suction pressure must never exceed  $-0.2$  bar.

When using the suction variants of the tracheostomy tubes, the secretions that have accumulated above the cuff must be suctioned off via the suction opening on the tracheostomy tube directly before deflating the cuff in order to prevent aspiration of secretions. In this way, additional suctioning with a suction catheter while deflating the cuff may be rendered unnecessary, which facilitates handling significantly.

#### **CAUTION!**

**When using any of the tracheostomy tube versions with a suctioning device (SUCTION), take special care during the suctioning procedure to generate a vacuum for as short a time as possible; drying out of the subglottic space as a result of this must be avoided.**

#### **CAUTION!**

**Contraindicated in patients with increased bleeding tendency (e.g. during anticoagulant therapy). In such cases, the Suction tracheostomy tube with suction opening may not be used, since suctioning entails an increased risk.**

#### **4. Inner cannula**

The convenient Easy Lock hexagonal lock serves to securely attach the inner cannula and compatible accessories.

The inner cannulas are, depending on specification, either equipped with a standard attachment or permanently connected to particular adapters/connectors, or they can be connected to detachable accessories such as for instance speaking valves.

The inner cannulas can easily be removed from the outer cannula, thus allowing the air supply to be increased quickly if necessary (for instance in case of dyspnoea).

The connection between inner and outer cannula must first be released by slight lifting (unbuttoning with thumb and index finger). This causes holding prongs/engaging elements of the inner cannula to be pulled off the Easy Lock hexagonal connection of the outer cannula.

Inner cannulas must never be used without outer cannula but must always be affixed to the outer cannula.

#### **4.1 Speaking valves**

Tracheostomy tubes with speaking valve (LINGO/PHON) are used after tracheotomy with a complete or partially retained larynx and enable the user to speak.

In the case of tracheostomy tubes with silicone speaking valve, the valve is directly attached to the inner cannula and is included in delivery of these tube versions.

The HUMIDOPHON® speaking valve or the COMBIPHON® speaking valve are also available. These can be used by tracheotomised patients who for instance use a perforated tracheostomy tube (with inner cannula) with 22 mm combi-adapter.

In tracheostomy tubes with silver speaking valve, the speaking valve can be detached from the inner cannula by pushing it out.

In tracheostomy tubes with silicone speaking valve, the speaking valve can be detached from the inner cannula by pulling it off.

## 5. Decannulation plug

The decannulation plug is included in delivery of the tracheostomy tubes with speaking valve and may only be used by tracheotomised patients who have retained their larynx. It may only be inserted under supervision by a doctor/physician. It allows the tracheostomy tube to be sealed and the air supply to be interrupted for a short time and helps the patient train to regain control of breathing via mouth/nose.

### CAUTION!

**Contraindicated in laryngectomised patients and patients with chronic obstructive pulmonary disease (COPD)! In such cases, the decannulation plug must not under any circumstances be inserted!**

**The decannulation plug must never under any circumstances be used when the tracheostomy tube is in sealed condition! The decannulation plug may only be used with perforated outer cannula without inner cannula.**

### CAUTION!

**A decannulation plug is enclosed with the tracheostomy tubes with speaking valve for the purpose of preparing for possible decannulation following temporary tracheotomy. The plug can be used to block the air supply via the tracheostomy tube for a short time so that the patient can get used to breathing via mouth/nose again. Decannulation may only be carried out under medical supervision. The plug may only be inserted upon instruction by the doctor/physician. There is a risk of suffocation! Please also make sure you observe the description of indications provided with the respective product variants/specifications!**

## 6. DURACUFF® CLIP

Use of the DURACUFF® tracheostomy tube version with inner cannula in patients undergoing mechanical ventilation was hitherto out of the question due to the risk of disconnection. The DURACUFF® CLIP (REF 13299) for locking the inner cannula enables fixation of the inner cannula with 15 mm swivel connector, so that DURACUFF® tracheostomy tubes with these inner cannulas can also be used during mechanical ventilation of a patient.

The DURACUFF® CLIP is a single-patient product.

It is made of medical-grade synthetic material.

The DURACUFF® CLIP is an adapter for locking DURACUFF® inner cannulas with 15 mm swivel connector in place in DURACUFF® outer cannulas.

Use of the following DURACUFF® tracheostomy tubes with inner cannula with 15 mm swivel connector is thus allowed in mechanically ventilated patients, provided that the tubes are used in combination with the clip:

- DURACUFF® UNI tracheostomy tubes with low-pressure cuff, during mechanical ventilation only in combination with the inner cannulas with 15 mm swivel connector (versions REF 13002, 13102, 13202)
- DURACUFF® SUCTION tracheostomy tube with low-pressure cuff, during mechanical ventilation only in combination with the inner cannulas with 15 mm swivel connector (versions REF 13032, 13132, 13232)
- DURACUFF® LINGO tracheostomy tube with low-pressure cuff, during mechanical ventilation only in combination with the closed inner cannulas with 15 mm swivel connector (versions REF 13022, 13122, 13222)
- DURACUFF® PHON tracheostomy tube with low-pressure cuff, during mechanical ventilation only in combination with the closed inner cannulas with 15 mm swivel connector (versions REF 13012, 13112, 13212)
- DURACUFF® KOMBI tracheostomy tube with low-pressure cuff, during mechanical ventilation only in combination with the closed inner cannulas with 15 mm swivel connector (versions REF 13072, 13172, 13272)
- DURACUFF® KOMBI VOICE tracheostomy tube with low-pressure cuff, during mechanical ventilation only in combination with the closed inner cannulas with 15 mm swivel connector (versions REF 13082, 13182, 13282)
- DURACUFF® VARIO-KOMBI LINGO-PHON SUCTION tracheostomy tube with low-pressure cuff, during mechanical ventilation only in combination with the closed inner annulas with 15 mm swivel connector (versions REF 13042, 13142, 13242)

The DURACUFF® CLIP must always be used for locking the inner cannula in position when DURACUFF® tracheostomy tubes with 15 mm swivel connector are used for patients who are being mechanically ventilated! The DURACUFF® CLIP is attached to the DURACUFF® tracheostomy tube before starting mechanical ventilation. It secures the inner cannula via the lock with the outer cannula and so prevents the inner cannula from accidentally detaching.

### CAUTION!

**Make sure - especially in mechanically ventilated patients - that the Easy Lock hexagonal connection is closed firmly and that the DURACUFF® CLIP is securely affixed to the tracheostomy tube.**

**CAUTION!**

If the DURACUFF® CLIP is not used, only the following tracheostomy tube versions (without inner cannula) may be used for mechanically ventilated patients:

- DURACUFF® VARIO tracheostomy tube without inner cannula (versions REF 13000, 13100, 13200)
- DURACUFF® VARIO SUCTION with universal attachment without inner cannula (versions REF 13040, 13140, 13240)

**CAUTION!**

Never use the DURACUFF® CLIP for locking the inner cannula in place if the functionality of the clip is impaired or if the clip shows damages such as sharp edges or cracks, since this could lead to injuries or disconnection. Defective clips must not be used and must be disposed of. Danger can be avoided by regular checking!

**6.1 Preparing for affixing the DURACUFF® CLIP on the tracheostomy tube before the patient is connected to a mechanical ventilation system:**

The DURACUFF® inner cannula with 15 mm swivel connector can be quickly and easily inserted into the DURACUFF® outer cannula; however, the inner cannula is only firmly attached to the outer cannula when the holding prongs of the inner cannula firmly engage the Easy Lock hexagonal connection of the outer cannula.

For inserting the tracheostomy tube, please observe the directions in these instructions for use.

Please note that an obturator can no longer be used if the DURACUFF® CLIP was already affixed before the tracheostomy tube was inserted and the inner cannula has thus already been inserted into the outer cannula. For complete change of tracheostomy tube, we recommend to slide the clip onto the neck flange before inserting the tracheostomy tube. Before the tracheostomy tube is used on the patient, the clip should be subjected to a function check and the handling procedure should at the same time be practiced several times. Please also observe the following points when attaching the clip:

When attaching the clip to the neck flange from above, avoid pinching the supply tubes on the tracheostomy tube (suctioning device and supply tube for inflating the cuff).

- Make sure that the inner cannula is correctly fixed in place and that the clip firmly encompasses the tracheostomy tube.

The DURACUFF® CLIP can also be attached or removed with the tracheostomy tube lying in place - In this way, it is possible to change only the inner cannula without having to remove the outer cannula. If you apply the clip when the tracheostomy tube is already inserted, please observe the following in addition to the points named above:

- Attach the clip carefully to avoid damaging the skin or the mucous membrane!
- Attach the clip before the compress which is already in place!

**6.2 Affixing and removing the DURACUFF® CLIP**

The DURACUFF® CLIP adapter for locking the inner cannula with 15 mm swivel connector in place is attached to the outer cannula. The clip is a semicircular attachment that is slid onto the neck flange of the tracheostomy tube from above (see picture 7d). This is done when the inner cannula is already in the outer cannula. When the clip is seated correctly, the flat, semicircular rear wall of the clip is located behind the neck flange (see picture 7e) and the front part of the clip firmly encompasses the ring of the inner cannula. In this way, the inner cannula is correctly fixed in position (see picture 7f).

To remove the inner cannula, detach the DURACUFF® CLIP from the connector (15 mm swivel connector) by applying upward pressure at both clip ends. The clip is under tension to achieve firm fixation of the inner cannula. Therefore, hold you index finger against the edge of the clip from above when removing the clip in order to prevent uncontrolled movements.

Cleaning and disinfection of the DURACUFF® CLIP can be performed as described in Chapter IX on cleaning and disinfection.

The clip should be replaced after not more than 29 days (reckoned from opening of the sterile package). As regards period of use and replacement, please also heed the corresponding chapter, Chapter XI on period of use.

**CAUTION!**

**Make sure to read all accompanying product information, instructions for use, indications and contraindications. Discuss use of the product with your doctor/physician before first use.**

**VIII. DIRECTIONS FOR TUBE INSERTION AND REMOVAL****For the Doctor/Physician**

The appropriate tracheostomy tube must be selected by a doctor/physician or trained medical professionals.

Select a tube that fits the patient's anatomy to optimise comfort and ventilation (breathing in and out).

The inner cannula can be removed at any time to increase air supply or for cleaning. This can for instance be necessary if the cannula is clogged up with secretion residues which cannot be removed by coughing or because no equipment for suctioning off the secretions is available.

## For the Patient

### CAUTION!

Only insert the tracheostomy tube with the cuff in fully deflated condition (see picture 7a)!

### CAUTION!

Carefully examine the sterile packaging to ensure that it has not been tampered with or damaged. Do not use the product if the packaging has been damaged.

Check the use-by or expiry date. Do not use after this date.

It is advisable to use sterile disposable gloves.

Carefully examine the tube before first use to make sure that it is not damaged and that there are no loose parts.

Should you notice any anomaly or anything unusual, DO NOT use the tube. Return the tube to the manufacturer for inspection.

If secretion collects in the lumen of the Fahl® tracheostomy tube or stoma button and cannot be removed by coughing or aspiration, the tube should be removed and cleaned.

After cleaning and/or disinfection, carefully examine the Fahl® tracheostomy tube for sharp edges, cracks, or other signs of damage, since these may impair function and/or injure the mucus membranes in the airways.

Never under any circumstances continue using damaged tracheostomy tubes.

### WARNING

Patients must be briefed by trained medical professionals in the safe use of the Fahl® tracheostomy tubes.

#### 1. Insertion of the tube

##### Step-by-step instructions to insert Fahl® tracheostomy tubes.

Before application, users should clean their hands (see picture 3).

Remove tube from the package (see picture 4).

If an obturator is to be used, this must first be fully inserted into the cannula tube so that the collar on the gripping piece of the obturator comes into contact with the outer edge of the 15 mm connector and the tip of the olive projects beyond the tip of the cannula (proximal end of cannula). The obturator must be held in this position during the entire procedure.

##### When using tracheostomy tubes with low-pressure cuff, pay special attention to the following points:

Before inserting the tracheostomy tube, check the cuff (balloon) as well - it must be free of damages of any kind and must be leakproof so as to ensure tight sealing as required. We therefore recommend to perform a leakproofness test prior to every insertion (see Section VII, No. 3.1.1). The balloon must be emptied completely prior to insertion of the tracheostomy tube (see picture 7b)! When using an aid for dilating the tracheostoma, take care to ensure that the tracheostomy tube, and especially the cuff, is not damaged by friction.

Next, push a tracheal compress, e.g. SENSOTRACH® DUO (REF 30608) or SENSOTRACH® 3-Plus (REF 30780), onto the tracheostomy tube.

To facilitate insertion of the tracheostomy tube, it is advisable to lubricate the outer tube by wiping it with an OPTIFLUID® Stoma oil wipe (REF 31550) which allows the stoma oil to be applied evenly to the whole surface of the tube (see picture 4a and 4b).

If you are inserting the tube yourself, use a mirror to make insertion of the Fahl® tracheostomy tube easier.

When inserting the Fahl® tracheostomy tube, hold it by the neck flange with one hand (see picture 5).

Pull the tracheostoma slightly apart with your free hand to allow the tip of the tube to fit into the tracheostoma more easily.

Special aids for dilating the tracheostoma are also available (tracheal dilator, REF 35500). These allow the tracheostoma to be dilated gently and evenly, for instance also in emergency situations with collapsing tracheostoma (see picture 6).

Make sure that the tube is not damaged by friction when using an instrument for assistance.

Now carefully insert the tube into the tracheostoma during the inspiration phase (while breathing in) while tilting your head slightly back (see picture 7).

Advance the tube into the trachea.

Straighten your head once the tube has been inserted further into the trachea.

If an obturator is used, this must then immediately be removed from the tracheostomy tube.

The tracheostomy tubes should always be attached with a special tube holder. This stabilises the tube and thus ensures that the tracheostomy tube is securely seated in the tracheostoma (see picture 1).

### 1.1 Inflating the low-pressure cuff (if present)

The low-pressure cuff is inflated by applying a defined pressure to the cuff via the Luer connection (standardised conical connection) of the inflation tube by means of a cuff pressure gauge (e.g. MUCOPROTECT®, REF 19500). If not instructed otherwise by the doctor/physician, we recommend a cuff pressure of at least 15 mmHg (20 cmH<sub>2</sub>O) to 22 mmHg (30 cmH<sub>2</sub>O). The cuff pressure should never under any circumstances exceed 22 mmHg (approx. 30 cmH<sub>2</sub>O).

Inflate the low pressure cuff at most to this target pressure and check to make sure that sufficient air is supplied via the tracheostomy tube.

Always make sure that the low-pressure cuff is undamaged and in perfect working order.

If the desired sealing is not achieved even after trying repeatedly with the specified limit volume, a tracheostomy tube with larger diameter may be indicated.

**The correct cuff pressure must be checked regularly, i.e. at least every 2 hours.**

#### CAUTION!

All instruments used for inflating the cuff must be clean and free of foreign particles! Detach the instruments from the Luer connection of the inflation tube as soon as the cuff has been inflated and close the connection with the cap.

#### CAUTION!

If the maximum pressure is exceeded for longer periods of time, the blood circulation in the mucus membrane can be impaired (risk of ischaemic necrosis, pressure ulcers, tracheomalacia, tracheal stenosis, pneumothorax). In patients undergoing artificial respiration, the cuff pressure should not be allowed to drop below the cuff pressure value specified by the doctor/physician in order to prevent unnoticed aspiration. Hissing noises in the region of the balloon, especially during expiration, indicate that the trachea is insufficiently sealed by the balloon. If the trachea cannot be sealed with the pressure values specified by the doctor/physician, the entire air should be withdrawn again from the balloon and the sealing process should be repeated. If this does not lead to success, we recommend to use the next larger tracheostomy tube with balloon. Due to the permeability of the balloon wall for gases, it is normal for the pressure in the balloon to decline slightly over time, but it can on the other hand also rise unintentionally during gas anaesthesia. Regular pressure monitoring is therefore urgently recommended.

The cuff must never under any circumstances be inflated with excessive amounts of air, since this can lead to damage of the tracheal wall, tears in the low pressure cuff with subsequent deflation, or distortion of the cuff, in which case airway obstruction cannot be ruled out.

#### CAUTION!

During anaesthesia, the cuff pressure can rise/fall due to nitrous oxide (laughing gas).

### 2. Removing the tube

#### CAUTION!

Accessories such as a tracheostoma valve or HME (Heat Moisture Exchanger) must be removed first before proceeding to remove the Fahl® tracheostomy tube.

#### CAUTION!

If the tracheostoma is unstable, or in emergency situations (puncture/dilation tracheostomy), the tracheostoma can collapse after withdrawal of the tracheostomy tube, thereby impairing air supply. A fresh tracheostomy tube must be kept ready for use in such cases and must be quickly inserted if necessary. A tracheal dilator (REF 35500) can be used for temporarily securing the air supply.

The cuff must be emptied before removing the tracheostomy tube. The head should be tilted back slightly for removal of the tube.

#### CAUTION!

Never use a cuff pressure gauge to empty the low-pressure cuff. Always use a syringe for this.

Before the air is removed from the balloon by means of a syringe and the tracheostomy tube is withdrawn, the region of the trachea above the balloon must first be cleaned by suctioning off secretions and mucus. If the patient is responsive and reflexes are intact, it is recommended that the patient be suctioned while at the same time unblocking the tracheostomy tube. Suctioning is performed by inserting a suction catheter through the cannula tube into the trachea. In this way, suctioning can be performed without any problems and gently for the patient and cough stimulus and the risk of aspiration are minimised.

**Next, deflate the low pressure cuff while suctioning off at same time.**

If secretions are present, these are now taken up by the suction tube and can no longer be aspirated. Please note that the tracheostomy tube must in every case be cleaned, if necessary disinfected, and lubricated with stoma oil as specified below prior to reinsertion.

Proceed very carefully to avoid injury to the mucus membranes.

The tube must always be cleaned and, if necessary, disinfected as follows before re-inserting according to the instructions provided below.

### Step-by-step instructions to remove the Fah!® tracheostomy tubes:

The tracheostomy tubes should be removed with the head slightly tilted back. Grip the tube or button at the side by the neck flange or the housing (see picture 7).

Carefully remove the tracheostomy tubes.

The outer cannula (if an inflated low-pressure cuff is present) remains in the tracheostoma.

For this purpose, the connection between inner and outer cannula must be released by slight lifting (unbuttoning with thumb and index finger). This causes holding prongs/engaging elements of the inner cannula to be pulled off the Easy Lock hexagonal connection of the outer cannula.

For reinsertion, proceed in reverse sequence as described above to reinsert the inner cannula and attach it to the Easy Lock hexagonal connection of the outer cannula.

## IX. CLEANING AND DISINFECTION

### CAUTION!

**For reasons of hygiene and to avoid the risk of infection, you should thoroughly clean the Fah!® tracheostomy tube at least twice a day, more often during heavier secretion production.**

**In case of an unstable tracheostoma, the airways must always be secured before removing the tracheostomy tube and a replacement tube must be kept ready for insertion. The replacement tube must be inserted immediately, even before cleaning and disinfecting the changed tube.**

### CAUTION!

**For reasons of hygiene and to avoid the risk of infection, you should thoroughly clean the Fah!® tracheostomy tube at least twice a day, more often during heavier secretion production.**

Keep this in mind, your personal cleaning schedule, including additional disinfection if required, must always be determined in consultation with your doctor/physician and according to your individual needs.

**Regular disinfection is only required if this is indicated from the medical point of view as ordered by the doctor/physician. The reason for this is that the upper airways are not free of microbes in healthy patients either.**

General cleaning may be insufficient to meet the special hygiene requirements to prevent infections for patients with specific problems (e.g. MRSA, ORSA etc.), who are subject to an increased danger of re-infection. We recommend chemical disinfection of the tubes according to the instructions provided below. Please consult your doctor/physician.

### CAUTION!

**Cleaning agent and disinfectant residues on the tracheostomy tube can lead to irritations of the mucus membranes or other health impairments.**

Tracheostomy tubes are considered as belonging to the category of instruments with hollow cavities. Special care must therefore be taken when performing cleaning or disinfection that the tube is fully wetted with the solution used and is unobstructed (see picture 8).

The solutions used for cleaning and disinfection must be freshly prepared prior to each use.

### 1. Cleaning

Fah!® tracheostomy tubes must be cleaned / replaced regularly in accordance with the individual patient's needs.

Only use cleaning agents if the tube is outside the tracheostoma.

A mild, pH-neutral washing lotion can be used for cleaning the tube. We recommend using the special cannula-cleaning powder (REF 31110) as directed by the manufacturer.

Never clean Fah!® tracheostomy tubes with cleaning agents that are not approved by the manufacturer. Never use aggressive household cleaners, detergents, high-concentration alcohol or denture cleaners.

This may be very hazardous to health. The tube may also be destroyed or damaged.

Tubes can also be cleaned by thermal disinfection at not more than 65 °C. Use clean water with a temperature of not more than 65 °C for this purpose. Make sure that the temperature is kept constant (use a thermometer to monitor the temperature) and never under any circumstances expose the tubes to boiling water for disinfection. This could seriously damage the tracheostomy tube.

### Cleaning Steps

Remove any additional items before cleaning.

First, rinse the tube thoroughly under running water (see picture 9).

Use only lukewarm water to prepare the cleaning solution and observe the instructions for use of the cleaning agent.

We recommend the use of the cannula-cleaning tub with sieve insert (REF 31200) to facilitate the cleaning procedure.

To do this, hold the sieve insert by the upper edge to prevent contact with and contamination of the cleaning solution (see picture 10).

Place only one tracheostomy tube at a time in the sieve insert of the cannula-cleaning tub. If more than one tube is cleaned at a time, they may be compressed too strongly and damaged by excessive pressure.

You can place the inner and outer cannula side by side for this procedure.

Immerse the sieve insert with the tracheostomy tube components in the prepared cleaning solution.

After the immersion time has expired (see instructions for use of the cannula cleaning powder), wash the tracheostomy tube thoroughly several times with lukewarm, clear water (see picture 9). There must be no residues of the cleaning agent on the tube when it is inserted into the tracheostoma.

If necessary, for example, if persistent secretion residues cannot be removed by the cleaning bath, additional cleaning with a special cannula-cleaning brush (OPTIBRUSH® REF 31850 or OPTIBRUSH® PLUS with fibre top, REF 31855) may help. Only use the cleaning brush, if the tube is removed and already outside the tracheostoma.

Always insert the cannula cleaning brush into the cannula from the cannula tip (see picture 11).

Use the brush as directed and proceed with great care to avoid damage to the soft tube material.

The outer cannula must only be cleaned with the cuff inflated, and the safety balloon must not come into contact with the cleaning solution in order to prevent intrusion of cleaning solution into the balloon, since this would lead to serious functional impairment and health risks for the user.

In tracheostomy tubes with speaking valve, the speaking valve must first be detached from the inner cannula. The valve itself must not be cleaned with the brush since it could otherwise be damaged or break off.

Thoroughly rinse the tube or button under lukewarm running water or sterile saline solution (0.9% NaCl solution).

Dry the tube thoroughly with a clean and lint-free cloth after wet cleaning.

Never use tracheostomy tubes with impaired functionality or with damages such as sharp edges or cracks, because this may lead to injuries of the mucus membranes in the trachea. If damages are found, the tube must not be used under any circumstances.

The obturator can be cleaned in the same way as the tracheostomy tube.

### **Cleaning the tracheostomy tubes with low-pressure cuff**

The indwelling time and thus the replacement intervals of the tracheostomy tube must be determined individually in consultation with the treating doctor/physician. The outer cannula must be cleaned or the tracheostomy tube replaced after 1 week at the latest though, since after that the risk of formation of granulation tissue or of tracheomalacia, etc. can increase significantly, depending on the patient's disease condition.

The outer cannula with cuff can be cleaned/rinsed with sterile saline solution.

A brush must not be used for tracheostomy tubes with cuff in order to avoid damaging the balloon!

Damage to the balloon and hence a defect of the tracheostomy tube can only be avoided by careful and gentle handling.

## **2. Chemical Disinfection Instructions**

### **2.1 Disinfecting the inner cannula / cleaning the outer cannula without low-pressure cuff**

The Fah® tracheostomy tube can be disinfected by cold disinfection with special chemical disinfectants.

Disinfection should always be done if ordered by a doctor/physician due to specific health concerns caused by disease, infection, or your specific situation.

Disinfecting is generally indicated to prevent cross-infection and in inpatient situations (e.g. hospitals, nursing homes, and/or other health care facilities) to limit infections.

#### **CAUTION!**

**Always clean according to the cleaning procedure described above before proceeding to disinfection (if applicable).**

**Disinfectants that release chlorine or that contain strong alkalis or phenol derivatives must under no circumstances be used. This could badly damage or even destroy the tracheostomy tube.**

### **2.2 Disinfecting the outer cannula with low-pressure cuff**

Tracheostomy tubes with low-pressure cuff may only be disinfected provided that this is performed and checked with utmost care. The balloon must always be inflated beforehand.

#### **Disinfection steps**

For this purpose, OPTICIT® tube disinfectant (REF 31180; not available in the USA) should be used in accordance with the manufacturer's instructions.

As alternative, we recommend a disinfectant based on glutaric aldehyde as active ingredient (also available in the USA). Always observe the area of application and spectrum of activity specified for the disinfectant by the manufacturer.

Follow the instructions for the disinfectant.

Dry the tube thoroughly with a clean and lint-free cloth after wet cleaning.

### 3. Sterilisation/Autoclaving

Resterilisation is not allowed.

#### CAUTION!

**Heating to over 65°C, disinfection with boiling water or steam sterilisation are not allowed and lead to damage of the tracheostomy tube.**

### X. STORAGE/CARE

Currently not used and cleaned tubes should be stored in a clean plastic container, dry environment away from dust, sunlight and heat.

Still sterile packed replacement tubes should be stored in a dry environment away from dust, sunlight and heat.

After cleaning tracheostomy tubes, make sure that the balloon has been deflated for storage.

**After the inner cannula has been cleaned, disinfected if necessary and then dried, the outer surface of the inner cannula should be lubricated by applying stoma oil, e.g. using an OPTIFLUID® Stoma oil wipe.**

**Only use stoma oil (REF 31525) as lubricant, or OPTIFLUID® Stoma oil wipes (REF 31550).**

We highly recommend keeping at least two replacement tubes on hand to ensure continuous supply.

### XI. SERVICE LIFE

These tracheostomy tubes are sterile products for single-patient use.

The maximum period of use should not exceed 29 days.

The service life of a tube is influenced by many factors. The composition of the secretions, the quality of cleaning and other aspects, for example, are very important.

Damaged tubes must be replaced immediately.

#### CAUTION!

**Modifications to the tracheostomy tube of any kind, in particular perforations as well as repairs to the tube, may only be carried out by the manufacturer himself or by companies expressly authorised for this in writing by the manufacturer! Inexpertly performed work on tracheostomy tubes can lead to serious injuries.**

### XII. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the tracheostomy tube, especially due to shortening of the same or perforations, if these modifications or repairs were not carried out by the manufacturer himself. This applies both to damages to the tubes themselves caused thereby and to any consequential damages caused thereby.

If the tracheostomy tube is used beyond the period of use specified under Section XI, and/or if the tracheostomy tube is used, maintained (cleaned, disinfected) or stored in non-compliance with the instructions and specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious adverse event occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, this is to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

Product specifications described herein are subject to change without notice.

DURACUFF® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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All measurements listed in the tables below are for standard cannulas.  
Customised cannulas are available upon request.

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#### Sizing information:

All lengths in mm | GR = Size | AD1 = Outer diameter at the neck flange | AD2 = Outer diameter at cannula tip | ID1 = Inner diameter at the cannula tip of the outer cannula | ID2 = Inner diameter at the cannula tip of the inner cannula | L1 = Length over the outer curve | L2 = Length over the middle curve | ø-CD = CUFF diameter | ⊖-BW = Bending angle

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