












SPIRAFLEX®

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
PICTROGRAM LEGEND

Pictogrammes listed below you'll find on the product packaging, if applicable.

	No inner cannula		Date of manufacture
	With two inner cannulas		Use by
	Outer cannula including 15 mm-connector (UNI)		Consult instructions for use
	Inner cannula with low profile		Batch code
	Inner cannula flexible		Catalogue number
	Cannula length adjustable		Sterilized using ethylene oxid
	With suction line (SUCTION)		Do not resterilize
	Cuff		Content (in pieces)
	Sieving		Single patient use
	Speaking valve (PHON)		Do not use if package is damaged
	15 mm connector (UNI)		Keep dry
	Decannulation plug		Keep away from sunlight
	Obturator (insertion aid)		Medical devices
	Tube holder included		
	Length short		
	Length XL		
	For ventilation		
	MRT suitable		
	Manufacturer		

SPIRAFLEX® TRACHEOSTOMY TUBES

EN

I. FOREWORD

These instructions for use are valid for Fahl® tracheostomy tubes. The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling of Fahl® tracheostomy tubes.

Please read the instructions for use carefully before using the product for the first time!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the tracheostomy tube. It contains important information on the product!

II. INTENDED USE

Fahl® tracheostomy tubes are intended to stabilise the tracheostoma following laryngectomy or tracheostomy.

The tracheostomy tube is designed to keep the tracheostoma open.

Tracheostomy tubes with cuff are always indicated after tracheostomy has been performed for any cause if the cannula needs to be sealed against the wall of the trachea.

Tracheostomy tubes with adjustable neck flange enable individual patient-specific adaptation of the tube, particularly in obese patients and in situations with post-operative oedemas or collapsed tracheostoma.

When used for the first time the products must only be selected, used and inserted by clinicians or medical professionals trained in their use.

Overlong tubes can be indicated particularly in cases of deep-seated tracheal stenoses.

Fahl® LINGO tube variants are intended for tracheotomised patients who retain their larynx or for laryngectomees who use a shunt valve (voice prosthesis users).

III. WARNINGS

Patients must be trained by medical professionals on how to use the Fahl® tracheostomy tubes safely.

Fahl® tracheostomy tubes must never be blocked, e.g. by mucus or encrustation. Blockage can cause asphyxiation!

Mucus in the trachea can be aspirated (removed by suction) through the tracheostomy tube with a tracheal suction catheter.

Damaged tracheostomy tubes must not be used and must be discarded immediately. Use of a damaged tube may result in airway compromise.

When the tracheostomy tubes are inserted or removed, this can lead to irritations, coughing or slight bleeding. If bleeding persists, consult your doctor immediately!

Do not use tracheostomy tubes during laser treatment (laser therapy) or with electrosurgical devices! Contact with a laser beam may damage the tube.

CAUTION!

Tracheostomy tubes with speech function are only recommended for tracheotomised patients with normal secretion and normal mucous tissue.

CAUTION!

Tracheostomy tubes that contain metal parts (SPIRAFLEX®/SPIRAFLEX® XL) must under no circumstances be used during a radiation therapy (radiotherapy), as this can lead to serious skin damage! If you need to wear a tracheostomy tube during a radiation therapy, then only use plastic tracheostomy tubes without any metal parts (SPIRAFLEX® MRT). In the case of plastic tracheostomy tubes with silver speaking valve, the valve can for instance be removed completely from the tube, including the safety chain, by removing the inner cannula from the outer cannula before starting a radiation therapy.

CAUTION!

If there is strong secretion, a tendency to form granulation tissue, if the patient is undergoing radiation therapy or if scabs have formed, a perforated cannula version is only recommended under regular surveillance by a doctor and provided that shorter replacement intervals are observed (as a rule once a week), because the perforations in the outer tube can increase the formation of granulation tissue.

IV. COMPLICATIONS

The following complications could emerge when using this device:

Contamination of the stoma may make it necessary to remove the tube. Contaminations can furthermore cause infections that necessitate the use of antibiotics.

Unintentional Inhalation of a tube that was not sized correctly may need to be removed by a doctor/physician. If mucus blocks the tube it should be removed and cleaned.

V. CONTRAINDICATIONS

Do not use if the patient is allergic to the material.

CAUTION!

Never under any circumstances use cannula variants without cuff for patients on mechanical ventilation!

CAUTION!

Only use perforated/fenestrated cannula variants during artificial respiration after consultation with the treating doctor/physician.

CAUTION!

Tracheostomy tubes with speaking valves must not be used on any account by laryngectomised patients, as this could lead to serious complications, even suffocation!

VI. PRECAUTIONS

The correct tracheostomy tube size should be selected by the treating physician or by trained medical professionals.

In order to rule out accidental detachment of auxiliary equipment or damage to the tube, insert only auxiliary equipment with a 15 mm connection into Fahl® UNI adapters.

We highly recommend always keeping at least two replacement tubes on hand to ensure continuous supply.

CAUTION!

When increased tensile forces occur during mechanical ventilation, e.g. due to stiff swivel connectors connected to the tracheostomy tube or due to uncontrolled movements of the patient, the inner cannula can be accidentally twisted out of the outer cannula. The patient must therefore be monitored, or the tracheostomy tube and/or the tube system must be replaced if necessary.

CAUTION!

Do not perform any type of adjustment, repair or alteration on a tracheostomy tube or compatible accessories. If a product is damaged, it must be discarded immediately.

VII. PRODUCT DESCRIPTION

The SPIRAFLEX®/SPIRAFLEX® XL/SPIRAFLEX® SHORT tracheostomy tubes are products that are made of PVC (cannula tube reinforced with a metal spiral).

The SPIRAFLEX® MRT tracheostomy tubes are products that are made of PVC (cannula tube reinforced with a plastic spiral).

The tracheostomy tubes are made of thermosensitive medical synthetics that develop their optimal product properties at body temperature.

We supply the Fahl® tracheostomy tubes in various sizes and lengths.

Please see the size table in the appendix.

Fahl® tracheostomy tubes are reusable medical products intended for single-patient use.

Fahl® tracheostomy tubes are only to be used by the same patient and not by another patient.

The package contains 1 tube, supplied sterile by ethylene oxide sterilization (EO).

The dimensional stability of the tube is ensured by the spiral in the wall of the soft tube material.

SPIRAFLEX®/SPIRAFLEX® XL/SPIRAFLEX® MRT:

The variable length adjustment of the external section of the tracheostomy tube by means of the height-adjustable neck flange with screw attachment as well as the anatomically optimised bend angle of the highly flexible, spiral-reinforced cannula tube facilitate adaptation to the shape of the trachea and guarantee high wearing comfort as well as tension-free seating in the tracheostoma.

The position of the neck flange can be easily documented and reproduced by means of an imprinted scale with 10 mm graduations extending over a length of 5 cm.

It may be advisable to use tubes of different lengths alternately so the tip of the tube does not always touch the same spot in the trachea and cause irritation. This prevents pressure points or the formation of granulation tissue. Always discuss the exact procedure with your doctor/physician.

MRT (MAGNETIC RESONANCE TOMOGRAPHY) INFORMATION

CAUTION!

Since the tracheostomy tubes with cuff have a small metal spring in the non-return valve of the pilot balloon with inflation tube, the tracheostomy tube with cuff must not be used while an MRI (magnetic resonance imaging) examination is being performed.

MRT is a diagnostic imaging technique in which inner organs, tissues and joints are depicted with the help of magnetic fields and radio waves. Metallic objects can be drawn into the magnetic field and can induce changes by their acceleration. Although the metal spring is extremely small and lightweight, interactions that could give rise to health impairments or to malfunction or damage of the technical equipment used or of the tracheostomy tube itself nevertheless cannot be ruled out in this context. If wearing a tracheostomy tube is indicated for keeping the tracheostoma open, we recommend use of a metal-free tracheostomy tube instead of the tracheostomy tube with cuff while an MRT is being performed.

CAUTION!

The SPIRAFLEX® MRT tracheostomy tubes can even be worn by the patient during magnetic resonance tomography (MRT).

1. Neck flange

A particular feature of the Fahl® tracheostomy tube is the specially shaped neck flange, which is designed to fit the anatomy of the neck.

The size specifications are printed on the neck flange.

The neck flange of the tracheostomy tubes has two side eyelets for inserting a tube holder.

All Fahl® tracheostomy tubes with fastening eyelets include a disposable tube holder. The tube holder secures the tracheostomy tube to the neck.

Please follow the instructions provided with the tube holder carefully when attaching and removing it from your tracheostomy tube.

Ensure that the Fahl® tracheostomy tube is placed in the tracheostoma without tension and the position is not changed when fastening the tube holder.

The insertion aid (obturator) included in delivery facilitates insertion of the tracheostomy tube.

2. Connectors and adapters

Connectors and adapters are used to connect compatible tube accessories.

The options available to a patient depend on the clinical condition, such as status post laryngectomy or tracheostomy.

Connectors/adapters are attached to the outer cannula by means of a rotary fastener.

A 15 mm standard connector enables secure connection to required accessories.

3. Cannula tube

The cannula tube borders directly on the neck flange and directs the airflow into the trachea.

The neck flange is attached to the cannula tube and can be variably affixed to the tube by means of a screw.

The tube tip is rounded to prevent irritation of the mucus membranes in the trachea.

3.1 Cuff

In the product versions with cuff, the very thin-walled, high-volume cuff adapts well to the trachea and ensures reliable sealing if inflated correctly. The cuff can be inflated like a balloon. The small pilot balloon on the inflation tube indicates whether the tracheostomy tube is in sealed (inflated) or non-sealed condition.

The cuff itself is inflated via a tube with one-way valve and pilot balloon.

3.1.1 Leakproofness test of the cannula and the cuff (if such exists)

The cannula and the cuff must be tested for leakages directly before and after every insertion and thereafter at regular intervals.

For this purpose, fill the cuff with 15 to 22 mmHg (1 mmHg corresponds to 1.35951 cmH₂O) and observe whether a spontaneous pressure drop occurs.

There should be no substantial drop in pressure in the cuff during the observation period.

This leakproofness test must also be performed prior to every renewed insertion (e.g. after cleaning the tracheostomy tube) (see picture 7c).

The following are possible signs indicating leakages in the cuff (balloon):

- Externally visible damage to the balloon (holes, fissures etc.)
- Audible and/or tactile perception of air escaping from the balloon
- Water in inlet tubes of the tracheostomy tube (after cleaning!)
- Water in the cuff (after cleaning!)
- Water in the pilot balloon (after cleaning!)
- No cough stimulus when pressure is applied to the pilot balloon

CAUTION!

When testing the balloon or when inserting, removing or cleaning the tracheostomy tube, never under any circumstances use sharp or pointed objects, such as for instance forceps or clamps, since these can damage or even destroy the balloon. Should one of the above-named signs of leakage be detectable, the tracheostomy tube must under no circumstances be used any more since it is no longer functioning properly!

3.2 Obturator

The insertion aid serves to stabilise the tracheostomy tube during insertion into the tracheostoma. Please check whether the obturator can easily be removed from the cannula before inserting the tracheostomy tube!

After checking the obturator for ease of movement, reinsert the obturator into the cannula for inserting the tracheostomy tube.

SPIRAFLEX®/SPIRAFLEX® XL/SPIRAFLEX® MRT:

The insertion aid can also be used to absorb secretion fluids and convey these to the secretion container. At the end of the insertion aid, there is a tube connector/Fingertip, by means of which the suction tube of the suction device is connected to the insertion aid (see picture 12)

The Fingertip/tube connector can be used to control the suction power of the suction device.

Please read the corresponding instructions for use of the suction device carefully.

SPIRAFLEX® SHORT:

The obturator enables changing the cannula by applying the Seldinger technique.

For this purpose, use a guide wire with a maximum diameter of 0.038 inches.

CAUTION!

Make sure to read all accompanying product information, instructions for use, indications and contraindications. Discuss use of the product with your doctor/physician before first use.

Insertion aids are single-patient products and are intended for single use only. They must not be cleaned or disinfected.

3.3 Suction opening (only in Suction tracheostomy tube versions)

Secretions that have accumulated above the inflated cuff can be removed via the suction opening in the outer tube of the Suction tracheostomy tubes.

The suction tube leading to the outside can be connected to a syringe or a suction device. The user must decide in consultation with the treating physician and after appropriate risk analysis, which of these suctioning variants to use. The individual disease condition of the patient should be considered in this context.

Suctioning using a suctioning device may in any case only be performed if the suctioning device is equipped with vacuum adjuster. The suction pressure must never exceed - 0.2 bar.

When using the suction variants of the tracheostomy tubes, the secretions that have accumulated above the cuff must be suctioned off via the suction opening on the tracheostomy tube directly before deflating the cuff in order to prevent aspiration of secretions. In this way, additional suctioning with a suction catheter while deflating the cuff may be rendered unnecessary, which facilitates handling significantly.

CAUTION!

When using any of the tracheostomy tube versions with a suctioning device („SUCTION“), take special care during the suctioning procedure to generate a vacuum for as short a time as possible; drying out of the subglottic space as a result of this must be avoided.

CAUTION!

Contraindicated in patients with increased bleeding tendency (e.g. during anticoagulant therapy). In such cases, the Suction tracheostomy tube with suction opening may not be used, since suctioning entails an increased risk.

4. Inner cannula

The inner cannulas can easily be removed from the outer cannula, thus allowing the air supply to be increased quickly if necessary (for instance in case of dyspnoea).

Inner cannulas must never be used without outer cannula but must always be affixed to the outer cannula.

4.1 Speaking valves

Tracheostomy tubes with speaking valve (LINGO/PHON) are used after tracheotomy with a complete or partially retained larynx and enable the user to speak.

5. Decannulation plug

The decannulation plug is included in delivery of the tracheostomy tubes with speaking valve and may only be used by tracheotomised patients who have retained their larynx. It may only be inserted under supervision by a doctor/physician. It allows the tracheostomy tube to be sealed and the air supply to be interrupted for a short time and helps the patient train to regain control of breathing via mouth/nose.

CAUTION!

Contraindicated in laryngectomised patients and patients with chronic obstructive pulmonary disease (COPD)! In such cases, the decannulation plug must not under any circumstances be inserted!

The decannulation plug must never under any circumstances be used when the tracheostomy tube is in sealed condition! The decannulation plug may only be used with perforated outer cannula without inner cannula.

CAUTION!

A decannulation plug is enclosed with the tracheostomy tubes with speaking valve for the purpose of preparing for possible decannulation following temporary tracheotomy. The plug can be used to block the air supply via the tracheostomy tube for a short time so that the patient can get used to breathing via mouth/nose again. Decannulation may only be carried out under medical supervision. The plug may only be inserted upon instruction by the doctor/physician. There is a risk of suffocation! Please also make sure you observe the description of indications provided with the respective product variants/specifications!

CAUTION!

Make sure to read all accompanying product information, instructions for use, indications and contraindications. Discuss use of the product with your doctor/physician before first use.

VIII. DIRECTIONS FOR TUBE INSERTION AND REMOVAL**For the Doctor/Physician**

The appropriate tracheostomy tube must be selected by a doctor/physician or trained medical professionals.

Select a tube that fits the patient's anatomy to optimise comfort and ventilation (breathing in and out).

The inner cannula can be removed at any time to increase air supply or for cleaning. This can for instance be necessary if the cannula is clogged up with secretion residues which cannot be removed by coughing or because no equipment for suctioning off the secretions is available.

For the Patient**CAUTION!**

Always insert the tracheostomy tube with the cuff in fully deflated condition (see picture 7a)!

CAUTION!

Carefully examine the sterile packaging to ensure that it has not been tampered with or damaged. Do not use the product if the packaging has been damaged.

Check the use-by or expiry date. Do not use after this date.

It is advisable to use sterile disposable gloves.

Carefully examine the tube before first use to make sure that it is not damaged and that there are no loose parts.

Should you notice any anomaly or anything unusual, DO NOT use the tube. Return the tube to the manufacturer for inspection.

The tube must always be cleaned and, if necessary, disinfected as follows before re-inserting according to the instructions provided below.

If secretion collects in the lumen of the Fah!® tracheostomy tube or stoma button and cannot be removed by coughing or aspiration, the tube should be removed and cleaned.

After cleaning and/or disinfection, carefully examine the Fah!® tracheostomy tube for sharp edges, cracks, or other signs of damage, since these may impair function and/or injure the mucus membranes in the airways.

Never under any circumstances continue using damaged tracheostomy tubes.

1. Insertion of the tube**Step-by-step instructions to insert Fah!® tracheostomy tubes.**

Before application, users should clean their hands (see picture 3).

Remove tube from the package (see picture 4).

If an obturator is used, please note the following instructions:

SPIRAFLEX®/SPIRAFLEX® XL/SPIRAFLEX® MRT:

Insert the obturator into the cannula tube.

Make sure that the tip of the insertion aid

- protrudes only a few millimetres beyond the tip of the tracheostomy tube (proximal end of tube) when inserting the tube into the tracheostoma (see picture 13),
- If a suction device is connected, make sure that the tip protrudes so far beyond the tip of the tracheostomy tube (proximal end) that at least the first hole and at most the second hole (see picture 14) is fully visible.

The insertion aid must be held in this position during the entire procedure.

SPIRAFLEX® SHORT:

First fully insert the obturator into the cannula tube so that the collar on the gripping piece of the obturator comes into contact with the outer edge of the 15 mm connector.

The tip of the olive projects beyond the tip of the cannula (proximal end of cannula).

The insertion aid must be held in this position during the entire procedure.

When using tracheostomy tubes with cuff, pay special attention to the following points:

Before inserting the tracheostomy tube, check the cuff (balloon) as well - the cuff must be free of damages of any kind and must be leakproof so as to ensure tight sealing as required. We therefore recommend to perform a leakproofness test prior to every insertion (see Section VII, 3.1.1). The balloon must be completely empty prior to insertion of the tracheostomy tube (see picture 7b)! When using an aid for dilating the tracheostoma, take care to ensure that the tracheostomy tube, and especially the cuff, is not damaged by friction.

Then a tracheal compress is pushed onto the tracheostomy tube.

To increase the lubricity of the tracheostomy tube and thus facilitate insertion into the trachea, it is recommended to wipe the outer tube with an OPTIFLUID® stoma oil wipe (REF 31550), which allows even distribution of the stoma oil on the tube (see Figs. 4a and 4b) or Fah!® OPTIFLUID® Lubricant Gel 20g tube (REF 36100) or Fah!® OPTIFLUID® Lubricant Gel 3g sachet (REF 36105).

If you are inserting the tube yourself, use a mirror to make insertion of the Fah!® tracheostomy tube easier.

SPIRAFLEX®/SPIRAFLEX® XL/SPIRAFLEX® MRT:

The locking screw of the neck flange is only slightly tightened to ensure that the neck canal can be adjusted for optimum fit after the tracheostomy tube has been inserted.

When inserting the Fah!® tracheostomy tube, hold it by the neck flange with one hand (see picture 5).

Pull the tracheostoma slightly apart with your free hand to allow the tip of the tube to fit into the tracheostoma more easily.

Special aids for dilating the tracheostoma are also available. These allow the tracheostoma to be dilated gently and evenly, for instance also in emergency situations with collapsing tracheostoma (see picture 6).

Make sure that the tube is not damaged by friction when using an instrument for assistance. Now carefully insert the tube into the tracheostoma during the inspiration phase (while breathing in) while tilting your head slightly back (see picture 7).

Advance the tube into the trachea.

Straighten your head once the tube has been inserted further into the trachea.

Once the tube has been brought into the correct position in the trachea, remove the obturator and tighten the locking screw.

The tracheostomy tubes should always be attached with a special tube holder. This stabilises the tube and thus ensures that the tracheostomy tube is securely seated in the tracheostoma (see picture 1).

1.1 Inflating the cuff (if such exists)

The cuff is inflated by applying a defined pressure to the cuff via the Luer connection (standardised conical connection) of the inflation tube by means of a cuff pressure gauge. If not instructed otherwise by the doctor/physician, we recommend a cuff pressure of at least 15 mmHg (20 cmH₂O) to 18mmHg (25 cmH₂O). The cuff pressure should never under any circumstances exceed 18mmHg (approx. 25 cmH₂O).

Inflate the cuff at most to this target pressure and check to make sure that sufficient air is supplied via the tracheostomy tube.

Always make sure that the cuff is undamaged and is functioning faultlessly.

If the desired sealing is not achieved even after repeated attempts with the specified limit volume, a tracheostomy tube with larger diameter may be indicated

The correct cuff pressure must be checked regularly, i.e. at least every 2 hours.

CAUTION!

All instruments used for inflating the cuff must be clean and free of foreign particles! Detach the instruments from the Luer connection of the inflation tube as soon as the cuff has been inflated.

CAUTION!

If the maximum pressure is exceeded for longer periods of time, the blood circulation in the mucus membrane can be impaired (risk of ischaemic necrosis, pressure ulcers, tracheomalacia, tracheal stenosis, pneumothorax). In patients undergoing artificial respiration, the cuff pressure should not be allowed to drop below the cuff pressure value specified by the doctor/physician in order to prevent unnoticed aspiration. Hissing noises in the region of the balloon, especially during expiration, indicate that the trachea is insufficiently sealed by the balloon. If the trachea cannot be sealed with the pressure values specified by the doctor/physician, the entire air should be withdrawn again from the balloon and the sealing process should be repeated. If this does not lead to success, we recommend to use the next larger tracheostomy tube with balloon. Due to the permeability of the balloon wall for gases, it is normal for the pressure in the balloon to decline slightly over time, but it can on the other hand also rise unintentionally during gas anaesthesia. Regular pressure monitoring is therefore urgently recommended.

The cuff must never under any circumstances be inflated with excessive amounts of air, since this can lead to damage of the tracheal wall, tears in the cuff with subsequent deflation, or deformation of the cuff, in which case airway obstruction cannot be ruled out.

CAUTION!

During anaesthesia, the cuff pressure can rise/fall due to nitrous oxide (laughing gas).

2. Removing the tube

CAUTION!

Accessories such as a tracheostoma valve or HME (Heat Moisture Exchanger) must be removed first before proceeding to remove the Fahi[®] tracheostomy tube.

CAUTION!

If the tracheostoma is unstable, or in emergency situations (puncture/dilation tracheostomy), the tracheostoma can collapse after withdrawal of the tracheostomy tube, thereby impairing air supply. A fresh tracheostomy tube must be kept ready for use in such cases and must be quickly inserted if necessary. A tracheal dilator can be used for temporarily securing the air supply.

The cuff must be emptied before removing the tracheostomy tube. The head should be tilted back slightly for removal of the tube.

CAUTION!

Never use a cuff pressure gauge to empty the cuff. Always use a syringe for this.

Before the air is removed from the balloon by means of a syringe and the tracheostomy tube is withdrawn, the region of the trachea above the balloon must first be cleaned by suctioning off secretions and mucus. If the patient is responsive and reflexes are intact, it is recommended that the patient be suctioned while at the same time unblocking the tracheostomy tube. Suctioning is performed by inserting a suction catheter through the cannula tube into the trachea. In this way, suctioning can be performed without any problems and gently for the patient and cough stimulus and the risk of aspiration are minimised.

Then deflate the cuff while suctioning at same time.

If secretions are present, these are now removed and can no longer be aspirated. Please note that prior to reinsertion, the tracheostomy tube must always be cleaned, disinfected if necessary and lubricated with stoma oil as specified below.

Proceed very carefully to avoid injury to the mucus membranes.

Step-by-step instructions to remove the Fah!® tracheostomy tubes:

The tracheostomy tubes should be removed with the head slightly tilted back. Grip the tube at the side by the neck flange or the housing (see picture 7).

Carefully remove the tracheostomy tubes.

The outer cannula (if a cuff with inflated balloon is present) remains in the tracheostoma.

When removing the inner cannula, the following must be observed:

First, the accessories (adapter/speaking valve/cough cap/decannulation stopper) must be removed by slight counterclockwise rotation (in inserted condition from the patient's point of view) (see picture 15).

The inner cannula can now be removed from the outer cannula by pulling it gently (see picture 16).

The inner cannula is reinserted by means of slight rotary movements (see picture 17).

Reinsertion of the inner cannula is performed in reverse sequence as described above.

IX. CLEANING AND DISINFECTION

CAUTION!

For reasons of hygiene and to avoid the risk of infection, you should thoroughly clean the Fah!® tracheostomy tube at least twice a day, more often during heavier secretion production.

CAUTION!

The tracheostomy tubes must not be cleaned using a dishwasher, a steam cooker, a microwave oven, a washing machine or similar appliances!

Keep this in mind, your personal cleaning schedule, including additional disinfection if required, must always be determined in consultation with your doctor/physician and according to your individual needs.

Regular disinfection is only required if this is indicated from the medical point of view as ordered by the doctor/physician. The reason for this is that the upper airways are not free of microbes in healthy patients either.

General cleaning may be insufficient to meet the special hygiene requirements to prevent infections for patients with specific problems (e.g. MRSA, ORSA etc.), who are subject to an increased danger of re-infection. We recommend chemical disinfection of the tubes according to the instructions provided below. Please consult your doctor/physician.

CAUTION!

Cleaning agent and disinfectant residues on the tracheostomy tube can lead to irritations of the mucus membranes or other health impairments.

Tracheostomy tubes are considered as belonging to the category of instruments with hollow cavities. Special care must therefore be taken when performing cleaning or disinfection that the tube is fully wetted with the solution used and is unobstructed (see picture 8).

The solutions used for cleaning and disinfection must be freshly prepared prior to each use.

1. Cleaning

Fahl® tracheostomy tubes must be cleaned/replaced regularly in accordance with the individual patient's needs.

Only use cleaning agents if the tube is outside the tracheostoma.

A mild, pH-neutral washing lotion can be used for cleaning the tube. We recommend using the special cannula-cleaning powder (REF 31110) as directed by the manufacturer.

Never clean Fahl® tracheostomy tubes with cleaning agents that are not approved by the manufacturer. Never use aggressive household cleaners, detergents, high-concentration alcohol or denture cleaners.

This may be very hazardous to health. The tube may also be destroyed or damaged.

Tubes can also be cleaned by thermal disinfection at not more than 65 °C. Use clean water with a temperature of not more than 65 °C for this purpose. Make sure that the temperature is kept constant (use a thermometer to monitor the temperature) and never under any circumstances expose the tubes to boiling water for disinfection. This could seriously damage the tracheostomy tube.

Cleaning Steps

Cleaning the inner cannula

Do not use cleaning powder or a cleaning brush for the inner cannula!

Cleaning is carried out under running water. Gentle and careful kneading movements can be used to remove secretion residues (see picture 18).

Cleaning the tracheostomy tubes without cuff

Remove any additional items before cleaning.

The inner cannula must also be removed from the outer cannula.

First, rinse the tube thoroughly under running water (see picture 9).

Use only lukewarm water to prepare the cleaning solution and observe the instructions for use of the cleaning agent.

We recommend the use of the cannula-cleaning tub with sieve insert (REF 31200) to facilitate the cleaning procedure.

To do this, hold the sieve insert by the upper edge to prevent contact with and contamination of the cleaning solution (see picture 10).

Place only one tracheostomy tube at a time in the sieve insert of the cannula-cleaning tub. If more than one tube is cleaned at a time, they may be compressed too strongly and damaged by excessive pressure.

Immerse the sieve insert with the cannula in the prepared cleaning solution.

After the immersion time has expired (see instructions for use of the cannula cleaning powder), wash the tracheostomy tube thoroughly several times with lukewarm, clear water (see picture 9). There must be no residues of the cleaning agent on the tube when it is inserted into the tracheostoma.

If necessary, for example, if persistent secretion residues cannot be removed by the cleaning bath, additional cleaning with a special cannula-cleaning brush (OPTIBRUSH®, REF 31850 or OPTIBRUSH® PLUS with fibre top, REF 31855) may help. Only use the cleaning brush, if the tube is removed and already outside the tracheostoma.

Always insert the cannula cleaning brush into the cannula from the cannula tip (see picture 11).

Use the brush as directed and proceed with great care to avoid damage to the soft tube material. Thoroughly rinse the tube or button under lukewarm running water or sterile saline solution (0.9% NaCl solution).

Dry the tube thoroughly with a clean and lint-free cloth after wet cleaning.

Never use tracheostomy tubes with impaired functionality or with damages such as sharp edges or cracks, because this may lead to injuries of the mucus membranes in the trachea. If damages are found, the tube must not be used under any circumstances.

The accessories (adapter/speaking valve/cough cap/decannulation stopper) can be cleaned in the same way as the tracheostomy tube.

Insertion aids are single-patient products and are intended for single use only. They must not be cleaned or disinfected.

Cleaning the tracheostomy tubes with cuff

The indwelling time and thus the replacement intervals of the tracheostomy tube must be determined individually in consultation with the treating doctor/physician.

The outer cannula must be cleaned or the tracheostomy tube replaced after 1 week at the latest though, since after that the risk of formation of granulation tissue or of tracheomalacia, etc. can increase significantly, depending on the patient's disease condition.

The outer cannula must only be cleaned with the cuff inflated, and the safety balloon must not come into contact with the cleaning solution in order to prevent intrusion of cleaning solution into the balloon, since this would lead to serious functional impairment and health risks for the user. The outer cannula with cuff can be cleaned/rinsed with sterile saline solution.

A brush must not be used for tracheostomy tubes with cuff in order to avoid damaging the balloon!

Damage to the balloon and hence a defect of the tracheostomy tube can only be avoided by careful and gentle handling.

2. Chemical Disinfection Instructions

2.1 Disinfecting the inner cannula/cleaning the outer cannula without cuff

The Fah!® tracheostomy tube can be disinfected by cold disinfection with special chemical disinfectants.

Disinfection should always be done if ordered by a doctor/physician due to specific health concerns caused by disease, infection, or your specific situation.

Disinfecting is generally indicated to prevent cross-infection and in inpatient situations (e.g. hospitals, nursing homes, and/or other health care facilities) to limit infections.

CAUTION

Always clean according to the cleaning procedure described above before proceeding to disinfection (if applicable).

Disinfectants that release chlorine or that contain strong alkalis or phenol derivatives must under no circumstances be used. This could badly damage or even destroy the tracheostomy tube.

2.2 Disinfecting the outer cannula with cuff

Tracheostomy tubes with cuff may only be disinfected provided that this is performed and checked with utmost care. The balloon must always be inflated beforehand.

Disinfection steps

For this purpose, OPTICIT® tube disinfectant (REF 31180) should be used in accordance with the manufacturer's instructions.

As alternative, we recommend a disinfectant based on glutaric aldehyde as active ingredient (also available in the USA). Always observe the area of application and spectrum of activity specified for the disinfectant by the manufacturer.

After disinfection the tube must be rinsed inside and out very thoroughly with sterile saline solution (NaCl 0.9%) and then dried.

SPIRAFLEX®/SPIRAFLEX® XL/SPIRAFLEX® MRT:

When doing this, the neck flange should be de-arrested and moved to and fro several times to ensure thorough cleaning underneath the locking screw of the neck flange as well

Dry the tube thoroughly with a clean and lint-free cloth after wet cleaning.

3. Sterilisation/Autoclaving

Resterilisation is not allowed.

CAUTION!

Heating to over 65°C, disinfection with boiling water or steam sterilisation are not allowed and lead to damage of the tracheostomy tube.

X. STORAGE/CARE

Currently not used and cleaned tubes should be stored in a clean plastic container, dry environment away from dust, sunlight and heat.

Still sterile packed replacement tubes should be stored in a dry environment away from dust, sunlight and heat.

After cleaning tracheostomy tubes, make sure that the balloon has been deflated for storage.

After the inner cannula has been cleaned, disinfected if necessary and then dried, the outer surface of the inner cannula should be lubricated by applying stoma oil (OPTIFLUID® Stoma Oil, 25 ml bottle REF 31525/Stoma oil wipe REF 31550) or lubricant gel (Fah!® OPTIFLUID® Lubricant Gel 20 g tube REF 36100/3 g sachet REF 36105).

We highly recommend keeping at least two replacement tubes on hand to ensure continuous supply.

EN

XI. SERVICE LIFE

These tracheostomy tubes are sterile products for single-patient use.

The maximum period of use is 29 days (reckoned from the date on which the sterile package was opened). If a pathogenic microbe (e.g. MRSA) is present, the period of use/indwelling time is reduced accordingly.

The service life of a tube is influenced by many factors. The composition of the secretions, the quality of cleaning and other aspects, for example, are very important.

The maximum service life is not extended by using the tube intermittently (e.g. alternately with other tubes at normal cleaning intervals).

Damaged tubes must be replaced immediately.

CAUTION!

Modifications to the tracheostomy tube of any kind as well as repairs on the tube must only be carried out by the manufacturer himself or by companies expressly authorised for this in writing by the manufacturer! Inexpertly performed work on tracheostomy tubes can lead to serious injuries.

XII. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the tracheostomy tube, especially due to shortening of the same or perforations, if these modifications or repairs were not carried out by the manufacturer himself. This applies both to damages to the tubes themselves caused thereby and to any consequential damages caused thereby.

If the tracheostomy tube is used beyond the period of use specified under Section XII., and/or if the tracheostomy tube is used, maintained (cleaned, disinfected) or stored in non-compliance with the instructions and specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious adverse event occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, this is to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.


Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

Product specifications described herein are subject to change without notice.

SPIRAFLEX® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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
ANHANG APPENDIX

SPIRAFLEX® UNI/SPIRAFLEX® MRT UNI/SPIRAFLEX® UNI CUFF XL/SPIRAFLEX® UNI CUFF/
SPIRAFLEX® UNI CUFF SUCTION/SPIRAFLEX® MRT UNI CUFF/SPIRAFLEX® SHORT UNI/
SPIRAFLEX® SHORT UNI CUFF/SPIRAFLEX® SHORT UNI LINGO PHON CUFF/SPIRAFLEX®
SHORT UNI CUFF SUCTION/SPIRAFLEX® SHORT UNI LINGO PHON

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All measurements listed in the tables below are for standard cannulas.
Customised cannulas are available upon request.

GRÖSSENTABELLE - SIZE TABLE							
GR	AD1	AD2	ID1	ID2	L1-AB	L2-MB	∅ BW
06	8,4	8,4	6,0	4,5	60,0	52,0	80°-110°
07	10,2	10,2	7,0	5,5	80,0	68,0	80°-110°
08	11,2	11,2	8,0	6,5	100,0	88,0	80°-110°
09	12,2	12,2	9,0	7,5	130,0	117,0	80°-110°
10	13,2	13,2	10,0	8,5	135,0	118,0	80°-110°
11	14,2	14,2	11,0	9,5	140,0	123,0	80°-110°

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Sizing information:

All lengths in mm | GR = Size | AD1 = Outer diameter at the neck flange | AD2 = Outer diameter at cannula tip | ID1 = Inner diameter at the cannula tip of the outer cannula | ID2 = Inner diameter at the cannula tip of the inner cannula | L1-AB = Length over the outer curve | L2-MB = Length over the middle curve | ø-CD = CUFF diameter | Ø-BW = Bending angle

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Sizing information:

All lengths in mm | GR = Size | AD1 = Outer diameter at the neck flange | AD2 = Outer diameter at cannula tip | ID1 = Inner diameter at the cannula tip of the outer cannula | ID2 = Inner diameter at the cannula tip of the inner cannula | L1-AB = Length over the outer curve | L2-MB = Length over the middle curve | Ø-CD = CUFF diameter | Ø-BW = Bending angle

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Sizing information:

All lengths in mm | GR = Size | AD1 = Outer diameter at the neck flange | AD2 = Outer diameter at cannula tip | ID1 = Inner diameter at the cannula tip of the outer cannula | ID2 = Inner diameter at the cannula tip of the inner cannula | L1-AB = Length over the outer curve | L2-MB = Length over the middle curve | ∅-CD = CUFF diameter | ∅-BW = Bending angle

SPIRAFLEX® SHORT UNI CUFF/SPIRAFLEX® SHORT UNI LINGO PHON CUFF

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Sizing information:

All lengths in mm | GR = Size | AD1 = Outer diameter at the neck flange | AD2 = Outer diameter at cannula tip | ID1 = Inner diameter at the cannula tip of the outer cannula | ID2 = Inner diameter at the cannula tip of the inner cannula | L1-AB = Length over the outer curve | L2-MB = Length over the middle curve | ø-CD = CUFF diameter | θ-BW = Bending angle

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