







FAHL® DECANNULATION TAPE

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I. FOREWORD

These instructions for use are valid for FAHL® DECANNULATION TAPE. The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling. Please read the instructions for use carefully before using the product!

II. INTENDED USE

FAHL® DECANNULATION TAPE with cuff is used for the temporary and secure sealing of the tracheostoma in tracheostomised patients.

III. WARNINGS

FAHL® DECANNULATION TAPE may only be inserted upon instruction by the doctor/physician. There is a risk of suffocation!

The FAHL® DECANNULATION TAPE with cuff may ONLY be used for tracheostomised patients!

IV. CONTRAINDICATIONS

The use of FAHL® DECANNULATION TAPE with cuff is not permitted in laryngectomised patients!

In case of increased secretion formation or pronounced vocal cord paralysis (recurrent paresis), the FAHL® DECANNULATION TAPE with cuff must NOT be used!

V. COMPLICATIONS

The following complications can occur prior to, during or after decannulation and can impede decannulation or render decannulation impossible: granulations, aspiration, infections as well as strong hypersalivation and atelectases, temporary sensation of shortness of breath after applying the decannulation patch and inflating the cuff.

VI. INSTRUCTIONS

1. Clean the skin around the tracheostoma in front of a mirror using warm water and pH-neutral soap or OPTIFAHL® stoma cleaning cloths (REF 33200) . Dry the skin. Before cleaning, remove residues of adhesive from the skin if necessary, using OPTICLEAR®(REF 33500)
2. To protect the skin from irritation, apply OPTIGARD® to the affected area of the skin.
3. Before applying the FAHL DECANNULATION TAPE to the tracheostoma, it is recommended to check the cuff for leakage.
To do this, the cuff is filled with 10 ml of air using a disposable syringe (volume: 20 ml) via the filling tube ❶.
The cuff takes on the shape of a hemisphere ❷.
The shape of the cuff must remain stable and no air must escape
If the cuff loses air, this FAHL® DECANNULATION TAPE is not to be used and must be replaced.
Once the cuff is stable, the air is then extracted from the cuff again via the filling tube with non-return valve with the aid of the disposable syringe and the FAHL® DECANNULATION TAPE can be used.

4. After removing the protective foil backing ❶, place the FAHL® DECANNULATION TAPE with the adhesive side centred over the tracheostoma.
5. To ensure optimal fixation and to avoid air enclosures, carefully stroke the entire tape surface with your fingers for 5 minutes after attaching the tape to the skin.
After placement of the patch, the cuff serves to seal the tracheostoma.
Due to the low pressure level which needs to be achieved for sealing, pressure necrosis can be avoided
The cuff can be filled with air using a disposable syringe (volume: 20 ml) via the filling tube with non-return valve ❷.
6. First fill the cuff with 10 ml of air.
Subsequent vocalisation by the patient can be used to check whether the cuff is sufficiently tight on the tracheostoma.
If there is no background noise, the cuff is sufficiently filled with air.
If background noise occurs during vocalisation, the cuff must be further filled with air until the tracheostoma is tight and vocalisation is performed without background noise.
Make sure that a total filling quantity of the cuff of 20 ml maximum is NOT exceeded!
7. In case of increased shortness of breath, relief can be provided by briefly unblocking the cuff or by reducing pressure via the filling valve using a disposable syringe.
8. To remove the FAHL® DECANNULATION TAPE, carefully grip an adhesive tab.
Carefully peel this off from the skin on one side. If necessary use OPTICLEAR® (REF 33500).

CAUTION!

Do not use creams that contain fat! Otherwise, secure attachment of the tape is no longer ensured.

VII. SERVICE LIFE

FAHL® DECANNULATION TAPE is non-sterile and intended for single use only.

Replacement should be performed after 5 days at the latest, unless an earlier replacement is indicated by the clinical picture.

When skin irritations occur, the LARYVOX® Tape must be changed immediately. Deviations from the expected wearing period occur as a result of individual wearing conditions and skin properties.

CAUTION!

Cleaning, disinfection, re-sterilisation and re-use can impair the safety and function of the product and are therefore not permissible!

VIII. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or other adverse events caused by unauthorised product alterations or improper use, care, and/or handling. If the FAHL® DECANNULATION TAPE neck tapes are used beyond the period of use specified under Section VII, and/or if they are used or maintained (cleaned) in non-compliance with the instructions and specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious adverse event occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, this is to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

Product specifications described herein are subject to change without notice.

FAHL® DECANNULATION TAPE is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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