# LARYVOX® INSERTER

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
# PICTROGRAM LEGEND

Pictogrammes listed below you'll find on the product packaging, if applicable.

EN

	Catalogue number
	Batch code
	Use by
	Date of manufacture
	Content (in pieces)
	Consult instructions for use
	Keep away from sunlight
	Keep dry
	Do not reuse
	Medical devices
	Do not re-sterilize
	Do not use if package is damaged
	CE marking
	Manufacturer

## APPLICATION TABLE

25900	25910	25915	25920	25925
<b>LARYVOX® INSERTER</b>	<b>LARYVOX® INSERTER EXPERT</b>	<b>LARYVOX® INSERTER EXPERT SMALL</b>	<b>LARYVOX® INSERTER SPECIAL</b>	<b>LARYVOX® INSERTER UNIEK</b>
<b>Colour</b>				
turquoise	black	yellow	blue	orange
				
<b>Assignment</b>				
Voice prostheses 20 Fr. or 22.5 Fr. ≥ 5 mm shaft length	Voice prostheses 16 Fr. or with dual valve ≥ 5 mm shaft length	Voice prostheses 16 Fr. or 20 Fr. with 4 mm shaft length	Changing a voice prosthesis with the gel capsule technique	LARYVOX® Uniek Ultra Low Resistance voice prosthesis 21 Fr. or 24 Fr.

# LARYVOX® INSERTER

## 1. FOREWORD

EN

These instructions apply to **LARYVOX® INSERTER (-EXPERT, -EXPERT SMALL, -UNIEK, -SPECIAL)**.

These instructions for use serve to inform the physician and medical professional staff in order to ensure technically correct handling of **LARYVOX® INSERTERS**.

Please read the instructions for use carefully before using the product for the first time!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the **LARYVOX® INSERTER**. It contains important information on the product!

## 2. SAFETY INFORMATION

### ⚠ WARNING!

The **LARYVOX® INSERTER (-EXPERT, -EXPERT SMALL, -UNIEK, -SPECIAL)** is a single-patient product and is intended for one single use only.

### ⚠ WARNING!

Never use the **LARYVOX® INSERTER (-EXPERT, -EXPERT SMALL, -UNIEK, -SPECIAL)** if the device exhibits visible deficiencies or defects.

Contact the manufacturer in such cases.

### ⚠ WARNING!

A damaged voice prosthesis must not be used.

### ⚠ WARNING!

The voice prosthesis must always be placed with care in the shunt channel. Pay attention to the correct position of the voice prosthesis.

### ⚠ WARNING!

A **LARYVOX® INSERTER (-EXPERT, -EXPERT SMALL, -UNIEK, -SPECIAL)** must not be cleaned or disinfected.

Sterilisation and reuse can impair the function and is therefore not permissible either.

### CAUTION!

Please also observe the notes / warnings in the instructions for use provided by the voice prosthesis manufacturer.

## 3. INTENDED PURPOSE

**LARYVOX® INSERTER (-EXPERT, -EXPERT SMALL, -UNIEK, -SPECIAL)** devices are used for the antegrade placement of a voice prosthesis when replacing the voice prosthesis.

## 4. INDICATION

A voice prosthesis must not remain permanently in the shunt channel, but needs to be replaced at regular intervals, depending on its proper functioning.

The leakage of fluids with the transfer of saliva or food particles from the oesophagus to the trachea is a direct indication for an examination by the physician and, where necessary, replacement of the voice prosthesis used.

The **LARYVOX® INSERTER (-EXPERT, -EXPERT SMALL, -UNIEK, -SPECIAL)** replacement is performed in patients who have been equipped with a voice valve following removal of the larynx (laryngectomy).

## 5. CONTRAINDICATION

The use of the **LARYVOX® INSERTER (-EXPERT, -EXPERT SMALL, -UNIEK, -SPECIAL)** for voice prosthesis placement in patients with increased bleeding tendency, e.g. due to treatment with anticoagulants, must always be weighed carefully on account of the risk.

Special caution is required in patients following radiation therapy and/or chemotherapy, as there is an increased risk of puncture-induced complications due to possible tissue damage (e.g. scar tissue, fibroses and the like).

## 6. COMPLICATIONS AND ADVERSE REACTIONS

### ⚠ WARNING!

Inappropriate handling of the **LARYVOX® INSERTER** can cause material damage to the voice prosthesis.

### CAUTION!

When the voice prosthesis is inserted, this can lead to irritations, coughing or slight bleeding.

### ⚠ WARNING!

Should the voice prosthesis come loose from the shunt channel or during insertion of the placement instrument, there is a risk of aspiration.

In the event of aspiration, it must be attempted to expel the voice prosthesis through strong coughing.

If the voice prosthesis cannot be coughed out, it must immediately be removed surgically.

### Further complications:

- General discomfort due to manipulation on the shunt channel

- Expansion/dilatation of the shunt channel
- Inflammatory reactions and formation of granulations in the region of the shunt channel

## 7. PRODUCT DESCRIPTION

The **LARYVOX® INSERTER (-EXPERT, -EXPERT SMALL, -UNIEK, -SPECIAL)** consists of a placement instrument ① for the secure fixation of the voice prosthesis for insertion with the applicator\* (\*with the exception of **LARYVOX® INSERTER SPECIAL**) ② into the shunt channel.

The **LARYVOX® INSERTER** is available in the variants **-EXPERT, -EXPERT SMALL, -UNIEK, -SPECIAL**.

Each variant has a distinct colour (see list) so that it can be clearly distinguished.

The inserter enables users of voice prostheses to apply a uniform placement technique for the insertion of different types of voice prostheses.

### NOTE:

In addition, also heed the instructions for use of the voice prosthesis used.

## 8. APPLICATION

### 8.1 BASIC PREPARATION

1. Prior to removing the voice prosthesis to be replaced, check its seating in the shunt channel.  
The correct length of the voice prosthesis is of essential importance for avoiding complications such as peri-prosthetic leakage, pressure necroses or the like.
2. It is recommended to use the **BLOM-SINGER® SIZER DILATOR (REF BE2051, BE2052)** for checking the length of the shunt channel in order to determine the correct length of the shunt channel prior to placement.
3. Check the packaging for damages and check the expiry date. Do not use the product after expiry of this date, or if the packaging and/or the product exhibit visible damages.

### 8.2 APPLICATION OF THE VARIANTS

The application of the various **LARYVOX® INSERTERS** is detailed for each variant in the following. Please follow the relevant instructions:

#### 8.2.1 LARYVOX® INSERTER REF 25900 (turquoise)

##### PLEASE NOTE!

The **LARYVOX® INSERTER** is used only and exclusively for the antero-grad placement / replacement of a 20 or 22.5 French voice prosthesis, with a shaft length of 5 mm or longer.

The **LARYVOX® INSERTER** must not be used for placement of a **BLOM-SINGER® DUAL-VALVE INDWELLING** device (REF DV20XX-NS, REF DV20XX-LF).

##### ⚠ WARNING!

**A voice prosthesis with a shaft length of less than 5 mm or with a dual valve flap cannot be securely fixed in place with the LARYVOX® INSERTER.**

**There is a risk that such a voice prosthesis may become detached from the inserter and may cause acute respiratory distress.**

#### 8.2.2 LARYVOX® INSERTER EXPERT (REF 25910) (black)

##### PLEASE NOTE!

The **LARYVOX® INSERTER EXPERT** is used only and exclusively for the antero-grad placement of the following voice prostheses with a shaft length of 5 mm or longer:

Voice prostheses with a diameter of 16 French, also **BLOM-SINGER® DUAL VALVE INDWELLING** (with dual valve flap) - **BLOM-SINGER® DUAL VALVE INDWELLING 20 French**

##### ⚠ WARNING!

**A voice prosthesis with a shaft length of less than 5 mm cannot be securely fixed in place with the LARYVOX® INSERTER EXPERT. There is a risk that such a voice prosthesis may become detached from the inserter and may cause acute respiratory distress.**

#### 8.2.3 LARYVOX® INSERTER EXPERT SMALL (REF 25915) (yellow)

##### PLEASE NOTE!

The **LARYVOX® INSERTER SMALL** is used only and exclusively for the antero-grad placement / replacement of voice prostheses with a shaft length of 4 mm, 16 French or 20 French.

##### ⚠ WARNING!

**A voice prosthesis with a shaft length of more than 4 mm cannot be securely fixed in place with the LARYVOX® INSERTER SMALL.**

**There is a risk that such a voice prosthesis may become detached from the inserter and may cause acute respiratory distress.**

**⚠ WARNING!**

**EN**

Should the voice prosthesis come loose from the shunt channel or during insertion of the placement instrument, there is a risk of aspiration.

In the event of aspiration, it must be attempted to expel the voice prosthesis through strong coughing.

If the voice prosthesis cannot be coughed out, it must immediately be removed surgically.

**INSTRUCTIONS FOR**

**LARYVOX® INSERTER (TURQUOISE, REF 25900)**

**LARYVOX® INSERTER EXPERT (BLACK, REF 25910) AND**

**LARYVOX® INSERTER EXPERT SMALL (YELLOW, REF 25915):**

1. Mount the voice prosthesis with the tracheal flange on the head of the placement instrument ③.
2. Fix the safety thread of the voice prosthesis in place with the opening at the intended position:

**LARYVOX® INSERTER** (turquoise) position B ⑤.

**LARYVOX® INSERTER EXPERT** (black) - if the safety thread cannot be securely fixed in place at position B ⑤, then choose position A ④.

**LARYVOX® INSERTER EXPERT SMALL** (yellow) position A ④.

3. Check to ensure that the voice prosthesis is securely seated on the placement instrument.
4. To achieve better lubrication of the applicator for advancing the voice prosthesis, you can use the **FAHL® OPTIFLUID® LUBRICANT GEL (REF 36105)** included in delivery.
5. Press the oesophageal flange together in forward direction with the thumb and index finger of one hand ⑦.  
Push the oesophageal flange forward by approx. 1 cm into the slotted opening of the applicator using the placement instrument ⑧.
6. While doing this, press on the oesophageal flange with the thumb of the other hand ⑨.  
At the same time, push the voice prosthesis further into the applicator with your other hand ⑩.
7. Marked positions 1 and 2 are shown on the placement instrument.  
Now guide the placement instrument with the voice prosthesis through the applicator up to the line of marked position 1 ⑪.
8. The voice prosthesis can now be introduced into the shunt channel using the applicator.
9. Cautiously introduce the tip of the applicator into the shunt channel ⑫.
10. Now slowly advance the voice prosthesis through the applicator up to the line of marked position 2.  
The oesophageal flange will open in the lumen of the oesophagus ⑬.  
Now cautiously retract the applicator via the placement instrument ⑭.
11. The tracheal flange opens in the lumen of the trachea.  
The applicator can be retracted via the placement instrument
12. The correct position of the voice prosthesis is verified by rotating the voice prosthesis with the placement instrument and by pulling slightly ⑮.
13. If the voice prosthesis can be easily rotated, fix the safety thread in place using an atraumatic instrument such as e.g. **LARYVOX® TWEEZERS (REF 35000)** and cautiously remove the placement instrument ⑯.
14. Perform a phonation test and, where appropriate, a swallowing test with the patient in order to verify the leakproofness.
15. If the voice prosthesis is in the correct position, cut the safety thread flush with the tracheal flange.

**8.2.4 LARYVOX® INSERTER UNIEK REF 25925 (orange)**

**PLEASE NOTE!**

The **LARYVOX® INSERTER UNIEK** is used only and exclusively for the anterograde placement of the **LARYVOX® UNIEK** Ultra Low Resistance (**REF 25010 and REF 25020**) voice prosthesis, 21 or 24 French, in the scope of voice prosthesis replacement.

**⚠ WARNING!**

Should the voice prosthesis come loose from the shunt channel or during insertion of the placement instrument, there is a risk of aspiration.

In the event of aspiration, it must be attempted to expel the voice prosthesis through strong coughing.

If the voice prosthesis cannot be coughed out, it must immediately be removed by professional medical means.

**INSTRUCTIONS:**

1. Fit the voice prosthesis with its tracheal flange onto the head of the placement instrument ③.
2. Fix the safety thread of the voice prosthesis in place on the placement instrument at position A ④.
3. Check to ensure that the voice prosthesis is securely seated on the placement instrument.

4. To achieve better lubrication of the applicator for advancing the voice prosthesis, it is recommended to use the **FAHL® OPTIFLUID® LUBRICANT GEL** included in delivery.
5. Use the placement instrument to push the voice prosthesis forward by approx. 1 cm into the slotted opening of the applicator ⑤.
6. Marked positions 1 and 2 are shown on the placement instrument.  
Now guide the placement instrument with the voice prosthesis through the applicator up to marked position 1 ⑥.
7. The voice prosthesis can now be introduced into the shunt channel using the applicator.
8. Cautiously introduce the tip of the applicator into the shunt channel ⑦.
9. Now slowly advance the voice prosthesis through the applicator up to marked position 2.  
The oesophageal flange will open in the lumen of the oesophagus ⑧.
10. Now cautiously retract the applicator via the placement instrument.  
The tracheal flange opens in the lumen of the trachea.  
The applicator can be retracted via the placement instrument ⑨.
11. The correct position of the voice prosthesis is verified by rotating the voice prosthesis with the placement instrument or by pulling slightly ⑩.
12. Rotate the voice prosthesis into position such that the safety thread is in the 12 o'clock position.
13. Fix the safety thread in place with an atraumatic clamp and remove the placement instrument ⑪.
14. Perform a phonation test. Check the leakproofness of the voice prosthesis with a swallowing test.
15. If the voice prosthesis is in the correct position, cut the safety thread flush with the tracheal flange ⑫.

### 8.2.5 LARYVOX® INSERTER SPECIAL REF 25920 (blue)

The **LARYVOX® INSERTER SPECIAL** consists of a placement instrument with a curved base ① for introducing a voice prosthesis into the shunt channel which is to be introduced using the gel capsule technique.

#### FIELDS OF APPLICATION:

- Fistula in low position or with oblique course
- Constricted shunt channel/fistula
- Voice prosthesis replacement in particularly small tracheostomas to optimise viewing conditions

#### PLEASE NOTE!

The **LARYVOX® INSERTER SPECIAL** is used only and exclusively for the antegrade placement of a voice prosthesis using the gel capsule technique in the scope of voice prosthesis replacement.

Heed the instructions for use of the voice prosthesis as well as these instructions for use for inserting voice prostheses using the gel capsule technique!

#### INSTRUCTIONS:

1. After folding the voice prosthesis into the gel capsule, mount it with its tracheal flange on the head of the placement instrument ②.
2. Fix the safety thread of the voice prosthesis in place at the intended position ③ + ④.
3. Check to ensure that the voice prosthesis is securely seated on the placement instrument.

#### ⚠ WARNING!

**Should the voice prosthesis come loose from the shunt channel or during insertion of the placement instrument, there is a risk of aspiration.**

**In the event of aspiration, it must be attempted to expel the voice prosthesis through strong coughing.**

**If the voice prosthesis cannot be coughed out, it must immediately be removed surgically.**

4. The voice prosthesis can now be introduced into the shunt channel ⑤.
5. The correct position of the voice prosthesis is verified after dissolution of the gel capsule by rotating the voice prosthesis with the placement instrument or by pulling slightly ⑥.
6. If the voice prosthesis can be easily rotated, fix the safety thread in place using an atraumatic instrument such as e.g. **LARYVOX® TWEEZERS (REF 35000)** and cautiously remove the placement instrument ⑦.
7. Perform a phonation test and, where appropriate, a swallowing test with the patient in order to verify the leakproofness.
8. If the voice prosthesis is in the correct position, cut the safety thread flush with the tracheal flange ⑧.

## 9. CLEANING & DISINFECTION

#### PLEASE NOTE!

**LARYVOX® INSERTERS (-EXPERT, -EXPERT SMALL, -UNIEK, -SPECIAL)** must not be cleaned or disinfected. Sterilisation and reuse can impair the function and is therefore not permissible either.

## 10. STORAGE

EN

**LARYVOX® INSERTERs** (-EXPERT, -EXPERT SMALL, -UNIEK, -SPECIAL) should be stored in a dry environment and protected from sunlight and/or heat.

## 11. SERVICE LIFE

**LARYVOX® INSERTER** (-EXPERT, -EXPERT SMALL, -UNIEK, -SPECIAL) is a single-patient product and is intended for one single use only.

## 12. DISPOSAL

**LARYVOX® INSERTER** (-EXPERT, -EXPERT SMALL, -UNIEK, -SPECIAL) must only be disposed of in accordance with the applicable national regulations.

After use, the product entails a possible biological hazard potential.

When disposing of medical devices, medical practice as well as the applicable laws and regulations of the respective country must be observed.

## 13. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the **LARYVOX® INSERTERs** (-EXPERT, -EXPERT SMALL, -UNIEK, -SPECIAL) or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies both to damages to the **LARYVOX® INSERTERs** themselves caused thereby and to all consequential damages caused thereby.

If the **LARYVOX® INSERTER** is used beyond the period of use specified in Section 11, and/or if the **LARYVOX® INSERTER** is used, maintained (cleaned, disinfected) or stored in non-compliance with the specifications laid down in these instructions for use, then Andreas Fahl Medizintechnik-Vertrieb GmbH shall, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.


The manufacturer reserves the right to make changes to the product at any time without notice.

**LARYVOX®** is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the European Member States.

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