





FAHL

HUMIDOTWIN®

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PICTOGRAMM LEGEND

EN

 REF	Catalogue number		Use by
 LOT	Batch code		Do not reuse
	Consult instructions for use		Content (in pieces)
	Keep away from sunlight		Do not use if package is damaged
	Keep dry		Manufacturer

HUMIDOTWIN®

I. FOREWORD

These instructions for use are valid for HUMIDOTWIN® HME (Heat-Moisture Exchanger).

The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling. Please read the instructions for use carefully before using the product for the first time!

II. INTENDED USE

HMEs filter the breathing air, thereby reducing intrusion of particles into the patient's airways.

The filter cassettes reduce the formation of viscous secretions in the lungs and protect the patient against infections. They store the moisture and warmth of the exhaled air in the filter medium and emit these to the inhaled air again during inspiration.

III. PRODUCT DESCRIPTION

HMEs are heat and moisture exchangers, as well as particle filters, consisting of a foam filter and a plastic housing.

The plastic housing is equipped on the patient side with a central opening with 15 mm internal diameter, thereby allowing for connection to a 15 mm standard connector.

The external diameter of the connector is 22 mm, and is thus also compatible with all commonly used 22 mm receptacle systems.

A small rim on the side of the opening close to the patient enables secure fixation of the HME in combi-adapters.

The overpressure valve in the lid of the housing prevents an excessive pressure build-up and allows secretions to be coughed up without having to remove the HME.

It minimises the risk of the moistening medium sticking and allows larger quantities of secretion to be coughed up.

IV. WARNINGS

HMEs are single-patient products and are intended for single use only.

It must not be cleaned or disinfected, as use of cleaning agent can render the filter unusable.

Sterilisation and reuse impair the function and therefore are also not permissible.

HMEs must under no circumstances be inserted directly into the tracheostoma (risk of suffocation!).

The product is suitable only for patients capable of breathing spontaneously.

Care must be taken to ensure that the valve opening is not blocked by clothing or the like in order to enable unhindered breathing (risk of respiratory distress).

The Artificial Noses must not be used in conjunction with a nebuliser or atomiser.

Do not fill any water into the HME.

This would always cause the danger of increased breathing resistance and thus a breathing obstruction.

V. PRECAUTIONS

If the filter is clogged or if signs of clogging appear, the HME must immediately be removed and disposed of.

HUMIDOTWIN® cannot be used in conjunction with TRACHEOTEC®, TRACHEOTEC® PRO / SILC and TRACHEOSILC®!

VI. CONTRAINDICATIONS

Use in patients with strong exsiccosis (dehydration), or excessive secretion and water accumulation in the lungs and airways (left-ventricular heart failure, pulmonary oedema), or strong air leakage flow (exhaled air not flowing through the HME) is not permissible.

VII. INSTRUCTIONS

- Check whether the packaging and the product are in a faultless condition.
- Check the use-by or expiry date. Do not use after this date.
- Use with a 15 mm connector: Mount the HME with the central 15 mm opening by means of light pressure on the 15 mm standard connector of the tracheostomy tube / endotracheal tube ❶.
- Use with a 22 mm connector: Mount the HME by means of light pressure on the housing ring of the tracheostomy tube / endotracheal tube ❷, until the housing snaps in audibly through the outer fastening rim.
- Remove the HME by holding the tracheostomy tube / endotracheal tube firmly with two fingers in order to extract the filter cassette by tugging lightly ❸/❹.

VIII. SERVICE LIFE

The maximum period of use / wearing period is 24 hours ❺.

The filter must be changed even several times a day if required, for instance if the respiratory resistance is in danger of increasing due to accumulation of secretions.

IX. TECHNICAL DATA (DIN EN ISO 9360-2)

Tidal volume (VT):	min. 250 ml
Internal Volume:	3,3 ml
Moisture lost:	VT = 250 ml : 8,5 mg/l VT = 500 ml : 10,5 mg/l VT = 750 ml : 13,0 mg/l VT = 1000 ml : 17,1 mg/l
Pressure loss:	Inspiration 0,5 l/s: 0,49 cmH2O (~hPa) 1,0 l/s: 1,47 cmH2O (~hPa) Expiration 0,5 l/s: 0,55 cmH2O (~hPa) 1,0 l/s: 1,67 cmH2O (~hPa)
Weight:	1,9 g
Length:	15,8 mm

X. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages (especially not for functional deficiencies, injuries, infections, and/or other complications or adverse events) caused by unauthorised product alterations, by repairs not performed by the manufacturer or by improper handling, care (cleaning/disinfection) and/or storage of the products in violation of the provisions of these instructions for use. This applies – to the extent permitted by law – both to damages to the products themselves caused thereby and to any consequential damages caused thereby.

Product specifications described herein are subject to change without notice.

HUMIDOTWIN® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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