




FAHL

HUMIDUAL®

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I. FOREWORD

These instructions for use are valid for HUMIDUAL® HME (Heat-Moisture Exchanger)

The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

Please read the instructions for use carefully before using the product for the first time!

II. INTENDED USE

HUMIDUAL® HMEs filter the breathing air, thereby reducing intrusion of particles into the patient's airways.

HUMIDUAL® HMEs reduce the formation of viscous secretions in the lungs. They store the moisture and warmth of the exhaled air in the filter medium and emit these to the inhaled air again during inspiration.

Application: spontaneous breathing, tracheotomised adult patients in inpatient or home care.

When used for the first time, the products must only be selected, used and inserted by appropriately trained physicians or medical professionals.

III. PRODUCT DESCRIPTION

HUMIDUAL® heat-moisture exchangers (HME) consist of a foam filter and a plastic housing

The plastic housing is equipped on the patient side with a central opening with 15 mm internal diameter, thereby allowing for connection to a 15 mm standard connector.

The overpressure valve in the lid of the housing prevents an excessive pressure build-up and allows secretions to be coughed up without having to remove the HME. It minimises the risk of the moistening medium sticking and allows larger quantities of secretion to be coughed up.

HUMIDUAL® is additionally equipped with an oxygen connecting piece at the side of the housing (5 mm plug-in connection), to which a connecting tube can be connected up, thereby allowing for combination with an oxygen source.

The product was packaged sterile and sterilised with ethylene oxide (EO).

IV. WARNING

HUMIDUAL® HMEs are single-patient products and are intended for single use only. They must not be cleaned or disinfected, as use of disinfectant can render the HME unusable.

Cleaning, disinfection, re-sterilisation and re-use can impair the safety and function of the product and are therefore not permissible!

The HUMIDUAL® HME must under no circumstances be inserted directly into the tracheostoma!

Care must be taken to ensure that the valve opening is not blocked by clothing or the like in order to enable unhindered breathing (risk of respiratory distress).

The HMEs must not be used in conjunction with a nebuliser or atomiser. Do not fill any water into the HME. This would invariably entail the risk of increased breathing resistance and thus obstruction of breathing.

V. PRECAUTIONS

If the HUMIDUAL® HME is clogged, or if signs of clogging appear, then the HME must immediately be removed and disposed of.

VI. CONTRAINDICATIONS

The HME should be removed if respiratory resistance is too high.

Use in patients with strong exsiccosis (dehydration), or excessive secretion in the lungs and airways, or strong air leakage flow (exhaled air not flowing through the HME), is not permissible.

Not for use in patients with impaired consciousness and/or in persons unable to look after themselves and requiring the help of others, who are not able to remove the HME by themselves in emergency situations.

The latter does not apply if the patient is in an institution under nursing supervision for 24 hours or in home care.

VII. INSTRUCTIONS

- Carefully examine the sterile packaging to ensure that it has not been tampered with or damaged. Do not use the product if the packaging has been damaged.
- Check the use-by or expiry date. Do not use the product after this date.
- Fit the HME with the central 15 mm opening onto the 15 mm standard connector of the tracheostomy tube by exerting light pressure ❶.
- Check the respiratory resistance.
- Remove the HUMIDUAL® HME after a wearing time of not more than 24 hours ❷, or if the HME is contaminated with secretions for instance.
To do this, hold the tracheostomy tube with two fingers and remove the HME by tilting it slightly downward. ❸

VIII. SERVICE LIFE

The maximum period of use / wearing period is 24 hours.

The HUMIDUAL® HME must be changed even several times a day if required, for instance if the respiratory resistance is in danger of increasing due to accumulation of secretions.

IX. TECHNICAL DATA (DIN EN ISO 9360-2)

Tidal volume (VT):	min. 250 ml	
Internal Volume:	16 cm ³	
Moisture lost:	9,7 md/dl	
Decrease of pressure:	Inspiration:	0,5 l/s: 0,30 cmH ₂ O (~hPa)
		1,0 l/s: 0,69 cmH ₂ O (~hPa)
		1,5 l/s: 1,12 cmH ₂ O (~hPa)
	Expiration:	0,5 l/s: 0,34 cmH ₂ O (~hPa)
		1,0 l/s: 0,78 cmH ₂ O (~hPa)
		1,5 l/s: 1,28 cmH ₂ O (~hPa)
Weight:	4,5 g	
Length:	33 mm	

X. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages (especially not for functional deficiencies, injuries, infections, and/or other complications or adverse events) caused by unauthorised product alterations, by repairs not performed by the manufacturer or by improper handling, care (cleaning/disinfection) and/or storage of the products in violation of the provisions of these instructions for use.

This applies – to the extent permitted by law – both to damages to the products themselves caused thereby and to any consequential damages caused thereby. Product specifications described herein are subject to change without notice.

Should a serious adverse event occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, this is to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

HUMIDUAL® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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