
# FAHL® Nebulizer-Set with aerosolmask

GEBRUCHSANWEISUNG • INSTRUCTIONS FOR USE

DE VERNEBLER-SET MIT AEROSOLMASKE  
EN NEBULIZER-SET WITH AEROSOLMASK

## BILDER / PICTURES

## BESTELLMINWEIS / ORDER INFO

### FAHL® Nebulizer-Set with aerosolmask

Größe/Size L

REF 51210-L

LOT 20200415

(01)04051948094206  
(10)20200415

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EN

## I. FOREWORD

These brief instructions for use are valid for the FAHL® Nebulizer-Set with aerosolmask. The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling. **Please read the instructions for use carefully before using the product for the first time!**

## II. INTENDED USE

Humidification of breathing air through inhalation of steam and medications and other inhalants approved for use in inhalation therapy.

## III. WARNINGS

Do not perform inhalation in the airways until you have been briefed by specialist staff.

Unless otherwise prescribed by the physician, inhalation should be performed 2 - 3 times daily for approx. 10 - 15 minutes. Once phlegm has loosened, it should be coughed up, or in the case of tracheostomy patients, suctioned using a suctioning device and a sterile suction catheter. Draughts and cold temperatures should be avoided immediately after inhalation. FAHL® Nebulizer-Set with aerosolmask is a single-patient product and is intended only for one single use.

## CAUTION!

Cleaning, disinfection, re-sterilisation and re-use can impair the safety and function of the product and are therefore not permissible!

## IV. PRECAUTIONS

Prior to any use, a visual and functional inspection of all system components must be performed.

A damaged product must not be used!

## V. PRODUCT DESCRIPTION

The FAHL® Nebulizer-Set with aerosolmask is a single-patient product and is ideally suited for aerosol-inhalation. The nebulizer-set can be used with inhalation devices powered by compressed air.

## VI. INSTRUCTIONS

### 1. CONNECTING THE NEBULIZER

Pull off the top section of the nebulizer by twisting both parts (anti-clockwise) gently against each other.

In the lower section of the nebulizer, pour in the exact quantity of the medication/ inhalant in accordance with the instructions provided by the manufacturer or the attending physician.

Observe the maximum capacity of 10 ml.

Place the top section of the nebulizer back onto the lower section. Connect the nebulizer by twisting both parts (clockwise) gently against each other. Ensure that these fit together securely.

### 2. CONNECTING THE AIR TUBE AND AEROSOLMASK

Check whether the air tube is bent or kinked.

Connect the air tube with the lower opening of the nebulizer.

Place the aerosolmask on the nebulizer head extension.

The connection can be established with the 22 ID or 22 AD / 15 ID connectors.

### 3. NEBULIZER

Hold the nebulizer in an upright position where possible.

Then connect the air tube with the air connection of the device.

## VII. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling. In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the products or by repairs, if these modifications or repairs were not carried out by the manufacturer or by an authorised specialist company.

This applies both to damages to the products itself caused thereby and to any consequential damages caused thereby.

Should a serious adverse event occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, this is to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

The manufacturer reserves the right to make changes to the product at any time without notice.

FAHL® is a trademark of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne, which is registered in Germany and the European Member States.