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EN I. FOREWORD

These instructions for use are valid for LARYVOX® Foam.

The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

Please read the instructions for use carefully before using the product for the first time!

II. INTENDED USE

LARYVOX® Foam is used to align the skin level of an irregular tracheostoma that is situated below the skin level.

In combination with LARYVOX® Tape, the contact surface of the tape is optimised, thus enabling a longer wearing time.

The speaking option is improved by exclusion of undesired losses of air.

III. WARNINGS

LARYVOX® Foam are single-patient products and are intended for single use only. It must not be cleaned or disinfected, as use of cleaning agent can render LARYVOX® Foam unusable.

Sterilisation and reuse can impair the function and is therefore not permissible either. LARYVOX® Foam must only and exclusively be used in combination with LARYVOX® Tape.

IV. COMPLICATIONS

The circular opening can, under certain circumstances, cause the usable diameter of the tracheostoma to be reduced, thereby leading to a reduction of the air supply.

Due to the specific adhesive properties of LARYVOX® Foam, reversible skin irritations can occur in individual cases.

V. CONTRAINDICATIONS

The use of LARYVOX® Foam should be avoided during radiotherapy, due to the risk of skin irritations.

After undergoing radiotherapy, only use LARYVOX® Foam after consulting your treating physician. Start using the tape at the earliest 2 weeks after completing radiotherapy.

INSTRUCTIONS

Clean the skin around the tracheostoma in front of a mirror using warm water and pH-neutral soap or OPTIFIAHL® stoma cleaning cloths. Dry the skin. Before cleaning, remove residues of adhesive from the skin if necessary, using OPTICLEAR®.

To protect the skin from irritation, apply OPTIGARD® to the affected area of the skin.

After removal of the protective foil on the back **●**, place LARYVOX® Foam concentrically over the tracheostoma **●**.

The different size versions (S / M / L) serve for individual adaptation in order to achieve an even skin level to be covered with LARYVOX® Tape **●**.

To accomplish this, sizes S / M / L can be used in combination.

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VIII. LEGAL NOTICES
The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages (especially not for functional deficiencies, injuries, infections, and/or other complications or adverse events) caused by unauthorised product alterations, by repairs not performed by the manufacturer or by improper handling, care (cleaning/disinfection) and/or storage of the products in violation of the provisions of these instructions for use.

This applies to the extent permitted by law to damages to the products themselves caused thereby and to any consequential damages caused thereby. Product specifications described herein are subject to change without notice. LARYVOX® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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