LARYVOX® HME FILTERKASSETTE

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EN

I. FOREWORD

These instructions for use are valid for LARYVOX® HME filter cassettes. The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling. **Please read the instructions for use carefully before using the product for the first time!**

II. INTENDED USE

LARYVOX® HME filter cassettes are particle filters consisting of a foam filter and a plastic housing. The LARYVOX® HME filter cassettes are affixed over the tracheostoma either on a special base plate with holding ring or on a tracheostomy tube with combi-adapter or 22 mm connector. LARYVOX® HME filter cassettes filter the breathing air, thereby reducing intrusion of particles into the patient's airways. The filter cassettes reduce the formation of viscous secretions in the lungs. They store the moisture and warmth of the exhaled air in the filter medium and emit these to the inhaled air again during inspiration.

When used for the first time the products must only be selected, used and inserted by clinicians or medical professionals trained in their use.

III. WARNINGS

LARYVOX® HME filter cassettes are single-patient products and are intended for single use only. They must not be cleaned or disinfected, as use of disinfectant can render the filter unusable. Sterilisation and reuse can impair the function and is therefore not permissible either. The LARYVOX® HME filter cassette must under no circumstances be inserted directly into the tracheostoma!

IV. PRODUCT DESCRIPTION

The LARYVOX® HME filter cassette is available in three different versions:

1. LARYVOX® HME filter cassette (REF 49800)

- Normal respiratory resistance
- For patients directly after laryngectomy
- For situations with normal physical exertion
- Plastic housing with outer diameter of 22 mm
- Compatible with all commonly used 22 mm receptacle systems (such as DURAVENT® KOMBI tracheostomy tube or LARYVOX® Tape base plate)

2. LARYVOX® HME HighFlow filter cassette (REF 49810)

- Low respiratory resistance
- For new users of HMEs
- For situations with strong physical exertion
- Plastic housing with outer diameter of 22 mm
- Compatible with all commonly used 22 mm receptacle systems (such as DURA-VENT® KOMBI tracheostomy tube or LARYVOX® Tape base plate)

3. LARYVOX® O₂ HME filter cassettes (REF 49802)

- Normal respiratory resistance
- With O₂ connection (the oxygen connecting piece fitted at the side of the housing enables a connecting tube to be connected, thereby allowing combination with an oxygen source).
- Ideal for postoperative application
- Plastic housing with an outer diameter of 22 mm and an inner diameter of 15 mm
- Compatible with all commonly used 22 mm receptacle systems (such as DURA-VENT® KOMBI tracheostomy tube or LARYVOX® Tape base plate) as well as with a 15 mm standard connector

V. COMPLICATIONS

If the particle filter is clogged or if signs of clogging appear, the LARYVOX® HME filter cassette must immediately be removed and disposed of.

VI. PRECAUTIONS

Care must be taken to ensure that the valve opening is not blocked by clothing or the like in order to enable unhindered breathing (risk of respiratory distress).

VII. CONTRAINDICATIONS

- The LARYVOX® HME should be removed if respiratory resistance is too high.
- Use in patients with strong exsiccosis (dehydration), or excessive secretion in the lungs and airways, or strong air leakage flow (exhaled air not flowing through the filter cassette) is not permissible.
- Not for use in patients with impaired consciousness and/or in non-independent persons requiring the help of others, who are not able to remove the LARYVOX® HME filter cassette by themselves in emergency situations.

VIII. INSTRUCTIONS

- Check whether the packaging and the product are in a faultless condition.
- Check the use-by or expiry date. Do not use after this date.
- Fit ❶ the LARYVOX® HME filter cassettes vertically onto the 22- or 15-mm connection of the tracheostomy tube or the base plate, applying gentle pressure.
- Check the respiratory resistance.
- Remove the LARYVOX® HME filter cassettes after a ❷ wearing period of 24 hours at most, or if the filter cassette is contaminated with secretions for instance. To do this, hold the base plate or the tracheostomy tube in position with two fingers ❸ and remove the filter cassette by tilting it downward slightly.

IX. SERVICE LIFE

- The maximum period of use / wearing period is 24 hours.
- The filter must be changed even several times a day if required, for instance if the respiratory resistance is in danger of increasing due to accumulation of secretions.

X. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH accepts no liability for damage caused by unauthorised product modifications, repairs not carried out by the manufacturer or by improper handling, care (cleaning/disinfection) and/or storage of the products contrary to the provisions of these instructions for use.

Should a serious adverse event occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, this is to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Product specifications described herein are subject to change without notice.

LARYVOX® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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