- The cannula must be entirely wet with the disinfection solution. Place only one tracheostomy tube in the sieve insert of the cannula-cleaning tub at one time.
- Once the exposure time for the disinfection has ended, remove the sieve insert with the tracheostomy tube and rinse the tracheostomy tube very thoroughly on the inside and the outside with lukewarm water, then dry the tracheostomy tube with a clean, soft, lint-free, cloth / non-woven cloth.
- No disinfectant residues must be present on the tube when it is inserted into the tracheostoma.
- Store the cleaned and disinfected tracheostomy tube in accordance with the tracheostomy tube manufacturer's instructions.

#### PRECAUTIONS

- **Always hold the sieve insert by the upper edge to avoid coming into contact with the disinfection solution.**
- **The tracheostomy tube must be completely submerged in the disinfection solution.**
- **After disinfection, the tracheostomy tube must be thoroughly rinsed with lukewarm water.**

#### III. DISPOSAL OF THE PRODUCT

##### CAUTION

The product must only be disposed of in accordance with the respective applicable national regulations for waste products. For further information, contact the manufacturer.

#### IV. LEGAL NOTES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages (especially not for functional deficiencies, injuries, infections, and/or other complications or adverse events) caused by unauthorised product alterations or by improper handling, care (cleaning/disinfection) and/or storage of the products in violation of the provisions of these instructions for use.

This applies – to the extent permitted by law – both to damages to the products themselves caused thereby and to any consequential damages caused hereby.

Should a serious adverse event occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, this is to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Product specifications described herein are subject to change without notice.

OPTICIT® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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## OPTICIT® DISINFECTANT FOR TRACHEOSTOMY TUBES

### I. INTENDED USE

OPTICIT® disinfectant (concentrate) is used for the disinfection of tracheostomy tubes. **Use only in diluted form as specified in the dosage directions.**

#### WARNING

**Please observe the manufacturer's directions regarding the disinfection and care of the tracheostomy tube.**

**Never disinfect the tracheostomy tube with cleaning agents and disinfectants not recommended by the manufacturer.**

**Failure to observe the tracheostomy tube manufacturer's directions for disinfection will result in an acute health hazard. This may also damage or destroy the tracheostomy tube.**

#### CAUTION

Contains glutaric dialdehyde.

#### CAUTION

Ensure that the warning notices and precautionary measures on the bottle label are always observed.

#### CAUTION

The date of expiry of the concentrate is indicated on the package label.

The solutions used for cleaning and disinfection must be freshly prepared prior to each use.

#### Concentration:

To prepare a disinfectant solution, dilute 8 ml of concentrate in water to 400 ml (this is equivalent to a 2% concentration).

#### Spectrum of activity & Exposure time

(at a concentration of 2%):

bactericidal (incl. MRSA), tuberculocidal, yeasticidal	15 Min.
virucidal (incl. HBV/HIV/HCV) pursuant to RKI recommendation	30 Min.

### II. DOSAGE DIRECTIONS AND APPLICATION

Remove all accessories that may possibly still be inserted before disinfecting the tracheostomy tube. First, rinse the tube thoroughly under running water.

We recommend using a suitable plastic tub with the sieve insert in which the tube can be placed together with the antiseptic (e.g. REF 31200, OPTIBRUSH® CONT cannula-cleaning tub with sieve insert) for disinfecting the tracheostomy tube with OPTICIT®. Then proceed as follows:

- Remove the sieve insert from the cannula-cleaning tub and fill 8 ml of OPTICIT® disinfection solution into the cannula-cleaning tub with the aid of the enclosed measuring container.
- Dilute the concentrate with water to 400 ml
- Place the pre-cleaned outer and/or inner cannula in the sieve insert of the cannula-cleaning tub and immerse this in the prepared disinfection solution in the cannula-cleaning tub. Hold the sieve insert by the upper edge in order to prevent coming into contact with and contaminating the solution.