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## EN PICTOGRAM LEGEND

<span></span>	<span></span>	<span></span>	<span></span>
Catalogue number	Manufacturer	<span></span>	<span></span>
Batch code	Date of manufacture	<span></span>	<span></span>
Consult instructions for use	Keep away from sunlight	<span></span>	<span></span>
Use by	Keep dry	<span></span>	<span></span>
Do not reuse	Sterilized using ethylene oxid	<span></span>	<span></span>
Content (in pieces)	<span></span>	<span></span>	<span></span>
Do not use if package is damaged	<span></span>	<span></span>	<span></span>

## I. FOREWORD

These instructions for use are valid for SENSFOAM® nursing compresses made of foam material. The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling. **Please read the instructions for use carefully before using the product for the first time!**

## II. INTENDED USE

Tracheal compresses serve to absorb tracheal secretions and to cushion the tracheal neck flange against the skin.

## III. WARNINGS

**Should you experience skin irritation (e.g. redness, itchiness), discontinue the product immediately and consult a physician. This may be due to special skin sensitivity or allergy-related incompatibility. Carefully examine the sterile packaging to ensure that it has not been tampered with or damaged. Do not use the product if the packaging has been damaged.**
**Check the use-by or expiry date. Do not use after this date.**

## IV. PRODUCT DESCRIPTION

Different versions of tracheal compresses are available. The compresses have a circular hole through which the cannula tube can be inserted and is also available with an additional slit opening, for fast compress changing.

## V. INSTRUCTIONS

Before replacing the compress, ensure that you have all required material available. The tracheal cannula should be cleaned and if required, disinfected before application.

The tracheal cannula is inserted through the existing, circular opening of the compress. The front of the tracheal compress is light beige and is indicated by a white skin-side of the compress is white. An additional slit next to the circular opening allows fast compress replacement without prior removal of the tracheal cannula from the tracheostoma.

## CAUTION!

When inserting the tracheal cannula, please observe the application instructions of its manufacturer.

## VI. SERVICE LIFE

The compress supplied sterile and is intended for single use.

For hygienic reasons, the compress must be replaced at least every 12 hours or several times daily during severe secretion.

## VII. LEGAL NOTICES

The manufacturer Andreas FahI Medizintechnik-Vertrieb GmbH will not accept any liability for damages (especially not for product deficiencies, injuries, infections, and/or other complications or adverse events) caused by unauthorised product alterations, by repairs not performed by the manufacturer or by manufacturer or by patient/user in violation of the provisions of the products in violation of the provisions of these instructions for use. This applies – to the extent permitted by law – both to damages to the products themselves caused thereby and to any consequential damages caused thereby.

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## VII. LEGAL NOTICES

The manufacturer Andreas FahI Medizintechnik-Vertrieb GmbH will not accept any liability for damages (especially not for product deficiencies, injuries, infections, and/or other complications or adverse events) caused by unauthorised product alterations, by repairs not performed by the manufacturer or by manufacturer or by patient/user in violation of the provisions of the products in violation of the provisions of these instructions for use. This

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