

# TRACHEOSILC® VARIO CUFF-SUCTION



- ▶ Tracheostomy tube with CUFF and integrated suction line (SUCTION) to remove secretions above the CUFF
- ▶ Also available with CUFF without integrated suction line (SUCTION) (REF 18830) and without CUFF (REF 18820 on separate data sheet)
- ▶ Soft, flexible neck-flange with key tube identification data clearly marked
- ▶ Clear and transparent tube shaft for easy handling during cleaning
- ▶ With rotating 15 mm connector (VARIO) to buffer torque and tension caused by tractive forces from ventilation equipment
- ▶ Rounded cannula tip to prevent irritation of the mucosa in the trachea
- ▶ Set includes tube holder, obturator and product ID card

## Order information:

**TRACHEOSILC® VARIO CUFF-SUCTION**

SIZE 5-9.5



**REF 18840**

When ordering, please state the size after the item number!

Size	O.D. neck flange mm	O.D. Tip mm	Length mm	REF
5.0	7.3	7.3	57	<b>18840-05</b>
6.0	8.7	8.7	63	<b>18840-06</b>
7.0	10.0	10.0	71	<b>18840-07</b>
7.5	10.7	10.7	73	<b>18840-075</b>
8.0	11.0	11.0	75	<b>18840-08</b>
8.5	11.7	11.7	78	<b>18840-085</b>
9.0	12.3	12.3	80	<b>18840-09</b>
9.5	13.3	13.3	83	<b>18840-095</b>

O.D. neck flange = outer diameter behind the neck flange · O.D. Tip = outer diameter at the tube tip · Length = length over the outer curve