

**OPTICLEAR®:**


- ▶ Adhesive remover for laryngectomized patients when using base plates.  
Removes adhesives of all types
- ▶ Non-woven cleaning wipe, soaked in an alcohol-based cleaning solution.  
Adhesive residues are quickly and gently removed from the skin and absorbed by the wipe
- ▶ Large, touch-proof
- ▶ Supplied in yellow box to make it easily discernable from its complementary product  
OPTIGARD® Protective Dressing (see below)



CE

**OPTIGARD®:**

- ▶ For preparative skin care for laryngectomized patients when using base plates or skin adhesives of all types
- ▶ Helps to protect against skin irritations which can be caused by skin adhesives
- ▶ Increases the adhesive power of the skin adhesive
- ▶ Easy application on cleaned skin: The skin is covered with a thin protective film which dries after a few seconds.
- ▶ When changing a base plate the film does not necessarily have to be removed but it can be easily washed off with cleaning lotion and water
- ▶ Supplied in a green box to make it easily discernable from its complementary product:  
OPTICLEAR® Adhesive Remover (see above)



CE

**Order Information:**

OPTICLEAR®, Pack of 30

REF 33500

OPTIGARD®, Pack of 30

REF 33600