




# LARYVOX® TAPE HYDROSOFT

- ▶ Particularly skin-friendly
- ▶ Non-invasive
- ▶ High fixation comfort
- ▶ Ability to enclose fluids/secretions
- ▶ Not suitable for use in combination with shower guard!
- ▶ No friction inside the tracheostoma
- ▶ With 22 mm housing to be combined with any HME and speaking valve with 22 mm adapter (KOMBI)
- ▶ Can be combined with special silicone tubes (LARYNGOTEC® KOMBI CLIP or LARYNGOTEC® KOMBI CLIP LINGO) if additional tracheal stabilisation is needed
- ▶ Available in 3 different shapes: round, oval, XL-oval



 Consult instructions for use  
 Do not reuse  
 CE-Kennzeichnung  
**Product made in Germany**

## Order information:

**LARYVOX® Tape HYDROSOFT round**, pack of 15

**REF 48130**

**LARYVOX® Tape HYDROSOFT oval**, pack of 15

**REF 48230**

**LARYVOX® Tape HYDROSOFT XL oval**, pack of 15

**REF 48330**