






- ▶ Adhesive foam filter for laryngectomized patients to provide cover and protection for the tracheostoma
- ▶ LARYNGOFIX® further humidifies, filters and warms the inhaled air
- ▶ Its rounded shape assures an excellent fit on the tracheostoma for proper cover and protection from dust particles or insects
- ▶ Simple to apply by adhesive stripe on the top side of the tracheostoma; adhesive strength may further be enhanced by use of the OPTIGARD® protective dressing wipe (cf separate data sheet)
- ▶ Comes in different sizes (large L and small S) and different colours (beige or white)
- ▶ **Special models:**
 - LARYNGOFIX® PLUS - with two adhesive stripes, one at the top side, one at the bottom side, to assure a tight fit on the tracheostoma
 - LARYNGOFIX® PRO - material thickness 5 mm
 - LARYNGOFIX® TOUCH L - with push button to facilitate finger occlusion of the tracheostoma when speaking with a voice prosthesis



-  Keep away from sunlight
-  Keep dry
-  Do not reuse
-  Consult instructions for use
-  CE marking

Order Information:

LARYNGOFIX® S, beige, 4 x 5.5 cm, Pack of 10	REF 45900
LARYNGOFIX® L, beige, 4.8 x 6.7 cm, Pack of 10	REF 45950
LARYNGOFIX® S, white, 4 x 5.5 cm, Pack of 10	REF 45905
LARYNGOFIX® L, white, 4.8 x 6.7 cm, Pack of 10	REF 45955
LARYNGOFIX® PLUS S, beige, 4 x 5.5 cm, Pack of 10	REF 45902
LARYNGOFIX® PLUS L, beige, 4.8 x 6.7 cm, Pack of 10	REF 45952
LARYNGOFIX® PRO S, beige, 4 x 5.5 cm - 5 mm, Pack of 10	REF 45901
LARYNGOFIX® PRO L, beige, 4.8 x 6.7 cm - 5 mm, Pack of 10	REF 45956
LARYNGOFIX® TOUCH L, beige, 4.8 x 6.7 cm, Pack of 10	REF 45911