






# TRACHEOTEX® BIB

- ▶ To cover and protect the tracheostoma of laryngectomized and tracheostomized patients
- ▶ Available in multiple layers (3-layers, 5L, 8L, 12L) - enables optimal coordination with each seasonal climate
- ▶ The fine mesh fabric (Diolen tulle) assures a high degree of air permeability and at the same time ensures that the inspired air is warmed, humidified and filtered
- ▶ Washable and reusable
- ▶ With adherent fastener
- ▶ The neck width can be individually adjusted to a size between 34 cm and 46 cm in the adult variants.
- ▶ The fastener is located at the side of the bib and is therefore easily accessible.
- ▶ Adult and pediatric variants (PED / 24 - 33 cm)
- ▶ Many colour variants  
(cf colour sample card – available as extra brochure)



-  Keep away from sunlight
-  Keep dry
-  Single patient use
-  Consult instructions for use
-  CE marking
-  Hand wash
-  Do not dry in the drum dryer

## Order information:

<b>TRACHEOTEX® BIB 3L, 34 - 46 cm</b>	<b>REF 40100</b>
<b>TRACHEOTEX® BIB 5L, 34 - 46 cm</b>	<b>REF 40200</b>
<b>TRACHEOTEX® BIB 8L, 34 - 46 cm</b>	<b>REF 40300</b>
<b>TRACHEOTEX® BIB 12L, 34 - 46 cm</b>	<b>REF 40400</b>
<b>TRACHEOTEX® BIB 3L PED, 24 - 33 cm</b>	<b>REF 40110</b>
<b>TRACHEOTEX® BIB 5L PED, 24 - 33 cm</b>	<b>REF 40210</b>

Please specify color number on the color sample card when ordering