





# TRACHEOTEX® EASY

- ▶ To cover and protect the tracheostoma of laryngectomized and tracheostomized patients
- ▶ Cut in the shape of a fashionable, triangular neckerchief
- ▶ Easily fastened by knotting the tips at the back of the neck
- ▶ An insert, made of multiple layers of fine mesh fabric warms, filters and humidifies the inspired air. At the same time escaping tracheal secretions are reliably absorbed
- ▶ Washable at 30°C and reusable
- ▶ Available in many different colors and patterns matching the wardrobe (cf color sample card – available as extra brochure)



-  Keep away from sunlight
-  Keep dry
-  Single patient use
-  Consult instructions for use
-  CE marking
-  Hand wash
-  Do not dry in the drum dryer

## Order information:

**TRACHEOTEX® EASY**, max. 55 cm

Please specify color number on the color sample card when ordering

**REF 42002**