

LARYNGOTEC® PRO KOMBI / LARYNGOTEC® PRO KOMBI LINGO



- ▶ Tracheostomy tube/silicone tube for long-term use, can be used up to 3-6 months
- ▶ The proximal cone optimally seals the tracheostoma in combination with a cannula support strap
- ▶ Medical grade silicone: soft, pliable, lightweight
- ▶ Transparent, inconspicuous design
- ▶ The ergonomically designed cannula shield offers a pleasing wearing comfort
- ▶ With 22 mm adapter (KOMBI)
- ▶ The cannula is available in the sizes: 8,9,10 and 12 and in two lengths (36 mm / 55 mm)
- ▶ The size designation is printed on the distal cone of the 22 mm KOMBI connector. This means that it is not necessary to remove a filter cassette from the cannula to identify the size used
- ▶ With (LARYNGOTEC® PRO KOMBI LINGO, REF 14722) or without (LARYNGOTEC® PRO KOMBI, REF 14922) sieve fenestration
- ▶ Can be tailor-made to customer requirements (length, sieve, window) in consultation with our sales representatives
- ▶ Including broad adjustable tube holder and a product ID card

Order information:

LARYNGOTEC® PRO KOMBI, Size 8-10, 12



REF 14922

LARYNGOTEC® PRO KOMBI LINGO LINGO, Size 8-10, 12



REF 14722

When ordering, please state the size after the item number!

Size	I.D. tip (mm)	O.D. tip (mm)	Length over middle curve (mm)	Length over outer curve (mm)	REF KOMBI	REF KOMBI LINGO
8	9,5	12,0	36,0	42,0	14922-0836	14722-0836
8	9,5	12,0	55,0	65,5	14922-0855	14722-0855
9	10,5	13,5	36,0	42,0	14922-0936	14722-0936
9	10,5	13,5	55,0	65,5	14922-0955	14722-0955
10	12,0	15,0	36,0	42,0	14922-1036	14722-1036
10	12,0	15,0	55,0	65,5	14922-1055	14722-1055
12	13,5	17,0	36,0	42,0	14922-1236	14722-1236
12	13,5	17,0	55,0	65,5	14922-1255	14722-1255

O.D. Tip = Outer diameter of the outer tube at the tip of the cannula · I.D. = Inner diameter at the tip of the cannula
Length = Length over the outer curve · Length of mean radius = Length on midline, from flange to tip