



NECKFIX®

REF 32620



REF 32610



REF 32612



REF 32622



- ▶ Tracheostomy tube holder
 - “K” variant: with adherent fastener - white
 - “Blue K” variant: with adherent fastener - blue
 - “Blue PED K” variant: with adherent fastener, for pediatric use – blue
 - “Beige K PED” variante: with adherent fastener, for pediatric use – beige
 - “Beige K” variante: with adherent fastener – beige
 - “Beige K XL” variante: with adherent fastener – beige
- ▶ Assures the secure fixation of tracheostomy tubes
- ▶ Foam core surrounded by soft non-woven fabric
- ▶ Economic tube holder for routine use
- ▶ Slim design, one part
- ▶ Individually packed

SCOPE OF DELIVERY

▶ Instruction for use	
NECKFIX® K	REF 32620
1 x Tracheostomy Tube Holder	
oder	
NECKFIX® Blue K	REF 3261x
1 x Tracheostomy Tube Holder	
oder	
NECKFIX® Beige K	REF 3262x
1 x Tracheostomy Tube Holder	

Product	Length	Width	REF
NECKFIX® K	45 cm	3,5 cm	32620
NECKFIX® Blue K	45 cm	3,5 cm	32610
NECKFIX® Blue PED K	34 cm		32612
NECKFIX® Beige K PED	34 cm		32621
NECKFIX® Beige K	46 cm		32622
NECKFIX® Beige K XL	50 cm		32623

