


# DURACUFF® VARIO

- ▶ Tracheostomy tube with a rotating 15 mm swivel connector (VARIO) to buffer external forces like e.g tension and torque forces exerted by the ventilation tubing
- ▶ With low-pressure cuff
- ▶ Soft, pliable, lightweight
- ▶ Transparent, inconspicuous design
- ▶ Set includes broad and adjustable tube holder (OPTIFLAUSCH®), obturator and product ID card
- ▶ Three lengths: Standard length, medium (M) and extra-long (XL)



- Sterilized using ethylene oxide
- Do not resterilize
- Keep away from sunlight
- Single patient use
- Consult instructions for use

**Rx only**  
**CAUTION:**  
 United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner (in the US only).

**CE 0482** CE marking identification number of the notified body

**Product made in Germany**

## Order information:

**DURACUFF® VARIO, SIZE 7-12**

**DURACUFF® VARIO M, SIZE 7-12**

**DURACUFF® VARIO XL, SIZE 7-12**

When ordering, please state the size after the item number!

- REF 13000**
- REF 13100**
- REF 13200**

Size	O.D. neck flange mm	O.D. Tip mm	Length mm standard	Length mm medium (M)	Length mm extra long (XL)
7	10.2	8.8	65	75	100
8	11.2	9.8	70	80	105
9	12.4	11.2	75	85	105
10	13.4	12.2	85	95	110
11	14.4	13.2	87	97	115
12	15.4	14.2	90	100	120

O.D. Tip = Outer diameter of the outer tube at the tip of the cannula · O.D. neck flange = Outer diameter behind the neck flange · I.D. = Inner diameter at the tip of the cannula · Length = Length over the outer curve