

SENSOFOAM® PAD PROFESSIONAL

TRACHEAL COMPRESSES

- ▶ Made of a high-quality 3-layer foam material to prevent adhesion to clotted blood or dried secretions and to absorb secretions
- ▶ Transmissible to air and humidity to prevent skin maceration and inflammation
- ▶ With two punched holes on each side of the central hole to allow fixation – together with the neck-flange of the tracheostomy tube – to a tube holder
- ▶ SENSOFOAM® PAD Professional with the shape of the flange of a tracheostomy tube to cushion the pressure of the neck-flange to reduce pain and the risk of pressure sores
- ▶ With slit to allow easy exchange of the compress while the tracheostomy tube remains in situ



STERILE Sterile
Do not reuse
Consult instructions for use
CE CE marking

Order information:

SENSOFOAM® PAD Professional, 7.8 x 3.8 cm, pack of 10

REF 30864