

SENSOTRACH® ALU

Tracheal compress

- ▶ Tracheal compresses with aluminium coating to reduce adhesion of the compress to skin around the tracheostoma
- ▶ **SLIT** variant with slit for easy exchange while the tracheostomy tube is in situ



Do not reuse



Consult instructions for use



CE marking

Order Information:

SENSOTRACH® ALU, 8 x 10 cm, pack of 10

SENSOTRACH® ALU SLIT, 8 x 10 cm, pack of 10

SENSOTRACH® ALU SLIT, 8 x 9 cm, pack of 10

REF 30000

REF 30030

REF 30036