

LARYVOX® SOS MASK

- ▶ Ventilation Funnel for emergency use
- ▶ Allows safe mouth-to-tracheostoma ventilation in cases of emergency
- ▶ Consists of two parts: a soft silicone tracheostoma mask and an adapter tube which is attached to the mask
- ▶ Suitable for direct ventilation through the tracheostoma as well as for an indirect ventilation through a tracheostomy tube
- ▶ Easy to clean, lightweight, handy, small
- ▶ The easy to detach plug-in connection enables the ventilation funnel to be easily disassembled into its individual parts which can be easily stored, e.g. in a jacket pocket
- ▶ Additionally included in scope of supply: Emergency ID card and Emergency sticker



Single patient use



Consult instructions for use



CE marking

Order information:

LARYVOX® SOS MASK

REF 75010