


# LARYVOX® TAPE STANDARD

Base Plate

- ▶ Transparent appearance with water-repellent effect
- ▶ Non-invasive
- ▶ Good adhesive strength on the skin around the tracheostoma
- ▶ Suitable in situations with unproblematic tracheostoma anatomy
- ▶ With 22 mm housing to be combined with any HME and speaking valve with 22 mm adapter (KOMBI)
- ▶ Can be combined with special silicone tubes (LARYNGOTEC® KOMBI CLIP or LARYNGOTEC® KOMBI CLIP LINGO) if additional tracheal stabilisation is needed
- ▶ Available in 3 different shapes: round, oval, XL-oval



 Consult instructions for use

 CE marking

 CE-Kennzeichnung

Product made in Germany

## Order information:

**LARYVOX® Tape STANDARD round**

Pack of 15

**REF 48100**

**LARYVOX® Tape STANDARD oval**

Pack of 15

**REF 48200**

**LARYVOX® Tape STANDARD XL oval**

Pack of 15

**REF 48300**