

LARYVOX® ALARM

- ▶ Device to generate a loud alarm, through an acoustic and visual signal, in cases of emergency when calling for help is not possible
- ▶ The acoustic alarm is performed as a continuous signal by pulling out the contact pin with neck strap. This makes it possible to quickly alert other people
- ▶ A label reading: “...Attention: I am a laryngectomee and cannot breathe through nose and mouth...” is printed on the device to explain the emergency situation to bystanders or paramedics.
- ▶ This label is available in several languages and can be supplied in any desired language
- ▶ Battery operated



Consult instructions for use

Order information:

LARYVOX® ALARM

REF 76000