

# DURAPUR® IC

Inner Cannulas non-fenestrated



▶ Three types of non-fenestrated inner cannulas for DURAPUR® tracheostomy tubes are available:

1. With flat profile (IC) if no connection to any peripheral device (ventilation equipment, or HME) is desired (REF 23801 and REF 23701 XL)
2. With fixed 15 mm-UNI-connector (ICU) for attachment of an HME (REF 23811)
3. With 22 mm-KOMBI-adapter (ICK) to connect HMEs with 22 mm-adapter (REF 23881)

▶ Size of the inner cannulas to be selected according to the size of the DURAPUR® outer cannula



Keep away from sunlight



Single patient use



Consult instructions for use



CE marking identification number of the notified body

## Order information:

<b>DURAPUR® IC</b>	SIZE 7-12	 MRT	<b>REF 23801</b>
<b>DURAPUR® IC XL</b>	SIZE 8-10	 MRT XL	<b>REF 23701</b>
<b>DURAPUR® ICU</b>	SIZE 7-12	 15mm MRT	<b>REF 23811</b>
<b>DURAPUR® ICK</b>	SIZE 7-12	 22mm MRT	<b>REF 23881</b>

When ordering, please state the size after the item number!