

HUMIDOLIFE®

REF 46445



- ▶ Heat and Moisture Exchanger / artificial nose (HME)
- ▶ 15 mm connector (UNI) - fits all tracheostomy tubes with 15 mm connector
- ▶ Cellulose moisture exchanger that ensures efficient humidification of the respiratory tract for up to 24 hours
- ▶ An oxygen connecting piece fitted to the side of the housing (6 mm ISO plug-in connection) enables a connecting tube to be connected up, thereby allowing combination with an oxygen source
- ▶ Small, interior chamber for the collection of sputum.
- ▶ Free of latex and DEHP

SCOPE OF DELIVERY**HUMIDOLIFE**

1 x HME filter cassette

REF 46445**Technical data**

Tidal volume	150 - 1500 ml		
Weight	8.5 g	0.5 l/s:	0.55 Hpa
Length	45 mm	1.0 l/s:	1.49 Hpa
Dead space	16 ml	1.5 l/s:	2.83 Hpa
Moisture loss			
VT = 250 ml:	11.0 mgH ₂ O/L	Expiration	
VT = 500 ml:	15.4 mgH ₂ O/L	0.5 l/s:	0.42 Hpa
VT = 750 ml:	16.6 mgH ₂ O/L	1.0 l/s:	1.18 Hpa
VT = 1000 ml:	18.7 mgH ₂ O/L	1.5 l/s:	2.12 Hpa

Product	Unit	REF
HUMIDOLIFE®	50 pieces, packed singly	46445