

LARYVOX® SOS SAFE-MASK

- ▶ Emergency Ventilation Mask for laryngectomees
- ▶ Allows safe mouth-to-tracheostoma ventilation in cases of emergency
- ▶ A disposable valve protects the user from possible risk of infection: The exhaled air, as well as any secretions escape through the valve at the side of the mask and are not inhaled by the first-aider
- ▶ The funnel-shaped design and the flexible silicone material provide a reliable seal of the tracheostoma, even on an uneven skin surface
- ▶ Additionally included in scope of supply: Emergency ID card and Emergency sticker



Single patient use



Consult instructions for use



CE marking

Order information:

LARYVOX® SOS SAFE-MASK

REF 75000