






LARYVOX® PAD CONNECT

Sealing plate



- ▶ LARYVOX® PAD CONNECT is a self-adhesive sealing plate made of soft, skin-coloured medical silicone with an integrated 22 mm connector (KOMBI)
- ▶ Particularly skin-friendly due to the material used, so the skin around the tracheostoma is less irritated
- ▶ Used to align the skin level of an irregular tracheostoma that is situated below the skin level
- ▶ Due to the integrated 22 mm connector (KOMBI) a suitable HME filter/artificial nose can be inserted directly into the PAD without any problems and without a base plate
- ▶ Due to its material properties and ergonomic shape, the LARYVOX® PAD CONNECT adapts to the skin around the tracheostoma
- ▶ The plano-convex shape of the LARYVOX® PAD CONNECT ensures high quality sealing that makes voice expression easier
- ▶ In addition to reduced leakage of secretion, the lateral flow of breathing air can also be minimised
- ▶ An adhesive for fixing is not necessary, but if an increased adhesive strength is required (e.g. at high tempering pressure), FAHL® Silicone Glue (REF 29810) can be used additionally



Order information:

LARYVOX® PAD CONNECT

pack of 2

REF 48710