

HUMIDOBAC®

REF 46832



REF 46831



REF 46830



- ▶ The following variants are available:
 - Humidobac® HME, REF 46832
 - Humidobac®, REF 46831
 - Humidobac® HME (sterile), REF 46830
- ▶ HUMIDOBAC® is for the filtration and retention of viruses and bacteria
- ▶ To be used with artificial respiratory equipment
- ▶ HUMIDOBAC® HMEs additionally humidify the breathing air
- ▶ Bacterial/viral filtration: > 99,999%
- ▶ Filter should be placed between ETT/Tracheal canula and Y-piece of breathing circuit
- ▶ Housing transparent and non-conductive
- ▶ Very lightweight (< 35 g)
- ▶ Should be replaced every 24 hours or more

SCOPE OF DELIVERY

- ▶ Manual

HUMIDOBAC® HME **REF 46832**

1 x Bacterial/viral filter

oder

HUMIDOBAC® **REF 46831**

1 x Bacterial/viral filter

oder

HUMIDOBAC® HME **REF 46830**

1 x Bacterial/viral filter (sterile)

Please specify quantity when ordering.

Product	Unit	REF
HUMIDOBAC® HME	1 or 50 pieces, individually packed	46832
HUMIDOBAC®	1 or 50 pieces, individually packed	46831
HUMIDOBAC® HME	1 or 50 pieces, individually packed, (sterile)	46830

CE 0482