

DURATWIX® UNI-VARIO LINGO

1 ICV + 1 ICFU

- ▶ Tracheostomy tube with two inner cannulas, both with 15 mm connector, one with rotating 15 mm swivel connector (ICV), the other one fenestrated (ICFU)
- ▶ Soft double-swivel ROTATWIX® neckflange, rotating in 360° for enhanced neck mobility
- ▶ With ASPILO® sieve-fenestration - roof-tile shaped (horizontal)
- ▶ Medical-grade polyurethane, soft, atraumatic, thermosensible
- ▶ Modular: 6 different types of inner cannulas
- ▶ Set includes broad and adjustable OPTIFLAUSCH® tube holder, obturator and product ID card
- ▶ Size 7-10 available in two lengths
- ▶ Colour coding
Inner cannula: fenestrated/unfenestrated, standard length/short
Imprint outer cannula neckflange: sieve-fenestrated/non-fenestrated



SHORT

STERILE EO Sterilized using ethylene oxid
 Do not resterilize
 Keep away from sunlight
 Single patient use
 Consult instructions for use
Rx only CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner (in the US only).
 CE 0482 CE marking identification number of the notified body
Product made in Germany

Order information:

DURATWIX® UNI-VARIO LINGO 1 ICV + 1 ICFU, SIZE 6-10



REF 19261

DURATWIX® UNI-VARIO LINGO short 1 ICV + 1 ICFU, SIZE 7-10



REF 19262

When ordering, please state the size after the item number!

Size	O.D. neck flange mm	O.D. Tip mm	I.D. mm	Length mm standard	REF standard	Length mm short	REF short
6	11.9	10.3	6.0	82.0	19261-06	-	-
7	11.9	10.8	7.0	82.0	19261-07	72.0	19262-07
8	12.5	11.4	8.0	85.0	19261-08	75.0	19262-08
9	13.7	12.5	9.0	88.0	19261-09	78.0	19262-09
10	15.0	13.8	10.0	91.0	19261-10	81.0	19262-10

O.D. Tip = Outer diameter of the outer tube at the tip of the cannula · O.D. neck flange = Outer diameter behind the neck flange · I.D. = Inner diameter at the tip of the cannula · Length = Length over the outer curve