

HUMIDOPHONE® / HUMIDOPHONE® TYP 15

REF 46480



REF 46485



- ▶ Speaking Valve with integrated filter
- ▶ Low profile, lightweight
- ▶ Open-pore polyurethane filter
- ▶ Filter in proximal position (towards the patient)
- ▶ **HUMIDOPHONE®** (REF 46480) with 22 mm (KOMBI) combi-adapter and **HUMIDOPHONE® TYP 15** (REF 46485) with 15 mm (UNI) connector

SCOPE OF DELIVERY

- ▶ Instruction for use

HUMIDOPHONE® **REF 46480**
30 x Speaking Valve

or

HUMIDOPHONE® TYP 15 **REF 46485**
30 x Speaking Valve

The HUMIDOPHONE® speaking valve must be only used by patients with a partially or fully preserved larynx!

Product	Unit	REF
HUMIDOPHONE®	30 pieces, packed by 5 pcs.	46480
HUMIDOPHONE® TYP 15	30 pieces, packed by 5 pcs.	46485

