



TRACHLINE® HUMIDOSTOM® MINI

Filter cassette

- ▶ HME filter (Heat-Moisture Exchanger)
- ▶ Open pore polyurethane foam filter
- ▶ The plastic housing is equipped on the patient side with a bayonet closure and thereby allows connection to a TRACHLINE® tracheostomy tube
- ▶ With oxygen supply port
- ▶ Individually packed in convenient peel-packs



Order information:

TRACHLINE® HUMIDOSTOM® MINI

Pack of 30

REF 46860