



BLOM-SINGER® CLASSIC™ INDWELLING LEF / LF (LARGE ESOPHAGEAL FLANGE / LARGE FLANGE) VOICE PROSTHESIS

SPECIAL ORDER

- ▶ The Blom-Singer® Special Order Series offers a large assortment of special indwelling prostheses that cover specific requirements
 - ▶ All variants are made of medical grade silicone and are furnished with specific product features to suit the respective indications
 - ▶ Blom-Singer® CLASSIC™ Indwelling LEF (Large Esophageal Flange) is furnished with an enlarged oesophageal flange (19 mm diameter) for treatment of periprosthetic leakages as appropriate for the indication
 - ▶ Optionally, this version is also available as a „Large Flange“ version
 - ▶ In the Blom-Singer® CLASSIC™ Indwelling LF (Esophageal + Tracheal Flange) „Large Flange“ variant, both the oesophageal flange and the tracheal flange are enlarged
 - ▶ Both voice prosthesis are available as 16 and 20 Fr. version in sizes 4 - 14
 - ▶ With the help of the LARYVOX® INSERTER the voice prosthesis can be placed anterograde quickly and easily:
 LARYVOX® INSERTER (REF 25900) - for different types of voice prosthesis
 LARYVOX® INSERTER EXPERT (REF 25910) - for 16 Fr. voice prosthesis (or DUAL VALVE™ 20 Fr.)
 LARYVOX® INSERTER EXPERT SMALL (REF 25915) - for 4 mm voice prosthesis
- Application with gel capsule technique:
 LARYVOX® INSERTER SPECIAL - (REF 25920) - for difficult tracheostoma



Single patient use
 Consult instructions for use

Order information:

Blom-Singer® CLASSIC™ Indwelling LEF (Large Esophageal Flange)	16 Fr.	SIZE 4 -14	REF IN16XX-LEF
	20 Fr.	SIZE 4 -14	REF IN20XX-LEF
Blom-Singer® CLASSIC™ Indwelling LF (Large Esophageal + Tracheal Flange)	16 Fr.	SIZE 4 -14	REF IN16XX-LF
	20 Fr.	SIZE 4 -14	REF IN20XX-LF
LARYVOX® INSERTER			REF 25900
LARYVOX® INSERTER EXPERT			REF 25910
LARYVOX® INSERTER EXPERT SMALL			REF 25915
LARYVOX® INSERTER SPECIAL			REF 25920