

LARYVOX® TAPE HYPOALLERGEN

BASE PLATE

- ▶ Particularly for sensitive skin
- ▶ Non-invasive
- ▶ May reduce allergic skin reactions
- ▶ Adapts comfortably to movements
- ▶ Very good adhesive properties
- ▶ High flexibility and long wearing period
- ▶ No friction inside the tracheostoma
- ▶ With 22mm housing to be combined with any HME and speaking valve with 22mm adapter
- ▶ Can be combined with special silicone tubes (LARYNGOTEC® KOMBI CLIP or LARYNGOTEC® KOMBI CLIP LINGO) if additional tracheal stabilisation is needed
- ▶ Available in 3 different shapes: round, oval, XL-oval




ROUND




OVAL



XL OVAL

 Consult instructions for use

 CE marking

 CE-Kennzeichnung

Product made in Germany

Order information:

LARYVOX® Tape HYPOALLERGEN round, pack of 15

REF 48050

LARYVOX® Tape HYPOALLERGEN oval, pack of 15

REF 48150

LARYVOX® Tape HYPOALLERGEN XL oval, pack of 15

REF 48250