

MUCOPROTECT®

Cuff Pressure Gauge

- ▶ For cuff pressure monitoring of tracheostomy tubes and endotracheal tubes with cuff
- ▶ Large surface allows precise adjustment to ensure maximum control accuracy
- ▶ Optimal pressure range (20-30 cm H₂O) highlighted in blue colour
- ▶ Lateral pressure release button
- ▶ Can be used for all control balloons with Luer-Lock-connector (market standard)
- ▶ Connection tube included
- ▶ Connection tube also separately available (REF 19506)



Order information:

MUCOPROTECT® Cuff pressure gauge

REF 19500

MUCOPROTECT® Connection tube for cuff pressure gauge

REF 19506