








LARYVOX® FOAM

Size S, M, L

- ▶ Ring-shaped foam pad with a central opening
- ▶ Is used to align the skin level of an irregular tracheostoma that is situated below the skin level.
- ▶ With adhesive contact surface
- ▶ In combination with LARYVOX® Tape, the contact surface of the tape is optimised, thus allowing a longer wearing time
- ▶ Further to this the speaking option is improved by the exclusion of undesired loss of air.
- ▶ The different size versions (S / M / L) serve for individual adaptation in order to achieve an even skin level to be covered with LARYVOX® Tape.



-  Keep away from sunlight
-  Keep dry
-  Do not reuse
-  Consult instructions for use
-  CE marking

Order Information:

LARYVOX® Foam S, Pack of 10

REF 48420

LARYVOX® Foam M, Pack of 10

REF 48430

LARYVOX® Foam L, Pack of 10

REF 48440