


DURATWIX® JUNIOR PED / PDC

- ▶ Tracheostomy tube for children (paediatrics), available without CUFF (REF 17436) and with CUFF (REF 17438)
- ▶ Adaptive, supportive and lightweight material
- ▶ Free of latex and DEHP
- ▶ Soft, flexible neck flange with imprint of the essential specifications of the cannula
- ▶ Rounded tube tip to prevent irritation of the mucous membranes in the trachea
- ▶ With X-ray contrast strip
- ▶ Set includes a broad, stable tube holder, an obturator and cannula passport



Order information:

DURATWIX® JUNIOR PED

Size 2,5-6



REF 17436

DURATWIX® JUNIOR PDC

Size 2,5-6



REF 17438

When ordering, please state the size after the item number!

Size	O.D. mm	I.D. mm	Angle ca.	Cuff (Ø) mm	Length over outer curve mm	Length over middle curve mm
2,5	4,0	2,5	120°	8	34	32
3	4,5	3,0	120°	9	37	35
3,5	5,1	3,5	120°	10	42	39
4	5,8	4,0	120°	12	46	43
4,5	6,5	4,5	120°	12	50	47
5	7,1	5,0	120°	15	53	49
5,5	7,6	5,5	120°	15	58	54
6	8,3	6,0	120°	18	64	60

O.D. = Outer diameter · I.D. = Inner diameter