



SENSOTRACH® PED

Tracheal compress for Pediatric Care

- ▶ High-quality tracheal compresses with anatomical rounded shape for maximum wearing comfort with a circular hole through which the cannula tube can be passed
- ▶ Size adjusted to the anatomy of neonates and small infants
- ▶ **UNO** variants single-layered, for patients with little production of secretion
- ▶ **DUO** variants with two layers, for patients with higher production of secretion. The layer towards the skin – marked by a green thread - effectively absorbs secretions and transfers them to the peripheral layer. This keeps the skin dry and warm.
- ▶ Standard variants without slit
- ▶ **SLIT** variants with slit for exchange while the tracheostomy tube is in situ
- ▶ **ALU** variants with aluminium coating to reduce adhesion of the compress to skin around the tracheostoma



Order Information:

SENSOTRACH® UNO PED , 6 x 7 cm, pack of 10	REF 30522
SENSOTRACH® UNO PED SLIT , 6 x 7 cm, pack of 10	REF 30523
SENSOTRACH® UNO PED SLIT Junior , 6 x 7 cm, pack of 10	REF 30525
SENSOTRACH® UNO PED ALU SLIT , 6 x 7 cm, pack of 10	REF 30524
SENSOTRACH® DUO PED , 6 x 7 cm, pack of 10	REF 30622
SENSOTRACH® DUO PED SLIT , 6 x 7 cm, pack of 10	REF 30623