


# DURAVENT® UNI LINGO-PHON

1 ICU + 1 ICF



XL

- ▶ Tracheostomy tube for long-term use. Can be used up to 3-6 months
- ▶ With two inner cannulas (IC): One fenestrated with silver speaking valve, one closed with 15 mm connector (ICU / UNI) to connect an HME
- ▶ Easy-Lock has been designed for fixing the inner cannula and compatible accessories securely - no twisting
- ▶ Outer cannula (AK) with sieve-fenestration
- ▶ Standard length from size 3.5 to 13
- ▶ Extra-long (XL) variants from size 05 to 13
- ▶ Soft, pliable, lightweight
- ▶ Transparent, inconspicuous design
- ▶ Custom-made variants possible (length. sieve-fenestration. bending)
- ▶ Including broad adjustable OPTIFLAUSCH® tube holder and a product ID card

## Order information:

**DURAVENT® UNI LINGO-PHON 1 ICU + 1 ICF, SIZE 3.5-13**



**REF 21042**

**DURAVENT® UNI LINGO-PHON XL 1 ICU + 1 ICF, SIZE 5-13**



**REF 21242**

When ordering, please state the size after the item number!

☀ Keep away from sunlight

📖 Single patient use

📖 Consult instructions for use

⚠ CAUTION:  
United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner (in the US only).

📄 Rx only

CE 0482 CE marking identification number of the notified body

**Product made in Germany**

Size	O.D. neck flange mm	O.D. Tip mm	I.D. mm O-cannula	I.D. mm I-cannula	Length mm standard	Length mm extra-long (XL)
3.5	6.2	4.8	3.6	2.7	55	-
4	7.2	5.8	4.4	3.3	55	-
5	8.4	6.6	5.4	4.3	55	95
6	9.4	7.8	6.2	4.9	60	95
7	10.2	8.8	7.0	5.5	65	100
8	11.2	9.8	8.0	6.5	70	105
9	12.4	11.2	9.0	7.5	75	105
10	13.4	12.2	10.0	8.5	85	110
11	14.4	13.2	11.0	9.5	87	115
12	15.4	14.2	12.0	10.2	90	120
13	16.6	15.2	13.0	11.2	90	125

O.D. neck flange = Outer diameter behind the neck flange · O.D. Tip = Outer diameter of the outer cannula at the tip of the cannula · I.D. = Inner diameter at the tip of the cannula · O-cannula = Outer cannula · I-cannula = Inner cannula · Length = Length over the outer curve