


HUMIDOSTOM® JUNIOR

REF 46850

PED

UNI



- ▶ HME filter (Heat-Moisture Exchanger) for children
- ▶ 15 mm-connector (UNI) fits to all tracheostomy tubes with 15 mm standard connector
- ▶ Open pore polyurethane foam filter
- ▶ Filter made of blue colour for cosmetic/aesthetic reasons
- ▶ Lightweight, only 2 g weight, size adjusted for pediatric anatomy
- ▶ Individually packed in convenient peel-packs
- ▶ Pack size of 30 each

SCOPE OF DELIVERY

- ▶ Instruction for use

HUMIDOSTOM® JUNIOR**REF 46850**

1 x HME filter cassettes

Product	Unit	REF
HUMIDOSTOM® JUNIOR	30 pieces, individually packed	46850

